

# Axion



## ROTARY INTERFACE

NEW ZEALAND  
SCRIPTING AND ORDERING GUIDE

# Client Details

Client/User Name:	
Date:	
NDIS/Funding:	





## Notes

**Further Information**

## Part One: Select Axion with either Whitmyer Pad, Spex Pad or Savant Pad

# Mounting Axion to Whitmyer Pad

### Step One: Select Whitmyer Head Support Pad


Cradle	Plush	Adjust-a-Plush	Complex Options
			
<b>Small</b> 1480-4211-281	<b>10" Pad</b> 1480-4110-281	<b>Small</b> 1480-4310-281	Refer to your Medifab rep!
<b>Medium</b> 1480-4221-281	<b>14" Pad</b> 1480-4114-281	<b>Medium</b> 1480-4320-281	
<b>Large</b> 1480-4231-281	<b>19" Pad</b> 1480-4119-281	<b>Large</b> 1480-4330-281	

### Step Two: Select Axion to fit Whitmyer Head Support Pad


<b>Axion Interface 60° rotation (19mm ball)</b> 1280-4461-100		<b>Axion Interface 90° rotation (19mm ball)</b> 1280-4491-100	
--	---	--	---

# Mounting Axion to Spex Pad


## Step One: Select Spex Head Support Pad (Includes Ball)

Comfort	Circle	Contour	Standard and Extended Lateral	Adjustable Lateral
				
<b>250</b> 1275-0031-000	<b>150</b> 1275-0011-000	<b>250</b> 1275-0021-000	<b>Standard Paediatric</b> 1275-0041-000	<b>150</b> 1275-0061-000
<b>350</b> 1275-0032-000	<b>250</b> 1275-0012-000	<b>350</b> 1275-0022-000	<b>Standard</b> 1275-0043-000	<b>250</b> 1275-0062-000
			<b>Extended Paediatric</b> 1275-0042-000	
			<b>Extended</b> 1275-0044-000	

## Step Two: Select Axion to fit Spex Head Support Pad

<b>Axion Interface 60° rotation (25mm ball)</b> 1280-4432-100		<b>Axion Interface 90° rotation (25mm ball)</b> 1280-4492-100	
--	---	--	---




## Step 3: Add Clamp Ball link (required)

<b>Add a link to clamp ball to ball (25mm)</b> Must be selected when using Spex Pad in conjunction with Axion 1280-4685-100	
---	---


# Part One: Select Axion with either Whitmyer Pad, Spex Pad or Savant Pad

## Mounting Axion to Savant Pad

### Step One: Select Savant Head Support Pad

		
<b>Paediatric</b> Circumference at forehead: 14" - 18" (36-46cm) 1280-4101-000	<b>Small</b> Circumference at forehead: 18" - 22" (46-56cm) 1280-4102-000	<b>Adult</b> Circumference at forehead: 22" - 24" (56-61cm) 1280-4104-000

### Step Two: Select Axion to fit Savant Head Support Pad

<b>Axion Interface 60° rotation (25mm ball)</b> 1280-4432-100		<b>Axion Interface 90° rotation (25mm ball)</b> 1280-4492-100	
--	--	--	--





# Part Two: Select either Whitmyer or Spex Stylo Headrest Stem

## Mounting to Whitmyer Stem

### Step One: Select Whitmyer Stem

Cobra Flip Back	Linx <sup>2</sup>	Linx	Axys	Pro
				
1480-4927-100	<b>4"</b> 1480-4962-004	<b>4"</b> 1480-4961-004	1480-4913-100	<b>Paediatric</b> 1480-4934-101
	<b>8"</b> 1480-4962-008	<b>8"</b> 1480-4961-008		<b>Adult</b> 1480-4933-100
	<b>10"</b> 1480-4962-010	<b>10"</b> 1480-4961-010		
	<b>15"</b> 1480-4962-015	<b>15"</b> 1480-4961-015		

### Step Two: Select Detach Bracket for Whitmyer Stem

Cobra Flipback	Linx/Linx <sup>2</sup>	Axys	Pro
			
1480-4921-100	1481-M2300	1480-4911-100	1480-4931-100

### Step Three: Add Adaptor Ring (required)

#### Savant Adapter ring for Whitmyer 19mm ball mount

Must be selected when ordering Whitmyer Stem in conjunction with Axion.

1280-4680-100



## Mounting to Stylo Stem

### Step One: Select Stylo Stem (includes ball joint mount)

Stylo 160	Stylo 260	Stylo Flipback 160	Stylo Flipback 260
			
1274-1252-000-3	1274-1282-000-3	1274-2252-000-3	1274-2282-000-3

# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script Axion\_20240412\_NZ