

BINGO Evo



PAEDIATRIC ADAPTIVE STROLLER

Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

July 2024

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

	Size 1	Size 2	Size 3
BINGO Stroller Package	2881-0015-000	2882-0015-000	2883-0008-000
BINGO Switch Stroller Package	2884-0001-000	2884-0002-000	N/A

Mobility Base Specifications			
Push handle height (cm)	75-110	75-110	
Overall width (cm)	63	68	
Wheel diameter, front/swivel (cm)	19	19	
Wheel diameter, rear (cm)	29	29	
Tilt in space (Forward and rear facing)	45°	45° (Size 3 - Forward facing only)	
Folded size, mobility base (cm) (with rear wheels on)	80 x 63 x 38	80 x 68 x 38	
Folded size, mobility base (cm) (with rear wheels removed)	80 x 63 x 31	80 x 68 x 31	
Weight, mobility base (kg) (swivelling front wheels)	8.6	8.7	
Load capacity, mobility base (kg)	60	60	

Seat Unit Specifications			
Seat width (cm)	18-30	23-35	23-35
Seat depth (cm)	16-34	22-40	30-45
Lower leg length (cm)	19-40	19-40	19-40
Back support height (cm)	45-63	54-75	62-83
Height of slots for shoulder straps (cm)	28-46	34-55	43-64
Seat to floor height (cm) (dependant on tilt position of seat)	46	46	46
Back support recline	90-160°	90-160°	90-135°
Weight of seat (kg)	6.4	7.4	8.4
Load capacity of seat unit (kg)	50	50	50
Weight of seat and mobility base combined (kg)	15	16.1	17.1




Stroller packages include:

- Matte black stroller frame with black wheels (BINGO Switch frame features interchangeable front end capabilities for mounting the Swivelling Front Wheel Attachment - page 6)
- Sun canopy and rain cover
- Removable seat unit
- Tie-downs for transport
- Removable hip supports (width & depth adjustable)
- Fabric bag (max. 3kg)
- Adjustable suspension set
- Receiver for grab rail or therapy tray
- Elevating foot support hanger (angle adjustable 90-160°)

Note: Upholstery not included in package

Items outlined in blue are required as standard package

Upholstery Colour Options

Berry Soft		Turquoise Soft		Black	
Upholstery (required)					
Size 1 2881-8801-347		Size 1 2881-8801-349		Size 1 2881-8801-346	
Size 2 2882-8802-347		Size 2 2882-8802-349		Size 2 2882-8802-346	
Size 3 2883-8803-347		Size 3 2883-8803-349		Size 3 2883-8803-346	



Hip Support Upholstery (required)					
Size 1 2881-1075-347		Size 1 2881-1075-349		Size 1 2881-1075-346	
Size 2 2882-1076-347		Size 2 2882-1076-349		Size 2 2882-1076-346	
Size 3 2883-1077-347		Size 3 2883-1077-349		Size 3 2883-1077-346	

Tubie Fun Covers

Cover Sets

Cover set includes seat liner with waterproof layer along with covers for hip support, shoulder wing and spacer fabric head support pad.


(Note: Spacer fabric head support pads required. Must be used in conjunction with standard covers selected on page 2).

Llamas Cover Set - For Size 1 BINGO 2880-9561-121	
Llamas Cover Set - For Size 2 BINGO 2880-9562-121	
Llamas 5-Point Harness Cover (One size) 2880-6511-121	
Planets Cover Set - For Size 1 BINGO 2880-9561-036	
Planets Cover Set - For Size 2 BINGO 2880-9562-036	
Planets 5-Point Harness Cover (One size) 2880-6511-036	




Head Support Options

Spacer Fabric Head Support Pads

Size 1 11cm (D) x 16cm (H) 2880-4471-346	
Size 2 13cm (D) x 18cm (H) 2880-4472-346	

Head Support Speakers

Requires Spacer Fabric Head Support Pads.
2880-8330-000






Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
Merino/Spacer Fabric 2010-4110-099		Small (17cm Internal width)		Small (17cm Internal width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
		Wipe Down Cover 2010-4210-209		Wipe Down Cover 2010-4310-209	
Wipe Down Cover 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)	
		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	

Lateral Support Options

Standard Lateral Supports

Height and angle adjustable. All sizes are compatible with all BINGO seat units


Size 1 8cm (H) x 10cm (L) 2880-3010-346	Size 2 10cm (H) x 12cm (L) 2880-3020-346	Size 3 12cm (H) x 13cm (L) 2880-3030-346	
--	---	---	---

Narrow Straight Lateral Supports		Broad Lateral Supports		Teardrop Lateral Supports	
Size 0 8.5cm (H) x 11.5cm (D) 2880-3008-100		Size 1 10.5cm (H) x 13cm (D) 2880-3008-310		Size 0 15cm (H) x 10.5cm (D) 2880-3008-500	

Pelvis and Lower Limb Support Options

Sacral Support (including spacer)


22cm (W) x 9cm (H) x 1cm (D)
2880-2510-060



Pelvis and Lower Limb Support Options


Thigh Wedge

22.5cm (W) x 10.5cm (D)
Height: 2cm tapering to 1cm.
(Recommended for size 1 BINGO Evo only).
2006-0505-000



Trochanter Wedges

(Recommended for size 1 BINGO Evo only).
2006-0660-000



Swing-Down Abductor Pommel (Fits all sizes)

Size 1
9cm (H) x 6cm (D)
2880-1651-346

Size 2
1cm (H) x 8cm (D)
2880-1652-346

Size 3
12.5cm (H) x 9cm (D)
2880-1653-346



Footplate Spacer

Spacer block to reduce lower leg length.
4cm thickness.

Size 1
2880-5321-000

Size 2 & 3
2880-5322-000



Ankle Footcasts

Requires technician to install. Purchase only.
Inner circumference listed below:

Size 1
16-25cm
2880-6101-000

Size 2
21-29cm
2880-6102-000



Hip Belts & Harnesses

Spex 2-Point Hip Belt

X-small
2 point slim line, dual pull
with neoprene sleeve padding.
1405-6220-027



Wonderseat 4-Point Hip Belt

Stabilises and controls the
pelvis. Can be quickly removed
for washing.
1406-5420-027



Groin Strap

To fit Size 1
2880-6361-346

To fit Size 2 & 3
2880-6362-346



Wonderseat Padded Centrepoint Shoulder Harness

Size 0
Pad height 22cm
1409-5660-017

Size 1
Pad height 28cm
1409-5661-017



Hip Belts & Harnesses

5-Point Harness

Includes front pull adjusters and neck padding.

2880-6511-346



Spex 5-Point Harness

Includes tee pad for babies and very small children. (Recommended for size 1 BINGO Evo only).

2006-6551-000



Tray & Grab Rail Options

Therapy Tray

Cannot be used with Grab Rail

To fit Size 1

2880-7111-000

To fit Size 2 & 3

2880-7112-000



Therapy Tray Infill (Fits all sizes)

For smaller trunk cut out.
2880-7180-000



Grab Rail

Cannot be used with therapy tray

To fit Size 1

2880-7251-346

To fit Size 2 & 3

2880-7252-346



Critical Care Options

IV Pole (Fits all sizes)

With mobility bracket for attachment to mobility base.

2880-8275-000

Extra Bracket

For QR mounting of IV pole to second base

2880-8275-050



Oxygen Bottle Holder (Fits all sizes)

With mounting bracket for attachment to mobility base. 120mm diameter for 'C Type' bottle.

2880-8253-000




Chassis Mount Tray

Size 1

2880-8058-000

Size 2 & 3

2880-8059-000



Base Frame Options

Bingo Switch Swivelling Front Wheel Attachment

Suitable for use with Bingo Switch frame only.

Size 1
2884-9571-000

Size 2
2884-9572-000



Separate Seat Unit

Separate seat unit for use with the Cobra Base if a combined stroller and seat unit is not purchased. This is also a great option if two seat units are preferred to avoid swapping the seat unit between the stroller and Cobra bases. Must be purchased with upholstery.

Size 1
2881-0011-000

Size 2
2882-0012-000

Size 3
2883-0013-000



Cobra Indoor Base

Cobra Hi-Low Indoor Frame with 75mm Twin Castors

Size 1
2880-8711-000

Size 2 & 3
2880-8721-000



Twin Bingo Stroller

This twin stroller connection is only possible between strollers of the same size, or for connection of the mini and Size 1 together, it is not available for the size 3.

Please fill out an additional script form to order the second stroller unit with a unique configuration and accessories.

Note that this can only be fitted at the time of ordering both stroller units. (Medifab cannot retrofit this attachment).

2880-8630-000

Attendant Brake

Recommended for purchase with Twin Stroller configuration due to the increased combined weight and inertia when both strollers are connected.

2880-8452-000



medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: BINGO_20240726_NZ