

BINGO Mini



EARLY INTERVENTION ADAPTIVE STROLLER

Distributed by
medifab

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	

Get it Faster

9999-ASMB-002

Medifab Assembly Service for Easy Delivery

Medifab offers our popular assembly service on this product, meaning the product arrives ready to be delivered, straight from the box, saving time and speeding up the delivery process.

Product Details

Bingo Mini Package Options

Order Code, Size 0		2880-0008-000	
Conversion Kit from Mini to Size 1 (stroller and Cobra) Converts BINGO Mini frames to accept the BINGO Evo size 1 seat unit, extending the growth of the solution (size 1 seat unit must be purchased separately).		2880-8080-086	
Seat Unit		Mobility Base with swivelling wheels	
Seat width	13-25cm	Push handle height	75-110cm
Seat depth	12-27cm	Overall width	63cm
Lower leg length	16-36cm	Wheel diameter front (swivel)	19cm
Back support height	37-51cm	Wheel diameter rear	29cm
Height of slots for shoulder straps	25.5-39.5cm	Tilt in space (forward and rear facing)	45°
Seat to floor height (dependant on tilt position of seat)	50cm	Folded size mobility base (incl. rear wheels)	80 x 63 x 38cm
Backrest recline	80-180°	Folded size mobility base (excl. rear wheels)	80 x 63 x 31cm
Weight of seat unit*	7.8kg	Weight mobility base (swivelling front wheels)	9.8kg
Load capacity of seat unit	25kg	Load capacity mobility base	60kg
Weight of seat unit and mobility base combined*	17.6kg		

*Including upholstery and standard accessories




Stroller package includes:

- Matte black foldable stroller frame with black wheels
- Removable seat unit
- Tie-downs for transport
- Removable hip supports (width & depth adjustable)
- Fabric bag (max. 3kg)
- Adjustable suspension set
- Receiver for grab rail or therapy tray
- Elevating foot support hanger (angle adjustable 85-180°)
- **Note: Upholstery not included in package.**



Items outlined in blue are required as standard package

Upholstery Colour Options

Berry Soft		Turquoise Soft		Black	
Seat Upholstery (required)					
2880-8800-347		2880-8800-349		2880-8800-346	

Berry Soft		Turquoise Soft		Black	
Hip Support Upholstery (required)					
2880-1074-347		2880-1074-349		2880-1074-346	

Upholstery Upgrades

Upholstery Pressure Relief Foam Upgrades

Seat
2880-0391-000

Back
2880-2392-000

Head Supports

Head Support Pads

2880-4470-346



Head Supports (continued)

Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
Merino/Spacer Fabric 2010-4110-099		Small (17cm Internal width)		Small (17cm Internal width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
		Wipe Down Cover 2010-4210-209		Wipe Down Cover 2010-4310-209	
Wipe Down Cover 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)	
		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	


Lateral Supports

Standard Lateral Supports		Spex Narrow Straight Lateral Supports	
Size 0 80mm (H) x 100mm (D) 2880-3009-346		Size 0 85mm (H) x 115mm (D) 2880-3007-100	
Spex Broad Lateral Supports		Spex Teardrop Lateral Supports	
Size 1 105mm (H) x 130mm (D) 2880-3007-310		Size 0 150mm (H) x 105mm (D) 2880-3007-500	

Sacral Support

Sacral foam insert including Spacer	
On velcro attachment. Also includes foam spacer which adheres behind Sacral insert where required. 10cm (H) x 18cm (W) 2880-2510-000	

Seat Options

Thigh Wedge		Trochanter Wedges	
Velcro attached under upholstery. 22.5cm (W) x 10.5cm (D). Height: 2cm tapering down to 1cm 2006-0505-000		Option to contour standard flat cushion. Velcro attached under upholstery. 2006-0660-000	

Abductor Pommel

Swing down abductor pommel

Size 1

For increased abduction of the legs.
9cm (D) x 6.5cm (W) x 9cm (H)
2880-1651-346



Hip Belts & Harnesses

Spex 5-Point Harness

Must be worn when stroller is used for vehicle transportation.
2006-6551-000



Wonderseat Harness

Wonderseat Padded Centrepoint Shoulder Harness - Size 0.
Pad height 22cm
1409-5660-017



Wonderseat Hip Belt


Wonderseat 4-Point Padded Hip Belt.
Pad width 24cm.
1406-5420-027



Foot Positioning

Footplate Cover

Neoprene fabric protects feet from the hard footplate surface.
2880-5305-000



Footplate Spacer

Spacer block to reduce lower leg length by 4cm.
2880-5321-000



Ankle Footfasts

Strap inner circumference 15-23cm.
2880-6101-000



Tray & Grab Rail Options

Therapy Tray

Cannot be used at the same time as Grab Rail.
2880-7110-000




Therapy Tray Infill

For smaller trunk cut out.
2880-7180-000



Grab Rail

Cannot be used at the same time as Therapy Tray.
2880-7250-346



Critical Care Options

IV Pole

With mounting bracket for attachment to mobility base. Fits all sizes.

2880-8275-000



Extra Bracket

For QR mounting of IV pole to second base. (Cobra or other indoor base).

2880-8275-050

Oxygen Bottle Holder

With mounting bracket for attachment to mobility base. 120mm diameter for "C Type" bottle. Fits all sizes.


2880-8253-000



Chassis Mount Tray

Size 1

2880-8058-000



Swivel Lock

To hold castors straight when brisk walking

2880-8440-000



Hand Strap

Wrist safety tether strap

2005-6930-000



Canopy

Includes mini mounts and rain cover

2880-7501-346



Cobra Indoor Base

Cobra Hi-Low Indoor Frame with 75mm Twin Castors

The seat unit easily transfers out of the stroller base and into the Cobra indoor base for feeding and activity times. Lower down the floor for playing and activities, and raise up to a table for meal time! Tilt settings can still be adjusted in the same way as the stroller base.

Size 0
2880-8720-000



Separate Seat Unit

Use with Cobra Base

Separate seat unit for use with the Cobra Base if a combined stroller and seat unit is not purchased. This is also a great option if two seat units are preferred to avoid swapping the seat unit between the stroller and Cobra bases. Must be purchased with upholstery.

One Size
2880-0010-000



Twin Bingo Stroller

Connection mount that allows two strollers to be connected one in front of the other for an additional sibling. The strollers are able to be easily detached and used separately where required, providing flexibility when travelling with or without siblings.

This is not available for use with a size 3 Bingo Evo. This can be used to connect either: bases of the same size (and corresponding sized Bingo Spex strollers), or the Bingo Mini and Evo Size 1. It is not suitable for use between any other size combinations.

Please fill out an additional script form to order the 2nd stroller unit with a unique configuration and accessories.

Please note that this can only be fitted at the time of ordering both stroller units (Medifab cannot retrofit this attachment).

2880-8630-000



Attendant Brake

Recommended for purchase with Twin Stroller configuration due to the increased combined weight and inertia when both strollers are connected.

2880-8452-000

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: BINGO Mini_20241031_AU