

# Duro Stroller



## PAEDIATRIC ADAPTIVE STROLLER

Distributed by  
**medifab**

AUSTRALIAN  
SCRIPTING AND ORDERING GUIDE

# Client Details

Client/User Name:	
Date:	
NDIS/Funding:	

Get it Faster

9999-ASMB-002

## Medifab Assembly Service for Easy Delivery

Medifab offers our popular assembly service on this product, meaning the product arrives ready to be delivered, straight from the box, saving time and speeding up the delivery process.

# Package Options

Specifications	Size 1	Size 2
<b>Order Code</b>	2780-0001-052	2780-0002-052
<b>Seat Specifications</b>		
Seat width (cm)	23-34	28-39
Seat depth (cm)	26-40	31-45
Seat angle (cm)	16 or 20°	16 or 20°
Back height (cm)	45-63	61-83
Back angle	78-120°	78-120°
Knee angle	85-160°	85-160°
Lower leg length (cm)	20-44	20-44
Folded size (L x W x H) (cm)	84 x 63 x 32	84 x 68 x 38
<b>Base Specifications</b>		
Max weight (kg)	75	75
Front wheel diameter (cm)	20	20
Rear wheel diameter (cm)	30	30

### Stroller Package Includes:

- Swivelling Front Wheels
- Foot Support Bracket, Angle Adjustable, Flip Up
- Foot Support With Heel Stop
- Accessories Bag (Max Weight: 3kg)
- Tie Down Kit
- Retainer for Accessories
- Hip Guides (Low Version)
- **Note: Upholstery Not Included In Package**

# Upholstery Colour Options

Turquoise	Black
Seat Upholstery (required)	
<b>Size 1</b> 2780-8801-349	<b>Size 1</b> 2780-8801-346
<b>Size 2</b> 2780-8802-349	<b>Size 2</b> 2780-8802-346



# Head Supports

## Head Support Pads

### Size 1

11cm (D) x 16cm (H)  
2880-4471-346

### Size 2

13cm (D) x 18cm (H)  
2880-4472-346



Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
<b>Merino/Spacer Fabric</b> 2010-4110-099		<b>Small (17cm Internal width)</b>		<b>Small (17cm Internal width)</b>	
		<b>Merino/Spacer Fabric</b> 2010-4210-099		<b>Merino/Spacer Fabric</b> 2010-4310-099	
		<b>Wipe Down Cover</b> 2010-4210-209		<b>Wipe Down Cover</b> 2010-4310-209	
<b>Wipe Down Cover</b> 2010-4110-209		<b>Large (21cm Internal width)</b>		<b>Large (21cm Internal width)</b>	
		<b>Merino/Spacer Fabric</b> 2010-4211-099		<b>Merino/Spacer Fabric</b> 2010-4311-099	
		<b>Wipe Down Cover</b> 2010-4211-209		<b>Wipe Down Cover</b> 2010-4311-209	

# Lateral Supports

## Standard Lateral Supports

Height and angle adjustable. All sizes are compatible with both Duro sizes.

### Size 1

8cm (H) x 10cm (L)  
2880-3010-346

### Size 2

10cm (H) x 12cm (L)  
2880-3020-346

### Size 3

12cm (H) x 13cm (L)  
2880-3030-346



## Teardrop Lateral Supports

### Size 1

13.5cm (H) x 7.5cm (L)  
2006-3110-099





## Contoured Lateral Supports



Spex Discovery Axial Design  
9cm (H) x 13cm (L)  
2006-3030-099





# Harness Options


Spex Centrepoint Harness		5-Point Harness	
<p><b>X-Small</b> Recommended for size 1 Duro Total pad dimensions: 26cm (H) x 22cm (W) 2006-6660-017</p>		<p>With front pull adjusters and neck padding 2880-6511-346</p>	
<p><b>Small</b> Recommended for size 1 &amp; 2 Duro Total pad dimensions: 36cm (H) x 32cm (W) 1409-6661-017</p>			
<p><b>Medium</b> Recommended for size 2 Duro Total pad dimensions: 36cm (H) x 32cm (W) 1409-6662-017</p>			

# Hip Belt Options

Spex Push Button		Pen Popper	
<p><b>X-Small</b> Pad width 30cm 1405-6210-028</p>		<p><b>X-Small</b> Pad width 30cm 1405-6210-025</p>	
<p><b>Small</b> Pad width 36cm 1405-6211-028</p>		<p><b>Small</b> Pad width 36cm 1405-6211-025</p>	
<p><b>Medium</b> Pad width 46cm 1405-6212-028</p>		<p><b>Medium</b> Pad width 46cm 1405-6212-025</p>	

# Tray & Grab Rail Options

Therapy Tray		Therapy Tray Infill	
Cannot be used with Grab Rail			<p>For smaller trunk cut out. Fits both tray sizes. 2880-7180-000</p> 
<p><b>To fit Size 1</b> 2880-7111-000</p>	<p><b>To fit Size 2</b> 2880-7112-000</p>		

Grab Rail	
Cannot be used with therapy tray	
<p><b>To fit Size 1</b> 2880-7251-346</p>	<p><b>To fit Size 2</b> 2880-7252-346</p> 

# Hand Strap

Wrist safety tether strap

2005-6930-000



# Swivel Lock

To hold castors straight when brisk walking

2880-8440-000



# Canopy

Includes mounts and rain cover

**Size 1**  
2780-7501-346

**Size 2**  
2780-7502-346



# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Duro\_20241029\_AU