

Duro Stroller



PAEDIATRIC ADAPTIVE STROLLER

Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

October 2024

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	

Notes

Package Options

Specifications	Size 1	Size 2
Turquoise	2780-0011-349	2780-0012-349
Black	2780-0011-346	2780-0012-346

Seat Specifications		
Seat width (cm)	23-34	28-39
Seat depth (cm)	26-40	31-45
Seat angle (cm)	16 or 20°	16 or 20°
Back height (cm)	45-63	61-83
Back angle	78-120°	78-120°
Knee angle	85-160°	85-160°
Lower leg length (cm)	20-44	20-44
Folded size (L x W x H) (cm)	84 x 63 x 32	84 x 68 x 38

Base Specifications		
Max weight (kg)	75	75
Front wheel diameter (cm)	20	20
Rear wheel diameter (cm)	30	30

Stroller package includes:

- Hoggi Duro Rehab-Buggy
- Hip Support, Aluminium
- Seat and Backrest Upholsteries
- Small Spex Hip Belt - Side-Release Buckle
- Canopy
- PU Wheels with Black Rim
- Swivelling Front Wheels
- Supporting Brackets
- Foot Support Bracket, Angle Adjustable, Flip Up
- Foot Support with Heel Stop
- Accessories Bag (Max weight: 3kg)
- Tie Down Kit (According to ISO 7176-19)
- Retainer for Accessories
- User Manual



Head Supports

Head Support Pads

Size 1

11cm (D) x 16cm (H)
2880-4471-346

Size 2

13cm (D) x 18cm (H)
2880-4472-346



Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
Merino/Spacer Fabric 2010-4110-099		Small (17cm Internal width)		Small (17cm Internal width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
		Wipe Down Cover 2010-4210-209		Wipe Down Cover 2010-4310-209	
Wipe Down Cover 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)	
		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	

Lateral Supports

Standard Lateral Supports

Height and angle adjustable. All sizes are compatible with both Duro sizes.

Size 1

8cm (H) x 10cm (L)
2880-3010-346

Size 2

10cm (H) x 12cm (L)
2880-3020-346

Size 3

12cm (H) x 13cm (L)
2880-3030-346



Teardrop Lateral Supports

Size 1

13.5cm (H)
x 7.5cm (L)
2006-3110-099

Size 2

15.5cm (H)
x 10.5cm (L)
2006-3120-099





Contoured Lateral Supports



Spex Discovery Axial
design 9cm (H) x 13cm (L)
2006-3030-099





Harness Options


Spex Centrepoint Harness		5-Point Harness	
<p>X-Small Recommended for size 1 Duro Total pad dimensions: 26cm (H) x 22cm (W) 2006-6660-017</p>		<p>With front pull adjusters and neck padding 2880-6511-346</p>	
<p>Small Recommended for size 1 & 2 Duro Total pad dimensions: 36cm (H) x 32cm (W) 1409-6661-017</p>			
<p>Medium Recommended for size 2 Duro Total pad dimensions: 36cm (H) x 32cm (W) 1409-6662-017</p>			

Hip Belt Options

Spex Push Button		Pen Popper	
<p>X-Small Pad width 30cm 1405-6210-028</p>		<p>X-Small Pad width 30cm 1405-6210-025</p>	
<p>Small Pad width 36cm 1405-6211-028</p>		<p>Small Pad width 36cm 1405-6211-025</p>	
<p>Medium Pad width 46cm 1405-6212-028</p>		<p>Medium Pad width 46cm 1405-6212-025</p>	

Tray & Grab Rail Options

Therapy Tray		Therapy Tray Infill	
Cannot be used with Grab Rail			<p>For smaller trunk cut out. Fits both tray sizes. 2880-7180-000</p> 
<p>To fit Size 1 2880-7111-000</p>	<p>To fit Size 2 2880-7112-000</p>		

Grab Rail	
Cannot be used with therapy tray	
<p>To fit Size 1 2880-7251-346</p>	<p>To fit Size 2 2880-7252-346</p> 

Hand Strap

Wrist safety tether strap

2005-6930-000



Swivel Lock

To hold castors straight when brisk walking

2880-8440-000



medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Duro_20241029_NZ