

# Jenx Atom



## MODULAR SEATING SYSTEM

Distributed by  
**medifab**

SCRIPTING AND ORDERING GUIDE

October 2024

# Client Details


Client/User Name: \_\_\_\_\_

Date: \_\_\_\_\_

NDIS/Funding: \_\_\_\_\_

Measurements	User Measurements
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	

Notes



## Product Specifications

Average age	6 Months-5 Years
Weight capacity	25kg
Seat width	18cm-27cm
Seat depth	9cm-34cm
Chest support width	20cm-28cm
Seat to top of backpad height	27cm-38cm
Foot support angle	-20° to 15°
Foot support height	11cm-36cm
Back support angle	-75° to 105°
Tilt range	30° posterior, 10° anterior



Therapists Normally Choose

## Seat Unit/Colour Choice

## Base Options

Red	Pink	Y-Base
-----	------	--------

**Includes:**

- Padded seat with memory foam
- Shaped back pad
- Arm supports
- Nursery tray
- Adjustable flip away hip support pads
- Push handle
- Upper back pad
- Flip up angle adjustable footplate


3431-1001-022

**Includes:**

- Padded seat with memory foam
- Shaped back pad
- Adjustable flip away hip support pads
- Push handle
- Upper back pad
- Arm supports
- Nursery tray
- Shoulder protractors
- High rise flexi thoracic pads
- Thoracic pad covers
- 2 Link multigrip head support
- Flip up angle adjustable footplate

3431-0010-028


Black Y-base with lockable swivel castors  
 Floor to seat height: 30-75cm  
 Tilt range: 30° posterior 10° anterior  
 3615-8740-000





# Head Supports

1 Link Multigrip	2 Link Multigrip	
Single link wings for some degree of lateral support. 33cm (W) x 5cm (D) x 7.5cm (H) <b>For Red Atom:</b> 3431-4601-022 <b>For Pink Atom:</b> 3431-4601-028	2 Links in each wing for extra support and positioning. 47.5cm (W) x 5cm (D) x 12cm (H)  <b>For Red Atom:</b> 3431-4602-022 <b>For Pink Atom:</b> 3431-4602-028	
Wipeable Cover (Optional) 5356-4671-000	Wipeable Cover (Optional) 5356-4672-000	

# Back Pad Options





Lumbar Pad	Back Pad Infill
Provides additional lumbar support. 3431-2501-000	
	A 50mm back pad to reduce the seat depth for smaller children. 3431-2370-000
	

Washable Cover	
For Seat and Back Support. Comes in Red. 3431-2470-022	

# Shoulder & Trunk Supports

Shoulder Protractors	High-Rise Short Thoracic Pad	High-Rise Flexible Thoracic Pads
Flip away when not required. Independently adjustable for width and angle. 13cm (H) x 15.5cm (D) 3431-2750-022	Short, slim trunk support pads with single-handed flip away action. No strap. 3430-3151-000	Slimline and super flexible trunk support pads with straps and single-handed flip away action. 8cm (H) x 18cm (D)  3431-3251-022
		
	Cover (optional) 3615-3171-022	Cover (required) 3431-3271-022

# Harnesses & Hip Belt

Waistcoat Harness	Wonderseat Harness	3 Point Lap Strap	4-Point Hip Belt
<p>Provides additional trunk support to children with low tone. 3615-6701-000</p>	<p>Size 0, Pad height 22cm Padded Centrepoint Shoulder Harness. 1409-5660-017</p>	<p>Dual lapstrap with additional safety of a groin harness. 3431-6251-000</p>	<p>Wonderseat Padded Hip Belt, Pad width 24cm. 1406-5420-027</p>
			

# Leg Supports

Leg Positioning Pads	Long Leg Sitting Board	Abduction Block	Abduction Block
<p>Ideal for controlling windsweeping or creating an abducted seating position. Required fixed angle foot support. 3431-1751-000</p>	<p>For very young children or those needing hamstring stretch. Used in place of foot support. Comes in red. 3431-0281-022</p>	<p><b>Size 0</b> Length: 9cm Height above Seat: 8.5cm Overall height: 14.5cm Width Front: 9cm Width Rear: 7cm 3615-1650-000</p>	<p><b>Size 1</b> Length: 18cm Height above Seat: 7.5cm Overall height: 15.5cm Width Front: 10cm Width Rear: 8.5cm 3615-1651-000</p>
			

# Foot Supports

Foot Support Pad
<p>Attaches to foot support for softer surface. Can only be used on the flip-up foot support. 3cm thickness. 3431-5701-000</p>


# Foot Supports (cont.)

Sandals	Sandal Raising Blocks	Ankle Snugs
<p>Ideal for positioning feet where a firmer hold on the foot is required. 8cm (W) x 17cm (L) 5356-6081-000</p>	<p>Easily accommodate leg length discrepancies by giving sandals extra height where needed. 4 x 1.5cm Blocks 5356-6091-000</p>	<p>Features unique magnetic buckle fastening mechanism. Sandals required. Inner Circumference: 12cm-22.5cm 5356-6131-000</p>

# Arm Supports

## Covers for Arm Supports

(Red) Provide soft top for arm supports when tray is not used.  
3431-5510-022



# Tray & Accessories

## Nursery Tray Infill


Designed to work with the Nursery Tray to reduce the cut-out for smaller children.  
3431-7182-000



## Nursery Tray Padding

Half Padding  
3431-7442-000

Full Padding  
3431-7452-000



## Elbow Blocks

Comes as a pair. Mounts to Nursery Tray corners for elbow support.  
9.5cm (H) x 10cm (W)  
3430-7350-000



## Activity Frame

Great for encouraging head control and hand/eye co-ordination. Requires tray.  
54cm (W) x 16cm (D) x 35cm (H)  
3615-7475-000



## Grab Bar

Nursery Tray required.  
Bar dimensions: 54cm (W) x 9.5cm (D) x 7cm (H)  
3615-7200-000



# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Jenx Atom\_20241003