

Jenx Multistander



**PRONE OR SUPINE
STANDING AID**


Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

December 2023

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	



Basic Packages

Multistander Size 1, Lime (Prone)	Multistander Size 1, Lime (Supine)
Order Code 5230-0011-054	Order Code 5230-0013-054
<p>Option one, for prone set-up.</p> <p>Includes:</p> <ul style="list-style-type: none"> • Multistander size 1 • Manual frame • Size 1 sandals • Size 1 knee blocks • Small hip support pads • Small trunk support pads and tray • MUS01-P1-03 (add headrest of choice to convert to supine) 	<p>Option one, for supine set-up.</p> <p>Includes:</p> <ul style="list-style-type: none"> • Multistander size 1 • Manual frame • Size 1 sandals • Size 1 knee blocks • Small hip support pads • Small trunk support pads • Tray and 1 link head support
Multistander Size 2, Lime (Prone)	Multistander Size 2, Lime (Supine)
Order Code 5230-0022-054	Order Code 5230-0021-054
<p>Option two, for prone set-up.</p> <p>Includes:</p> <ul style="list-style-type: none"> • Multistander size 2 • Manual frame • Size 2 sandals • Size 2 knee blocks • Medium hip support • Medium trunk support and tray 	<p>Option two, for supine set-up.</p> <p>Includes:</p> <ul style="list-style-type: none"> • Multistander size 2 • Manual frame • Size 2 sandals • Size 2 knee blocks • Slim thoracic • Medium hips supports • Moulded head rest • Multigrip body support and tray • MUS02-P1-03

Product Details

Product Specifications	Size 1	Size 1 - With Abduction	Size 2
Average age	9 Months – 6 Years	9 Months – 6 Years	3 – 13 Years
Weight capacity (kg)	35	35	60
Hip width (cm)	14 – 27.5	14 – 27.5	20 – 32
Hip height (cm)	25 – 44	25 – 44	37.5 – 60
Chest width (cm)	14 – 27.5	14 – 27.5	20 – 32
Chest height (cm)	49 – 85	50 – 89	75 – 120
Height range (cm)	62 – 116	60 – 116	110 – 155
Tilt range (°) (90 = upright)	10 – 90	10 – 90	0 – 85 (Powered) / 10 – 85 (Manual)
Frame width (cm)	61.5	61.5	67.5
Frame length (cm)	76	76	107
Abduction range (°)	-	30 (Either Leg)	-
Height of surface when horizontal (cm)	56	56	77.5

*Max ranges are achieved with additional accessories

Key

Blue	Suitable for Prone and Supine
Green	Suitable for Supine only
Yellow	Suitable for Prone only





Items outlined in blue are recommended as part of standard package


Power & Size Options

Manual		Battery Powered	
Size 1 5230-0001-000		Size 2 5230-0008-000	
Size 2 5230-0002-000			


Colour Options

Blue	Lime
-066	-054
	

Abduction Accessory

Only Available on Size 1	
Supports 0-30 degrees in each leg in 5 degree increments. Fitted prior to shipping. 5230-5000-000	



Support Boards

Size 1	Size 2	
PU chest & hip support boards. 450mm Wide 5230-2201-XXX	PU chest & hip support boards. 520mm Wide 5230-2202-XXX	

Head Supports



Moulded Head Support	Oval Head Support
Moulded polyurethane headrest for children who require minimal to moderate head support. 280mm wide x 210mm high 5230-4101-XXX	A simple support for children who require minimal head positioning. 190mm wide x 100mm high 5356-4611-XXX
	

Key Blue Prone and Supine Green Supine Only Yellow Prone Only

Multigrip - 1 Link	Multigrip - 2 Link	Multigrip - 3 Link	
Single link wings for some degree of lateral support. 33cm W x 5cm D x 7.5cm H 5356-4601-000	2 Links in each wing for extra support and positioning. 47.5cm W x 5cm D x 12cm H 5356-4602-000 	3 adjustable links in each wing for ultimate lateral headrest control and comfort. 62cm W x 5cm D x 12cm H 5356-4603-000	
Wipeable Cover (Optional) 5356-4671-000	Wipeable Cover (Optional) 5356-4672-000	Wipeable Cover (Optional) 5356-4673-000	

Pelvic Supports

Support Pad & Straps Small	Support Pad & Straps Medium	
A 100mm B 90mm C 90mm D 120mm E 25mm 5356-3221-XXX	A 100mm B 130mm C 90mm D 140mm E 25mm 5356-3222-XXX	


Flexi Thoracics (Pair) Size 1	Flexi Thoracics (Pair) Size 1 Long	Flexi Thoracics (Pair) Size 2	
A 180mm B 80mm 5356-3333-000	A 215mm B 80mm 5356-3334-000 	A 220mm B 110mm 5356-3335-000	
Covers Required 5356-3363-XXX	Covers Required 5356-3364-XXX	Covers Required 5356-3365-XXX	

Hi-Rise Thoracic Laterals

Adjustable for height, width and angle.

A 90mm B 90mm
C 200mm D 120mm
E 25mm

5230-3251-XXX



Trunk Supports

Support Pad & Straps Small	Support Pad & Straps Medium	
A 100mm B 90mm C 90mm D 120mm E 25mm 5356-3221-XXX	A 100mm B 130mm C 90mm D 140mm E 25mm 5356-3222-XXX	

Trunk Supports (cont.)

Flexi Thoracics (Pair) Size 1	Flexi Thoracics (Pair) Size 1 Long	Flexi Thoracics (Pair) Size 2	
A 180mm B 80mm 5356-3333-000	A 215mm B 80mm 5356-3334-000	A 220mm B 110mm 5356-3335-000	
Covers Required 5356-3363-XXX	Covers Required 5356-3364-XXX	Covers Required 5356-3365-XXX	

Hi-Rise Thoracic Laterals	Multigrip Body Support
----------------------------------	-------------------------------


Adjustable for height, width and angle. A 90mm B 90mm C 200mm D 120mm E 25mm 5230-3251-XXX		For Size 2 only. Provides extra support on shoulders or lumbar area. 420mm Wide 5230-2721-XXX	
---	--	--	--

Sternal Prompt

For Prone 5230-2910-000	
Cover required 5230-8801-XXX	

Waistcoat Harness

Size 2	Size 3	Size 4	Mounting Kit
Provides additional trunk and shoulder stability			Mounting bracket for Harness (Required for Size 1 Multistander with Abduction Only). 5230-part-056
A 5cm B 32cm C 14cm D 15.5cm 5230-6712-000	A 6.5cm B 37cm C 17.5cm D 18cm 5230-6713-000	A 8.5cm B 39.5cm C 20.5cm D 18.5cm 5230-6714-000	



Knee Supports

Size 1	Size 2	
Easy to attach and adjust for height, angle & depth		
A 100mm / B 45mm / C 75mm 5356-1751-XXX	A 120mm / B 60mm / C 95mm 5356-1752-XXX	


Tray Options

Elbow Blocks	
<p>Suits all sizes. 10cm W x 14cm H x 8.5cm D 5356-7351-000</p>	

Adjustable Tray - Size 1	Adjustable Tray - Size 2	
<p>A 650mm B 450mm C 370mm D (With Insert) 330mm / (Without Insert) 290mm 5230-7111-000</p>	<p>A 700mm B 530mm C 460mm D (With Insert) 410mm / (Without Insert) 360mm 5230-7112-000</p>	

Foot Positioning

Sandals			
<p>Size 1 170mm long x 80mm wide 5356-6081-000</p>	<p>Size 2 220mm long x 95mm wide 5356-6082-000</p>	<p>Size 3 275mm long x 105mm wide 5356-6083-000</p>	

Sandal Raiser Blocks			
<p>Size 1 Kit of 4x15mm blocks 5356-6091-000</p>	<p>Size 2 Kit of 4x15mm blocks 5356-6092-000</p>	<p>Size 3 Kit of 4x15mm blocks 5356-6093-000</p>	

Ankle Snugs		
<p>Size 1 120mm - 225mm Circumference 120mm High 5356-6131-000</p>	<p>Size 2 165mm - 295mm Circumference 120mm High 5356-6132-000</p>	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Multistander_20231208_NZ