

Jenx Multistander



**PRONE OR SUPINE
STANDING AID**

Distributed by
medifab

**AUSTRALIAN
SCRIPTING AND ORDERING GUIDE**

October 2024

Client Details

Client/User Name:

Date:

NDIS/Funding:

User Measurements

A Hip width	<input type="text"/>
B Seat depth	<input type="text"/>
C Seat to top of shoulder	<input type="text"/>
D Seat to top of head	<input type="text"/>
E Chest width	<input type="text"/>
F Seat to axilla	<input type="text"/>
G Leg length	<input type="text"/>

Get it Faster


Size 1: 9999-ASMB-001

Size 1 with Abduction: 9999-ASMB-003

Medifab Assembly Service for Easy Delivery

Medifab offers our popular assembly service on this product, meaning the product arrives ready to be delivered, straight from the box, saving time and speeding up the delivery process.

(Does not include Size 2)



Notes

Product Details

Product Specifications	Size 1	Size 1 - With Abduction	Size 2
Average age	9 Months – 6 Years	9 Months – 6 Years	3 – 13 Years
Weight capacity (kg)	35	35	60
Hip width (cm)	14 – 27.5	14 – 27.5	20 – 32
Hip height (cm)	25 – 44	25 – 44	37.5 – 60
Chest width (cm)	14 – 27.5	14 – 27.5	20 – 32
Chest height (cm)	49 – 85	50 – 89	75 – 120
Height range (cm)	62 – 116	60 – 116	110 – 155
Tilt range (°) (90 = upright)	10 – 90	10 – 90	0 – 85 (Powered) / 10 – 85 (Manual)
Frame width (cm)	61.5	61.5	67.5
Frame length (cm)	76	76	107
Abduction range (°)	-	30 (Either Leg)	-
Height of surface when horizontal (cm)	56	56	77.5

*Max ranges are achieved with additional accessories

Jenx Multistander Package Includes:

- Multistander manual or powered frame
- Base with 4 lockable swivel castors
- Knee block mounting brackets
- Angle adjustable footplate

Key

Blue	Suitable for Prone and Supine
Green	Suitable for Supine only
Yellow	Suitable for Prone only





Items outlined in blue are recommended as part of standard package


Power & Size Options

Manual		Battery Powered	
Size 1 5230-0001-000		Size 2 5230-0008-000	
Size 2 5230-0002-000			


Colour Options

Blue		Lime	
-066		-054	



Abduction Accessory



Only Available on Size 1	
Supports 0-30° in each leg in 5° increments. Fitted prior to shipping. 5230-5000-000	

Support Boards

Size 1	Size 2	
PU chest & hip support boards. 450mm (W) 5230-2201-XXX	PU chest & hip support boards. 520mm (W) 5230-2202-XXX	



Head Supports

Moulded Head Support		Oval Head Support	
Moulded polyurethane headrest for children who require minimal to moderate head support. 280mm (W) x 210mm (H) 5230-4101-XXX		A simple support for children who require minimal head positioning. 190mm (W) x 100mm (H) 5356-4611-XXX	

Multigrip - 1 Link	Multigrip - 2 Link	Multigrip - 3 Link	
Single link wings for some degree of lateral support. 33cm (W) x 5cm (D) x 7.5cm (H) 5356-4601-000	2 Links in each wing for extra support and positioning. 47.5cm (W) x 5cm (D) x 12cm (H) 5356-4602-000 	3 adjustable links in each wing for ultimate lateral headrest control and comfort. 62cm (W) x 5cm (D) x 12cm (H) 5356-4603-000	
Wipeable Cover (optional) 5356-4671-000	Wipeable Cover (optional) 5356-4672-000	Wipeable Cover (optional) 5356-4673-000	

Pelvic Supports

Support Pad & Straps, Small	Support Pad & Straps, Medium	
A 100mm B 90mm C 90mm D 120mm E 25mm 5356-3221-XXX	A 100mm B 130mm C 90mm D 140mm E 25mm 5356-3222-XXX	



Flexi Thoracics (Pair)			
Size 1 A 180mm B 80mm 5356-3333-000	Size 1 Long  A 215mm B 80mm 5356-3334-000	Size 2 A 220mm B 110mm 5356-3335-000	
Covers (required) 5356-3363-XXX	Covers (required) 5356-3364-XXX	Covers (required) 5356-3365-XXX	



Hi-Rise Thoracic Laterals


Adjustable for height, width and angle. A 90mm B 90mm C 200mm D 120mm E 25mm 5230-3251-XXX	
--	---

Trunk Supports

Support Pad & Straps, Small	Support Pad & Straps, Medium	
A 100mm B 90mm C 90mm D 120mm E 25mm 5356-3221-XXX	A 100mm B 130mm C 90mm D 140mm E 25mm 5356-3222-XXX	

Flexi Thoracics (Pair)			
Size 1 A 180mm B 80mm 5356-3333-000	 Size 1 Long A 215mm B 80mm 5356-3334-000	Size 2 A 220mm B 110mm 5356-3335-000	
Covers (required) 5356-3363-XXX	Covers (required) 5356-3364-XXX	Covers (required) 5356-3365-XXX	

Hi-Rise Thoracic Laterals	Multigrip Body Support
Adjustable for height, width and angle. A 90mm B 90mm C 200mm D 120mm E 25mm 5230-3251-XXX	 For Size 2 only. Provides extra support on shoulders or lumbar area. 420mm (W) 5230-2721-XXX
	

Sternal Prompt
(For Prone). 5230-2910-000
Cover (required) 5230-8801-XXX


Waistcoat Harness

Size 2	Size 3	Size 4	Mounting Kit	
Provides additional trunk and shoulder stability.				
A 5cm B 32cm C 14cm D 15.5cm 5230-6712-000	A 6.5cm B 37cm C 17.5cm D 18cm 5230-6713-000	A 8.5cm B 39.5cm C 20.5cm D 18.5cm 5230-6714-000	Mounting bracket for Harness (Required for Size 1 Multistander with Abduction only). 5230-PART-056	

Knee Supports


Size 1	Size 2	
Easy to attach and adjust for height, angle & depth.		
A 100mm / B 45mm / C 75mm 5356-1751-XXX	A 120mm / B 60mm / C 95mm 5356-1752-XXX	


Tray Options

Elbow Blocks	
Suits all sizes. 10cm (W) x 14cm (H) x 8.5cm (D) 5356-7351-000	

Adjustable Tray - Size 1	Adjustable Tray - Size 2	
A 650mm B 450mm C 370mm D (With Insert) 330mm / (Without Insert) 290mm 5230-7111-000	A 700mm B 530mm C 460mm D (With Insert) 410mm / (Without Insert) 360mm 5230-7112-000	

Foot Positioning

Sandals			
Size 1 170mm (L) x 80mm (W) 5356-6081-000	Size 2 220mm (L) x 95mm (W) 5356-6082-000	Size 3 275mm (L) x 105mm (W) 5356-6083-000	

Sandal Raiser Blocks			
Size 1 Kit of 4x15mm blocks 5356-6091-000	Size 2 Kit of 4x15mm blocks 5356-6092-000	Size 3 Kit of 4x15mm blocks 5356-6093-000	

Ankle Snugs		
Size 1 120mm - 225mm Circumference 120mm (H) 5356-6131-000	Size 2 165mm - 295mm Circumference 120mm (H) 5356-6132-000	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Multistander_20241031_AU