

# Jenx Standz



## ABDUCTION STANDER

AUSTRALIAN  
SCRIPTING AND ORDERING GUIDE

Distributed by  
**medifab**

# Client Details


Client/User Name:

Date:

NDIS/Funding:

### User Measurements

A Hip width	<input type="text"/>
B Axilla height	<input type="text"/>
C Chest width	<input type="text"/>
D Overall user height	<input type="text"/>
E Heel to knee	<input type="text"/>
F Hip height	<input type="text"/>
User weight	<input type="text"/>



# Medifab Notes

---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---

# Product Details

## Size 1 - Manual Tilt

Average Age	1 - 9 Years
Weight Capacity	45 kg
Hip Width	17.5cm - 30.5cm
Hip Height	45cm - 70cm
Chest Width	17.5cm - 30.5cm
Chest Height	58cm - 100cm
Max Height	138cm
Tilt Range (90° = Upright)	5° to 90°
Frame Width	64.5cm
Frame Length	89cm
Height of surface when horizontal	81cm
Abduction Range	30° Either Leg

## Key

**Blue** Suitable for **Prone and Supine**

**Green** Suitable for **Supine only**

**Yellow** Suitable for **Prone only**

Items outlined in blue are required and an option must be selected

## Base unit



### Standz Basic Package



Tool free adjustment standing frame includes base with castors, hydraulic assisted tilt, pelvic and trunk support pads with strap, abducting legs, foot positioning plates, tray mount.

5240-0005-044



# Knee support options

Standard Knee Support		Multi-Axis Knee Support	
Dual wall knee supports with accessories available for extra function (see below).		Width, height, depth and angle adjustable, with soft arched foam pads behind the knee and support straps to enable knee flexion.	
5240-1750-044		5240-1770-044	
A: 11.5cm, B: 12.5cm, C: 9.5cm & Min: 6.5cm		A: 13.5cm & B: 11.5cm	

Knee Blocks		Knee Block Mounting Bracket	
Flexible PU knee cups, complete with padded knee straps. Adjustable for height, width, angle and depth.			
<b>Knee Blocks Size 1</b> Internal Support Size: 7.5cm W x 10cm H		Mounting bracket <b>required</b> when using knee block option for knee support.	
5356-1751-011			
<b>Knee Blocks Size 2</b> Internal Support Size: 9.5cm W x 12cm H			
5356-1752-011		5240-1720-000	

# Accessories for Standard Knee Support

Standard Strap (Required for Supine use)	De-rotational Strap	Arched Foam Pads (Recommended for Supine use)	Qty
Holds knee in position and prevents knee flexion. Can also be used in Prone.	Adjustable strap for reducing internal or external rotation.	Ensures a comfortable standing position when accommodating knee flexion. Supplied individually. 13.5cm W x 13.5 cm H x 5.5cm D	
			
5240-1766-044	5240-1786-044	5353-1775-000	

# Foot support options

Size 1 sandals		Size 2 sandals	
Ensure solid foot placement		Ensure solid foot placement	
5356-6081-000		5356-6082-000	
8cm x 17cm long		9.5cm x 22cm long	

# Head support option

## Flat Pad

For children with reasonable head control. Can also be used as chest pad in prone.  
5240-4110-044

19cm W x 10cm H



## Multigrip - 1 Link

Single link wings for some degree of lateral support.  
33cm W x 5cm D x 7.5cm H  
5356-4601-000

Wipeable Cover (Optional)  
5356-4671-000

## Multigrip - 2 Link

2 Links in each wing for extra support and positioning.  
47.5cm W x 5cm D x 12cm H  
5356-4602-000

Wipeable Cover (Optional)  
5356-4672-000

## Multigrip - 3 Link

3 adjustable links in each wing for ultimate lateral headrest control and comfort.  
62cm W x 5cm D x 12cm H  
5356-4603-000

Wipeable Cover (Optional)  
5356-4673-000



# Shoulder support option

## Flat Shoulder Support

Shoulder support for supine use only (too wide for prone)

5240-2750-044

31.5cm W x 9.5cm H



## Multigrip body support (for prone and supine)

Adjustable contour shoulder support

5240-2721-044

41cm max width



# Pelvic support options

## Pelvic Support Board

Lower back support

5240-1901-044

19.5cm W x 15cm H



## Back Pad (Requires Pelvic support board)






Central back support

5240-2402-044

19.5cm W x 15cm H



# Accessories

Anti Tamper Conversion Kit		Multi-position tray (mounts to bracket in base package)	
Allows for all wing knob adjustments to be converted to a tool fixing – ideal for environments where easy adjustment can be problematic.		Includes extension wings with slot.	
5240-0170-000		5240-7112-044	
		A: 58.5cm & B: 39.5cm C: (Prone - 25cm, Supine - 45cm)	
Tray holder (useful for Supine)		Tray Elbow Blocks	
When in Supine, the tray needs to be removed to remove the child. Place the tray on the convenient tray holder to avoid tripping.		Aids in supporting the arms and promoting a functional arm and hand position. <b>Only usable in conjunction with Tray Extension Wings.</b>	
5240-7365-000		5356-7351-000	
		10cm W x 14cm H	

# Medifab Notes

---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---



---

# medifab

Distributed by Medifab Ltd

26 Pardoe Street,  
East Devonport,  
Tasmania, 7310

1300 543 343

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Jenx Standz\_20231025\_AU