

Juni Chair



INDOOR SEATING

SCRIPTING AND ORDERING GUIDE

Distributed by
medifab

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Size & Colour Options

	Size 1	Size 2	Size 3	Size 4
Black	3209-1001-099	3209-1002-099	3209-1003-099	3209-1004-099
Blue	3209-1001-060	3209-1002-060	3209-1003-060	3209-1004-060
Red	3209-1001-022	3209-1002-022	3209-1003-022	3209-1004-022
Lime	3209-1001-054	3209-1002-054	3209-1003-054	3209-1004-054

Seat Unit Measurements				
Seat height (cm)	16-22	20-28	24-34	30-40
Seat width (cm)	24	28	32.5	36.5
Seat depth (cm)	18-24	22-30	26-36	32-42
Back support height (cm)	26	29	32.5	40
Load capacity (kg)	30	40	50	60
Product weight (kg)	4	6	8	10

Base package includes:

- Seat Unit
- Hip Belt



Postural Support Accessories


Side Cushions - For Size 1 & 2 Juni	Side Cushions - For Size 3 & 4 Juni
-------------------------------------	-------------------------------------

Reduce seat width using side support (supplied as a pair).	
To reduce width by 5cm 3209-0001-025	To reduce width by 5cm 3209-0003-025
To reduce width by 10cm 3209-0001-050	To reduce width by 10cm 3209-0003-050
To reduce width by 15cm 3209-0011-050	To reduce width by 15cm 3209-0013-050



Lateral Supports - For Size 1 & 2 Juni	Lateral Supports - For Size 3 & 4 Juni
----------------------------------------	----------------------------------------

3209-3001-000	3209-3003-000
---------------	---------------



Medial Thigh Support

Padded support for increased medial thigh support and seat contouring.	
For size 1, 2 & 3 3209-1501-099	
For size 4 3209-1504-099	



Anterior Chest Support

Size 1	Size 2	Size 3	Size 4
Red 3201-6401-022	Red 3201-6402-022	Red 3201-6403-022	Red 3201-6404-022
Lime 3201-6401-054	Lime 3201-6402-054	Lime 3201-6403-054	Lime 3201-6404-054
Blue 3201-6401-060	Blue 3201-6402-060	Blue 3201-6403-060	Blue 3201-6404-060
Black 3209-6351-099	Black 3209-6352-099	Black 3209-6353-099	Black 3209-6354-099




Footsupport

Adjustable depth tray table, providing a work surface and an upper limb support.	
Size 1 3209-5001-000	
Size 2 3209-5002-000	
Size 3 3209-5003-000	
Size 4 3209-5004-000	



Base Options

Standard Skis		Standard Skis with Footboard	
Makes moving Juni chair in a straight line easier (under tables).		Standard skis with built in foot support, great for users who will push and could cause their chair to tip.	
Size 1 3209-5021-000		Size 1 3209-5061-000	
Size 2 3209-5022-000		Size 2 3209-5062-000	
Size 3 3209-5023-000		Size 3 3209-5063-000	
Size 4 3209-5024-000		Size 4 3209-5064-000	

Mobile Skis		
Makes the Juni chair easier to move, ideal for pushing under tables!		
Size 1 3209-5051-000		
Size 2 3209-5052-000		
Size 3 3209-5053-000		
Size 4 3209-5054-000		

Sliding Footboard (for use with mobile skis)		Height Adjustable Push Handle	
Add foot support to mobile skis to support users legs.		Great for use with standard or mobile skis.	
Size 1 3209-5071-000		Size 1 3209-8611-000	
Size 2 3209-5072-000		Size 2 3209-8612-000	
Size 3 3209-5073-000		Size 3 3209-8613-000	
Size 4 3209-5074-000		Size 4 3209-8614-000	

Growth Accessories

Backrest Extension (Increases back height as below)

Size 1: 11cm	Size 2: 11cm	Size 3: 13cm	Size 4: 16cm
Red 3209-2101-022	Red 3209-2102-022	Red 3209-2103-022	Red 3209-2104-022
Lime 3209-2101-054	Lime 3209-2102-054	Lime 3209-2103-054	Lime 3209-2104-054
Blue 3209-2101-060	Blue 3209-2102-060	Blue 3209-2103-060	Blue 3209-2104-060
Black 3209-2101-099	Black 3209-2102-099	Black 3209-2103-099	Black 3209-2104-099



Outriggers

Increases seat to floor height by 5cm.

Size 1 3209-2811-000	Size 2 3209-2812-000	Size 3 3209-2813-000
--------------------------------	--------------------------------	--------------------------------



Table Options

Adjustable & Foldable Table Quick Release Tray

Adjustable height work surface for children which can be easily folded for storage, a great companion to the Junii chair. This table comes in a beautiful 'natural' finish that complements the colour of the legs. (We don't currently offer the green finish shown in the image).

Size 1
Table height 31-45cm
3209-7191-000



Size 2
Table height 40-60cm
3209-7192-000



Adjustable depth tray table, providing a work surface and an upper limb support.

Size 1
3201-7020-000

Size 2
3201-7021-000

Size 3
3201-7022-000

Size 4
3201-7023-000



medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Juni_20240823