

Kidsflex





KIDSFLEX CAR SEAT

Distributed by
medifab

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE

Basic Packages

Size 2	Size 3
 <p><i>Laterals not included, select on next page</i></p>	 <p><i>Laterals not included, select on next page</i></p>
<p>Includes:</p> <ul style="list-style-type: none"> • Blue • 5 Point Harness • Recline Function • Tilt angle adjustment (0 to 20°) • Isofix attachment • Standard head support • Slim left turn swivel base • Seat belt guide on the left 	<p>Includes:</p> <ul style="list-style-type: none"> • Grey • 5 Point Harness • Recline Function • Tilt angle adjustment (0 to 20°) • Isofix attachment • Standard head support • Slim left turn swivel base • Seat belt guide on the left
4170-2022-060	4170-2023-090

Specifications


Size	Total Height	Seat to Head	Back Height	Seat Width	Seat Depth	Lateral Sizes	Distance of Lateral Supports	Weight	Max User Weight
2	80 cm	73 cm	41-57 cm	35 cm	34-39 cm	26-36 cm	19-32 cm	8.9 kg	36 kg
3	90 cm	82 cm	47 - 65cm	40 cm	40-46 cm	38-48 cm	24-37 cm	12.3 kg	75 kg

Tray

Upper Limb Support Tray

Size 2
4170-5532-000

Size 3
4170-5533-000



Foot Support

Foot support for use with either size Kidsflex

4170-5242-000



Lateral Supports

Small Lateral Supports (Recommended)

Fixed
4170-3301-000

Swing-Away
4170-3401-000



Large Lateral Supports

Fixed
4170-3302-000

Swing-Away
4170-3402-000



medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Kidsflex_20240624_AU