

# Klaxon Range



**KLAXON**<sup>®</sup>

Distributed by  
**medifab**

SCRIPTING AND ORDERING GUIDE

October 2023

# Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
<b>User Measurements</b>	
Chair Model	
Chair Width	
Client Weight	
Tube Diameter	

# Klaxon Klick Race

## Basic Package Includes:

- Battery (11Ah) with Battery Charger
- LCD Colour Display
- USB-Port Phone Charger
- Foldable Kick Stand
- Splashguard
- 14" Smooth Tyre
- 1500w Motor
- 11kg without Battery
- Cruise Control
- Electronic Brake
- Max Speed - 20km




## Standard Handlebar

1595-2022-KLICK-R-HUB



## Tetra Handlebar

1595-2021-KLICK-R-TETRA-HUB



# Klaxon Klick Power

## Basic Package Includes:

- Standard Handlebar
- Battery (11Ah) with Battery Charger
- LCD Colour Display with Automatic Night Mode
- USB-Port Phone Charger
- Foldable Kick Stand
- Splashguard
- 14" Smooth Tyre
- 1000w Motor
- 8.5kg without Battery



## Standard Handlebar

1595-2021-KLICK-P-HUB



## Tetra Handlebar

1595-2021-KLICK-P-TETRA-HUB





# Connectors

## Complete Connector (Front Connector + Rear Connector + Front Hook)

	Small 1595-LSTD-0002	Standard 1595-LSTD-0000	Extra Long 1595-LSTD-0001
--	-------------------------	----------------------------	------------------------------

## Complete Connector Tetra (Complete Connector + Quick-Release Lever Kit + Velcro Grip Kit)

	Standard 1595-LTET-0000	Extra Long 1595-LTET-0001
--	----------------------------	------------------------------

# Clamping Rings

## Round-shaped, Convergent Tube Frames (1 piece) (Tapered Seat Frame)

	Diameter 25,4mm 1595-AGCP125	Diameter 32,4mm 1595-AGCP132
--	---------------------------------	---------------------------------

## Round-Shaped, Parallel Tube Frames (1 piece)

	Diameter 25,4mm 1595-AGP25	Diameter 32,4mm 1595-AGP32	Diameter 34,925mm 1595-AGP35
--	-------------------------------	-------------------------------	---------------------------------

## Clamping Ring for Oval-shaped, Parallel Tube Frames (1 piece)

	For Progeo 25mm x 35mm (WxH) 1595-AGPOV	For Sunrise 27mm x 33mm (WxH) 1595-AGPOV SX
--	---	---

## Clamping Ring for Progeo Carbon Noir 2.0 and Duke

	(1 piece) 1595-AGPOVPC
--	---------------------------

# Reducers

## Reducer

	Converts 25, 4 into 18/18, 4	1595-RID25-18-0000
	Converts 25, 4 into 20/20, 4	1595-RID25-20-0000
	Converts 25,4 into 23	1595-RID25-23-0000
	Converts 25,4 into 22/22,4	1595-RID25-22-0000
	Converts 32,4 into 28/28,4	<b>\$150</b> 1595-RID32-28-0000
	Converts 32,4 into 30/30,4	1595-RID32-30-0000

# Crossbeam Options

## Complete Crossbeam Range (Includes Crossbeam + Fixing Clamp + Aluminium Axle)

	230-270mm (Small Crossbeam)	1595-TRAV2327
	271-290mm (Medium Crossbeam)	1595-TRAV2729
	291-360mm (Medium Crossbeam)	1595-TRAV2936
	361-480mm (Medium Crossbeam)	1595-TRAV3648
	480-540mm (Medium Crossbeam)	1595-TRAV4854


# Quick Release Locking Pin

## with Ring Handle



(1 piece)  
1595-DF-34-0001




# Handlebar Options

Tetra Handle - Vertical (1 Piece)		Quick Release Adjustable Handlebar Stem	
1595-T-HANDLE-02		(ES, PS) 1595-CC-23-22-24-0001	
Appendices for Handlebar		Handlebar Extension Standard	
1595-T-HORN-0000		1595-CC-25-0005	

# Tyre Options



Spare Tyres		Fat Wheels (Pair) with Complete Mounting Kit	
20 x 210 (HS)	1595-400730		12, 7mm 1595-FATWHEELS-01-12,7
FAT (FS)	1595-CM-3-2000		
Tread Tyre 14"	1595-CC-3-1400	<b>Spare Tyre Tubes</b>	
Tread Tyre 14" Smooth	1595-CC-3-1401		20" (FS) 1595-CM-87-2000


# Battery Options


Spare Battery		
	48V 11Ah (HS) - Range 50km Standard Spare Battery for Handcycle	1595-CH-1-11AH
	Fly Race 48V 2,9Ah (P, R, M) - Range 12km	1595-CM-1-2,9AH
	Race 48V 11,6Ah (P, R, M) Standard Spare Battery for Race, Power, Monster	1595-CM-1-11AH
	Race Lite 48V 5Ah (P, R, M) - Range 12km	1595-CM-1-5AH
Spare Battery Charger		
	48V 11,6Ah (RS, RT, FS, PSL, PTL, Fly version too)	1595-CC-88-0001
	48V 11,6Ah (HS)	1595-K500050

# Accessories


Shopping Bag Frame		Shopping Bag (Requires Shopping Bag Frame)			
1595-T-SHOPBAG		Left/Right Side (1 piece) 1595-SHOPBAG-0000		Front Side 1595-SHOPBAG-0001	

Basket		Basket Linking Kit	
	10L 1595-BASKET-0000		(ES, PSL, RS, FS, P, R, M) 1595-KIT-TEL-0000
	16L 1595-BASKET-0001		(PS, PT) 1595-KIT-TEL-0001

Flybag	
14" 1595-FLYBAG00	

Single Rear View Mirror	
(Excluding HS) 1595-SPE-0000	

LED Light (includes batteries and mount)		
	Front (White) (RS, RT, FS, PSL, PTL)	1595-CC-83-KIT
	Back (Red)	1595-LED-RED-0000

Additional Traction Weight Kit	
	(2 pieces) 1595-CC-93-0000-2



# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Klaxon\_20231012