


Flux Walker

The Flux is a lightweight but sturdy posterior walking aid which encourages upright posture and promotes natural walking movements.

Dependable support for children & adults



Height-adjustable grip bars grow with the user, making it easy to adjust to different body heights

Reverse-roll locks (included in package) can be deactivated so the walker can roll forward and backward

Opens out in walking direction and optically and physically free from interfering structural parts


Folding mechanism with safety lock eases transport and storage

Friction brakes and bumper wheels included in package


Which size?




Added extras



Swivel Lock
5640-8440-000



Carrier
Size 1-2 5640-8061-000
Size 3-4 5640-8063-000



Basket (can only be used if Carrier is installed)
5640-8051-000

ACCESSORIES

Fold-Up Seat

1-3 5640-0201-000
4 5640-0204-000



Multi Adjustable Hip/Bottom Pad

5640-1900-000



Dynamic Hip Control

For Orthoses
5640-1951
Soft Suspension (Yellow)
5640-1952
Hard Suspension (Red)
5640-1953

Add Hip Belt (Required)
Hip 53-83 cm 5640-6201
Hip 65-115 cm 5640-6202



Sling Seat

Hip 57-77 cm
5640-1921-000
Hip 73-103 cm
5640-1922-000

Add Sling Seat Rail
Size 1 5640-1931-000
Size 2 5640-1932-000
Size 3-4 5640-1933-000



Forearm Supports with Hand Grips

2 Axial 5640-5501-000
3 Axial 5640-5502-000



Driving Brake & Wheel Lock

Size 1-3 5640-8451-000
Size 4 5640-8454-000



Anti-Tipper

Size 1-3 5640-8421-000
Size 4 5640-8424-000

Anti-Tipper & Castor

Size 1-3 5640-8425-000
Size 4 5640-8428-000



Universal Handlebars

Size 1 5641-5751-000
Size 2 5642-5752-000
Size 3 5643-5753-000
Size 4 5644-5754-000

Not compatible with forearm supports



Size	Code	Colour	Standard Grip Height	Universal Grip Height	Frame Width	Frame Length	Frame Weight	Max User Weight
1	5641-0001-000	Luminous Red	48-61 cm	43-68 cm	64 cm	65 cm	6.5 kg	60 kg
2	5642-0001-000	Sky Blue	57-81 cm	51-88 cm	69 cm	81 cm	6.8 kg	75 kg
3	5643-0001-000	Gentian Blue	69-87 cm	62-94 cm	74 cm	88 cm	7.1 kg	75 kg
4	5644-0001-000	Grey Aluminium	82-102 cm	75-102 cm	78 cm	91 cm	9.1 kg	100 kg