


Flux Walker

The Flux walker is a walking aid that promotes natural walking. Opening out in walking direction and optically and physically free from interfering structural parts, the Flux Walker encourages upright posture and natural walking movements.

Lightweight posterior walker: Dependable support for children & adults



Height-adjustable grip bars grow with the user, making it easy to adjust to different body heights

Opens out in walking direction and optically and physically free from interfering structural parts

Folding mechanism with **safety lock** eases transport and storage


Reverse-roll locks (included in package) can be deactivated so the walker can roll forward and backward

Friction brakes and bumper wheels included in package


Which size?




Added extras



Swivel Lock
5640-8440-000



Carrier
Size 1-2 5640-8061-000
Size 3-4 5640-8063-000



Basket (can only be used if Carrier is installed)
5640-8051-000

ACCESSORIES

Fold-Up Seat

1-2 5640-0201-000

3-4 5640-0204-000



Multi Adjustable Hip/Bottom Pad

5640-1900-000



Dynamic Hip Control

Kit for Orthoses

5640-1951-000

Soft Suspension

5640-1952-000

Hard Suspension

5640-1953-000

Add Hip Belt (Required)

Small - Hip Size: 53-83 cm

5640-6201-000

Large - Hip Size: 65-115 cm

5640-6202-000



Sling Seat

Hip 57-77 cm

5640-1921-000

Hip 73-103 cm

5640-1922-000

Add Sling Seat Rail

Size 1 5640-1931-000

Size 2 5640-1932-000

Size 3-4 5640-1933-000



Forearm Supports with Hand Grips

2 Axial 5640-5501-000

3 Axial 5640-5502-000



Driving Brake & Wheel Lock

Size 1-3 5640-8451-000

Size 4 5640-8454-000



Anti-Tipper

Size 1-3 5640-8421-000

Size 4 5640-8424-000

Anti-Tipper & Castor

Size 1-3 5640-8425-000

Size 4 5640-8428-000



Universal Handlebars


Size 1 5641-5751-000

Size 2 5642-5752-000

Size 3 5643-5753-000

Size 4 5644-5754-000

Not compatible with forearm supports



Size	Code	Colour	Standard Grip Height	Universal Grip Height	Frame Width	Frame Length	Frame Weight	Max User Weight
1	5641-0001-000	Luminous Red	48-61 cm	43-68 cm	64 cm	65 cm	6.5 kg	60 kg
2	5642-0001-000	Sky Blue	57-81 cm	51-88 cm	69 cm	81 cm	6.8 kg	75 kg
3	5643-0001-000	Gentian Blue	69-87 cm	62-94 cm	74 cm	88 cm	7.1 kg	75 kg
4	5644-0001-000	Grey Aluminium	82-102 cm	75-102 cm	78 cm	91 cm	9.1 kg	100 kg