


Jenx Atom Chair

Modern, child-centred styling, combined with a comprehensive range of support options make Atom the ideal seating system to develop correct sitting posture and function.

A modular seating system for users 6 months to 5 years



Y base provides a **range of seated height options** to suit both floor level and table top activities

Colour choice of red or pink

Open sides ensures adjustments are easily made with the user remaining seated

Tool free adjustments, cleverly designed wing knobs offer a comfortable grip to enable effective clamping with minimal effort

Unrivalled seat depth growth to accommodate early intervention requirements

Memory foam ischial well designed to nestle the pelvis. Provides comfort, pressure distribution and stability to promote correct pelvic alignment

RED PACKAGE:

- Padded seat with memory foam
- Shaped back pad
- Arm supports
- Nursery tray
- Adjustable flip away hip support pads
- Push Handle
- Upper back pad
- Flip up angle adjustable footplate

3431-1001-022



Y-Base:

- Versatile and flexible, with lockable swivel castors
- High-low range of 25.8-75cm seat height from floor
- Tilt from 30° recline - 10° prone
- Base footprint of 65 x 78cm

3615-8740-000



PINK PACKAGE:

- Padded seat with memory foam
- Shaped back pad
- Adjustable flip away hip support pads
- Push handle
- Upper back pad
- Arm supports
- Tray table support
- Shoulder protractors
- High rise flexi thoracic pads
- Thoracic pad covers
- 2 Link Multigrip Head Support
- Flip up angle adjustable footplate

3431-0010-028



Functional benefits:

- **Head Positioning:** Through the innovative Multigrip Head Support, there are a number of solutions that adapt to varying head positioning needs.
- **Trunk Positioning:** All trunk positioning accessories are optional and are intended to offer just the right level of support.
- **Leg Positioning:** With a number of different control options windsweeping, abduction and adduction can be easily influenced or controlled using our standard leg positioning accessories.
- **Foot Positioning:** The child's feet can be precisely positioned using a number of accessories.

Colour	Code	Approx Age Range	Max User Weight	Tilt	Back Pad Height	Chest Width	Hip Width	Seat Depth	Lower Leg Length	Foot Support Angle
Red	3431-1001-022	6mths - 5yrs	25kg	15° - 15°	27-38cm	20-28cm	18-27cm	9-34cm	16.5-32.5cm	5° - 35°
Pink	3431-0010-028									

Jenx Atom Accessories

Multigrip Head Support
1 Link 3431-4601-022
2 Link 3431-4602-022




Multigrip Head Support
1 Link Wipeable Cover 5356-4671-000
2 Link Wipeable Cover 5356-4672-000



Shoulder Protractors
3431-2750-022



Sacral Block
3431-2501-000




Back Pad Infill
3431-2370-000



High-Rise Flexible Thoracic Support Pads
3431-3251-022




High-Rise Short Thoracic Support Pads
3430-3151-000



Nursery Tray Infill
3431-7182-000



Half Padding
3431-7442-000
Full Padding
3431-7452-000



Activity Frame
3615-7475-000
Grab Rail
3615-7200-000



3 Point Lap Strap
3431-6251-000
Spex 4 Point Lap Strap
1406-6420-027



Waistcoat Harness
3615-6701-000



Wonderseat Centrepoint Shoulder Harness:
Size 0 1409-5660-017



Wonderseat 4-Point Hip Belt:
1406-5420-027



Abduction Block
Size 0 3615-1650-000
Size 1 3615-1651-000



Long Leg Sitting Board
3431-0281-022



Leg Control Pads
3431-1751-000



Sandals
5356-6081-000



Ankle Snugs
5356-6131-000




Sandal Raising Blocks
5356-6091-000



Shoulder Protractors
3431-2750-022



Elbow Blocks
3730-7350-000



Foot Support Pad
3431-5701-000

