

Sharky Bathroom Chair

Sharky has been designed for multifunctional applications and ease of handling.

Multifunctional commode for for shower, bath and toilet use.



Removable rear seat section to assist with hygiene care

Removable seat with a range of adapters offers versatile support in different areas of the bathroom

Detachable foot support system

Easy to clean all surfaces

Slim, compact foldable frame



Sharky Size 1




Sharky Size 2

Seat Angle	Back Support Angle	Knee Angle	Foot Support Angle
-10° to 45°	72-120°	78-138°	80-100°

Size	Back Support Height			Seat Width	Seat Depth	Toilet Clearance	Armrest Height	Wheel Diameter	Weight Limit
	35cm	45cm	55cm						
1	6081-0010-100	6081-0010-000	N/A	24, 27, 30cm	27-36cm	52cm	14.5 - 20.5cm	75mm	60kg
2	N/A	6082-0015-100	6082-0015-000	29, 32, 35cm	31-40cm	52cm	14.5 - 20.5cm		

BACK SUPPORT RECLINE

Simple to use system for easy and secure handling.



45° TILT

The special construction of the tilt mechanism allows for 45° tilt of seat to ensure comfort and safety especially when showering.



HEIGHT ADJUSTMENT

Smooth movement with comfortable handles make adjustments easy and controlled.



COMPACT FOLDING

Fold down to compact height of just 28cm, even with seat attached, makes transport and storage easy.



REMOVABLE PADDING

Easy care removable padding for cleaning the Sharky.



Accessories


Removable Footrest System
 Size 1: 6085-1810-000
 Size 2: 6085-1820-000



Push Handles
 6085-8600-000



Lateral Trunk Supports
 Size 1: 6085-3310-000
 Size 2: 6085-3320-000



Headrest
 6085-4750-000



Neoprene Shoulder Anterior Chest Support
 Size 1: 6085-6610-000
 Size 2: 6085-6620-000



Hip Belt with Neoprene Pad Sizes
 Size 1: 6085-1910-000
 Size 2: 2880-1920-000



Removable & Height Adjustable Arm Supports
 6085-5500-000



Abduction Block
 Size 1: 6085-1500-000
 Size 2: 6085-1520-000



Depth Adjustable, removable tray table
 Size 1: 6085-7100-000
 Size 2: 6085-7200-000



Note: Arm Supports required

Neoprene Calf Strap
 Size 1: 6085-6160-000
 Size 2: 6085-6260-000



Neoprene Ankle Huggers
 Size 1: 6085-6101-000
 Size 2: 6085-6102-000



Removable Commode Pan with retainer
 6085-8210-000



Removable Splash Guard
 6085-7675-000



Seat Infill pad
 6085-0401-000




Bath subframe with seat angle adjustment
 Size 1: 6085-0500-000
 Size 2: 6085-0520-000



Note: Minimum width of bath 36 cm

Toilet Adapter
 Size 1: 6085-0600-000
 Size 2: 6085-0620-000



Seat width reducing hip supports - Size 1 Sharky
 24cm: 6085-1110-000
 27cm: 6085-1120-000
 30cm: 6085-1130-000



Seat width reducing hip supports - Size 2 Sharky
 29cm: 6085-1210-000
 32cm: 6085-1220-000
 35cm: 6085-1230-000



Upgrade Kit
 Upgrade from Sharky Size 1 to Size 2
 6085-8081-086

