



Back Support Hardware






















	QUAD MOUNT	QUAD MOUNT FIXED	QUAD MOUNT HEAVY-DUTY	DUAL MOUNT	DIRECT MOUNT	DROP MOUNT	XLELLA BARIATRIC
POINTS OF FIXATION	4	4	4	2	4	2	4
DESCRIPTION	Quick-release, compact and highly-adjustable hardware. Default for most Spex backs.	Fixed and compact. Best option for non-removable backs.	Compact and strong. Supports users who apply a lot of force and/or weight to the back support.	Ideal option when extra adjustment in angle and depth is needed.	The most minimal Spex option available. It is non-removable and only allows for width adjustment.	Easy drop-in hardware for active and lightweight back supports.	Bariatric hardware for heavier users.
HARDWARE WEIGHT	0.9kg	0.6kg	1.1kg	1.1kg	0.6kg	0.6kg - 0.8kg	1.1kg (Fixed) 1.3kg (Quick-Release)
FITS TO BACKS	Spex Backs (Excluding Adapta/Zygo).	Spex backs (Excluding Adapta/Zygo).	Spex backs (Excluding Adapta/Zygo). Use for backs > 20" wide or 23" high.	All Spex backs (Can't be used with wheelchairs with cross-brace).	Spex backs (Excluding Adapta).	Spex and Zygo backs up to 18" high.	XLella backs.
FITS CANE SIZES	1.9 - 2.6cm	1.9 - 2.6cm	1.9 - 2.6cm	1.9 - 2.6cm	1.9 - 2.6cm	1.9 - 3.2cm	1.9 - 2.6cm
BACK ANGLE ADJUSTMENT	+/- 10° (Add Depth Extension accessory to upper brackets for additional +10°)	Non-Adjustable	+/- 10° (Add Depth Extension accessory to upper brackets for additional +10°)	+/- 25°	Non-Adjustable	+/- 30°	+/- 10°
SEAT DEPTH ADJUSTMENT	4cm (Achieve 10cm with Depth Extension accessory)	4cm (Achieve 10cm with Depth Extension accessory)	4cm (Achieve 10cm with Depth Extension accessory)	4cm (Additional 4cm is available when Dual Mount is mounted in 'Forward' configuration)	Non-Adjustable	2.5cm (Small) 5cm (Paediatric, Medium) 7.5cm (Large)	3cm
OVERALL WIDTH ADJUSTMENT	5cm	25 mm	5cm	5cm	5cm	8cm	5cm
WEIGHT CAPACITY	10 – 15" Backs: 90kg 16" + Backs: 136kg	10 – 15" Backs: 90kg 16" + Backs: 136kg	10 – 15" Backs: 136kg 16" + Backs: 180kg	10 - 15" Backs: 90kg 16" + Backs: 136kg	10 - 15" Backs: 90kg 16" + Backs: 136kg	136kg	300kg
CODE	1259-2400-100	1259-2470-100	1259-2405-100 (Backs up to 23" wide) 1259-2406-100 (Backs 23" wide +)	1259-2200-100 (Standard) 1259-2205-100 (For Adapta) 1259-2210-100 (Forward) 1259-2215-100 (Forward-Extended)	1259-2460-100	1259-1400-00P (Paediatric) 1259-1400-00S (Small) 1259-1400-00M (Medium) 1259-1400-00L (Large)	1259-2500-100 (Fixed) 1259-2505-100 (Quick-Release)

Spex is as diligent as possible in compiling this information. Product availability may differ between regions, please contact your local Product Specialist to confirm availability. Information is correct at time of printing/release: 02/07/24 For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Spex. Product specifications are subject to change without notice. © Spex 2024