

Pelvi-Loc Hip Belt

Experience ultimate pelvis retraction with the 'Pelvi-Loc' hip belt! These superior high tensile strength 3D retraction belts are suitable for all wheelchair users.

ENHANCED PELVIC STABILITY FOR FUNCTION AND COMFORT



Optimal **3-DIMENSIONAL** contoured pelvic positioning pads provide safety with high comfort.

Contoured shape **REDUCES PELVIC PRESSURE**

Positive **TOOTHED RATCHET** tensioning system with easy to use ratchet levers

Always maintain **CORRECT PELVIC TENSION** if the user wears a different thickness of clothing (when the weather changes!)

ADD EXTRA PRESSURE RELIEF

Experience ultimate comfort with Medifab Artisan upgrades to the Pelvi-Loc hip belt. Combining enhanced stability with the pressure-reducing comfort of Pudgee foam.



12mm Pudgee foam provides increased shear resistance, optimal to allow increased movement while **REDUCING PRESSURE INJURIES.**



Open cell structured foam for **INCREASED BREATHABILITY.**

BUILD CODE FOR PELVI-LOC HIP BELT

1 4 3 0 - 6 3 X X - X 2 X

Shape/Mounting	Size	Comfort Option	Buckle Option
Basic/Solid 5	X-Small 0	Standard 0	Standard Compact 4
Basic/Webbing 6	Small 1	Integrated Pudgee 2	Easy 5
Active/Solid 7	Medium 2	Removeable pudgee padding also available	
Active/Webbing 8	Large 3		

SPECIFICATIONS

	X-Small	Small	Medium	Large
W/C WIDTH	6" - 9"	8" - 12"	12" - 16"	14" - 18"

SHAPE OPTIONS



BASIC STYLE:
EXPOSED BUCKLES,
SLIM PAD

Low profile pads which fit securely over the pelvis. Suited for all wheelchairs including manual and power chairs.



ACTIVE STYLE:
COVERED BUCKLES,
LARGE PAD

Provides more support and surface area over thighs and has front padded cover over buckles to avoid unintentional opening of the belt and prevent harm to users with uncontrolled arm movements e.g. Athetoid. Suited for sport players.

BUCKLE OPTIONS



COMPACT BUCKLE
Low profile use.
Most commonly chosen buckle.



EASY BUCKLE
Easily opened for those with limited hand function, e.g. tetraplegics, stroke victims.

MOUNT OPTIONS



WEBBING MOUNT
Loop through standard fixings or attach via mounting kit to the seat canes. Most common mount option; however not recommended for clients with strong extensor thrust.



SOLID MOUNT
Rigid mount which bolts to the chair. For heavy thrusting clients or active sports players to prevent movement.

ACCESSORIES



REMOVEABLE PADDING
Basic: 1430-636X-334
Active: 1430-637X-334
X-Small: **1** Medium: **3**
Small: **2** Large: **4**