

Pixi Stroller

An attractive, lightweight stroller with rigid seating shell construction.

Outstanding stability & quick, easy folding



Easy kerb navigation
Handle bar which extends over rear axle

Back support recline
allows the child to sleep or relax

Easy rear wheel brake system

Robust, sturdy design
for excellent reliability

Manoeuvrable!
Great stroller to push, with ball bearing castors to easily navigate shopping malls and more

Standard Package Includes:

- Foldable Frame
- Grey Upholstery (wipeclean seat surfaces)
- Flip-up Foot Plate
- 5 Point Harness

Configure for behaviours of concern

- Anti-escape hip belt 1405-6211-026
- Pen-release hip belt 1405-6211-025
- Anti-escape tee-pad anterior chest support configuration 2560-6556-026




Parents feel great with this smart design.


Quickly folds down to a compact size to place in a vehicle boot. Easy 'lift-to-fold' catches.

ACCESSORIES


Basket
1 2560-8011-017
2 2560-8012-017




Chest Strap complete with lateral pads
1 2560-6801-090
2 2560-6802-090




Standard Lateral Head Support
Sml 2010-4210-XXX
Lrg 2010-4211-XXX



Extended Lateral Head Support
Sml 2010-4310-XXX
Lrg 2010-4311-XXX



Occipital Head Support
2010-4110-XXX



Head Support Fabric Options:

- Merino/Spacer -099
- Wipe Down -209


Tray Table
1 2560-7120-013
2 2560-7130-013




Medial Knee Support
2560-1501-090




Rain Shield Cover
2560-7601-010



UV Shade Cover
2560-7701-010



Sun Canopy
1 2560-7543-090
2 2560-7544-090



Grab Rail
1 2560-7253-299
2 2560-7254-299



Tie-Downs
2560-8501-021



Size	Code	Approx Age Of Child	Approx Weight Of Child	Seat Depth	Seat Width	Shoulder Height	Foot Support Height	Approx Stroller Weight
1	2563-0006-191	3 - 7yrs	35kg	25 or 30cm	28cm	36 - 45cm	18 - 30cm	12.5kg
2	2566-0006-191	6 - 14yrs	50kg	35 or 40cm	32cm	45 - 55cm	18 - 35cm	15kg