

Pixi Stroller



**RIGID &
LIGHTWEIGHT**

Distributed by
medifab

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE


Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Measurements	Size 1	Size 2
Pixi Order Code	2563-0006-191	2566-0006-191
Max. User Weight	35kg	50kg
Max. User Height	100cm (120cm with extension)	125cm (150cm with extension)
A Seat Depth	25cm (30cm with extension)	35cm (40cm with extension)
B Seat Width	28cm	33cm
C Shoulder Height	36-45cm	45-55cm
D Footplate Height	13-35cm	13-39cm
E Overall Width	58.5cm	63.5cm
F Push Handle Height	57-110cm	62-120cm
Stroller Weight (approx)	12.5kg	15kg
Folded Dimensions	L: 73cm W: 59cm H: 40cm	L: 85cm W: 66cm H: 51cm
Back Angle	90-125°	90-123°
Basket Weight Limited	5kg	5kg

Stroller Package includes:

- Foldable Frame
- Grey Upholstery
- Flip-up Foot Plate
- 5 Point Harness









Lateral Head Support

Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
Merino/Spacer Fabric 2010-4110-099		Small (17cm Internal width)		Small (17cm Internal width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
		Wipe Down Cover 2010-4210-209		Wipe Down Cover 2010-4310-209	
Wipe Down Cover 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)	
		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	

Growth Adjustments

Seat Depth Extension Bracket		Foot Plate Spacer	
To increase seat depth by 5cm		Size 1 To decrease footplate height 2560-5301-299	
Size 1 2560-0281-199			
Size 2 2560-0282-199			


Harnesses & Hip Belts

Anti-Escape Tee Pad Harness		Chest Strap with Lateral Pads	
Fitted with anti-escape buckles to prevent user from escaping stroller. 2560-6556-026		Size 1 2560-6801-090	
		Size 2 2560-6802-090	
2-Point Padded Centre Pull Hip Belt			
Centre-Release Buckle 1405-6211-028	Anti-Escape Buckle 1405-6211-026	Pen Popper Buckle 1405-6211-025	
			

Abduction Pommel

To keep user's legs in the correct position

2560-1501-090



Tray Table

Activity and Feeding Tray

Size 1
2560-7120-013

Size 2
2560-7130-013



Grab Rail

Attaches to the Pixi frame

Size 1
2560-7253-299

Size 2
2560-7254-299



Weather Protection

Rain Shield Cover

For protection from the rain
2560-7601-010



UV Shade Cover


For protection from the sun
2560-7701-010



Sunroof

Size 1
2560-7543-090

Size 2
2560-7544-090



Storage Options

Mesh Basket with Plastic Bottom (weight limit 5kg)

Size 1
2560-8011-017

Size 2
2560-8012-017



Tie-Downs

ISO 7176-19

2560-8501-021



medifab

Distributed by Medifab Ltd

26 Pardoe Street
East Devonport
Tasmania 7310

+61 3 6427 0105

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Pixi_20231018_AU