

# Poz'In'Form



## CLINICAL LYING & POSITIONING SOLUTIONS

medifab


SCRIPTING AND ORDERING GUIDE

# Client Details

Client/User Name:	
Clients DOB:	
Date:	

Notes:	
--------	--

- Poz'In'Form® elements and modules are entirely removable and washable. They are delivered with a Pharmatex or Lenzing FR® cover.
- The Pharmatex bi-elastic cover is made of polyurethane and is recommended for use in healthcare facilities. Pharmatex is soft to the touch and very comfortable next to the skin, while still easy to clean.
- The Lenzing FR® cover is made of viscose, coated with a waterproof polyurethane film. Lenzing FR® is a major innovation in positioning of the long-term bedridden, preventing maceration. Hygroscopic properties of viscose allow moisture absorption which is particularly important for people with thermoregulation issues. Moreover elasticity and flexibility of Lenzing FR® are particularly suitable for problems regularly observed in classical positioning.
- Lenzing spare covers can be purchased to cover Pharmatex pillows for a more absorbent feel.



# Poz'In'Form Kit

<b>Pharmatex</b> 7009-0001-269	<b>Lenzing FR</b> 7009-0001-279
-----------------------------------	------------------------------------

- Transport bag
- Universal cushion (550mm x 400mm)
- Half-moon cushion
- Semi-lateral decubitus cushion
- Cylindrical cushion
- Ring cushion
- Heel cushion with strap
- Occiput head cushion



# Half-Moon

<b>Pharmatex</b> 7009-1851-269	<b>Lenzing FR</b> 7009-1851-279
	
<b>Width</b> 300mm <b>Length</b> 2000mm	
<b>Spare Cover</b> 7009-1851-269-C	<b>Spare Cover</b> 7009-1851-279-C

# Delta

<b>Pharmatex</b> 7009-0551-269	<b>Lenzing FR</b> 7009-0551-279
	
<b>Length</b> 560mm <b>Depth</b> 320mm <b>Height</b> 280mm	
<b>Spare Cover</b> 7009-0551-269-C	<b>Spare Cover</b> 7009-0551-279-C

# Semi-Lateral Decubitus

<b>Pharmatex</b> 7009-3001-269	<b>Lenzing FR</b> 7009-3001-279
	
<b>Width</b> 350mm <b>Length</b> 1850mm	
<b>Spare Cover</b> 7009-3001-269-C	<b>Spare Cover</b> 7009-3001-279-C


# Cylindrical

<b>Pharmatex</b> 7009-1501-269	<b>Lenzing FR</b> 7009-1501-279
	
<b>Length</b> 750mm <b>Diameter</b> 180mm	
<b>Spare Cover</b> 7009-1501-269-C	<b>Spare Cover</b> 7009-1501-279-C

# Wedge Abduction

<b>Pharmatex</b> 7009-1551-269	<b>Lenzing FR</b> 7009-1551-279
	
	<b>Length</b> 350mm <b>Height</b> 150mm <b>Small Width</b> 100mm <b>Large Width</b> 300mm
<b>Spare Cover</b> 7009-1551-269-C	<b>Spare Cover</b> 7009-1551-279-C

# Occiput

<b>Pharmatex</b> 7009-4751-269	<b>Lenzing FR</b> 7009-4751-279
	
<b>Length</b> 500mm <b>Width</b> 400mm <b>Height</b> 90mm	
<b>Spare Cover</b> 7009-4751-269-C	<b>Spare Cover</b> 7009-4751-279-C

# Ring & Half Ring

Pharmatex 7009-4081-269		Lenzing FR 7009-4081-279		Pharmatex 7009-4051-269		Lenzing FR 7009-4051-279	
<b>Ring</b>				<b>Half Ring</b>			
<b>Width</b>	200mm			<b>Width</b>	220mm		
<b>Length</b>	1350mm			<b>Length</b>	1050mm		
<b>Spare Cover</b>	7009-4081-269-C	<b>Spare Cover</b>	7009-4081-279-C	<b>Spare Cover</b>	7009-4051-269-C	<b>Spare Cover</b>	7009-4051-279-C

# Universal

Pharmatex		Lenzing FR	
<b>Code</b>	7009-1251-269	<b>Code</b>	7009-1251-279
<b>Length</b>	350mm	<b>Length</b>	350mm
<b>Depth</b>	250mm	<b>Depth</b>	250mm
<b>Spare cover:</b> 7009-1251-269-C		<b>Spare cover:</b> 7009-1251-279-C	
<b>Code</b>	7009-1253-269	<b>Code</b>	7009-1253-279
<b>Length</b>	400mm	<b>Length</b>	400mm
<b>Depth</b>	150mm	<b>Depth</b>	150mm
<b>Spare cover:</b> 7009-1253-269-C		<b>Spare cover:</b> 7009-1253-279-C	
<b>Code</b>	7009-1252-269	<b>Code</b>	7009-1252-279
<b>Length</b>	550mm	<b>Length</b>	550mm
<b>Depth</b>	400mm	<b>Depth</b>	400mm
<b>Spare cover:</b> 7009-1252-269-C		<b>Spare cover:</b> 7009-1252-279-C	

# Hand Wedge

Pharmatex 7009-5411-269		Lenzing FR 7009-5411-279	
<b>Length</b>	390mm		
<b>Height</b>	90mm		
<b>Width</b>	220mm		
<b>Spare Cover</b>	7009-5411-269-C	<b>Spare Cover</b>	7009-5411-279-C

# Heel Wedge

Pharmatex 7009-5401-269		Lenzing FR 7009-5401-279	
<b>Length</b>	480mm		
<b>Height</b>	120mm		
<b>Width</b>	300mm		
<b>Spare Cover</b>	7009-5401-269-C	<b>Spare Cover</b>	7009-5401-279-C

# Heel Cushion

Pharmatex 7009-5201-269		Lenzing FR 7009-5201-279	
<b>Length</b>	850mm		
<b>Diameter</b>	250mm		
<b>Spare Cover</b>	7009-5201-269-C	<b>Spare Cover</b>	7009-5201-279-C

# medifab

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY: