

Simple Stuff Works



SLEEPING & POSITIONING SYSTEM

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

Client details

Client/User name:			
DOB:		Date:	
Therapist:			

Notes

Further Information

CHOOSE A **PACKAGE** TO USE ON **FOAM MATTRESS**


Packages

Infant Package	Child Package	Adult Package
<p>For infants from 0 to 5 years of age</p> <ul style="list-style-type: none">• 2x cot Bed Mesh• 2x cot Toppers• 2x cot cotton sheets• 2x small Lateral Supports• 2x small Wedges• Small Horseshoe pillow• Simple Stabiliser Size 1• Travel bag <p>7150-7001-C00</p>	<p>For children from 5 to 12 years of age</p> <ul style="list-style-type: none">• 2x single Bed Mesh• 2x single Toppers• 2x single cotton sheets• 2x medium Lateral Supports• 2x small Wedges• 2x small Lateral Supports• Simple Stabiliser Size 2• Small Terry Towelling Pillow• Travel bag <p>7150-7002-S00</p>	<p>For adults from 12+ years of age</p> <ul style="list-style-type: none">• 2x single Bed Mesh• 2x single Toppers• 2x single cotton sheets• 2x medium Lateral Supports• 2x medium Wedges• 2x large Lateral Supports• Large Terry Towelling Pillow• Travel bag <p>7150-7003-S00</p>





OR REFER TO **PAGE 3** TO **CONFIGURE** YOUR OWN

Sheet Sets

Cot Bed	Single Bed	
<p>Set Includes:</p> <ul style="list-style-type: none"> Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top) 2 x cot temperate toppers 2 x cot fitted sheets <p>7150-7021-C00</p>	<p>Set Includes:</p> <ul style="list-style-type: none"> Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top) 2 x single temperate toppers 2 x single fitted sheets <p>7150-7022-S00</p>	

Sheets


Stabilising Mesh (Set of 2)		
<p>Cot Bed</p> <p>Base layer: 1500 mm (L) x 700 mm (W) Top layer: 1800 mm (L) x 1300 mm (W)</p> <p>7150-7011-C00</p>	<p>Single Bed</p> <p>Base layer: 1800 mm (L) x 1000 mm (W) Top layer: 1800 mm (L) x 1500 mm (W)</p> <p>7150-7012-S00</p>	

Temperate Topper (Set of 2)		
<p>Cot Bed</p> <p>1000 mm (L) x 590 mm (W)</p> <p>7150-7051-C00</p>	<p>Single Bed</p> <p>1000 mm (L) x 1000 mm (W)</p> <p>7150-7052-S00</p>	

Fitted Sheets			
<p>Cot Bed</p> <p>1400 mm (L) x 700 mm (W)</p> <p>7150-7121-C00</p>	<p>Single Bed/King Single</p> <p>1900 mm (L) x 900 mm (W)</p> <p>7150-7122-S00</p>	<p>Double/Queen Bed</p> <p>2030 mm (L) x 1530 mm (W)</p> <p>7150-7123-S00</p>	

Flat Sheet
<p>Single</p> <p>1620 mm (L) x 2600 mm (W)</p> <p>7150-7126-S00</p>

Wedges


Wedges (pair)		
With padded covers.		
<p>Small</p> <p>A: 220 mm B: 110 mm C: 110 mm</p> <p>7150-7351-000</p>	<p>Medium</p> <p>A: 360 mm B: 80 mm C: 90 mm</p> <p>7150-7352-000</p>	

Lateral Supports

Padded Lateral Supports




Also available in an option for air mattresses, designed to grip directly to the mattress surface.

Small (Foam) 7150-7201-000	Medium (Foam) 7150-7202-000	Large (Foam) 7150-7203-000
Small (Air) 7150-7211-000	Medium (Air) 7150-7212-000	Large (Air) 7150-7213-000
A: 180 mm B: 120 mm C: 180 mm	A: 220 mm B: 140 mm C: 220 mm	A: 280 mm B: 220 mm C: 280 mm



Simple Stabiliser

Supports hips and thighs in neutral supine alignment

Size 1	Size 2	Size 3	Size 4
H: 110mm D: 110mm W: 420mm	H: 190mm D: 350mm W: 460mm	H: 250mm D: 400mm W: 650mm	H: 280mm D: 400mm W: 650mm
To fit approx leg circumference of 250mm. Suggested age range 2 - 4yrs.	To fit approx leg circumference of 350mm. Suggested age range 2 - 10yrs.	To fit approx leg circumference of 550mm. Suggested age range 5 - 13yrs.	To fit approx leg circumference of 600mm. Suggested age range 13 - adult.
With breathable inner & terry cotton outer 7150-7341-000	With breathable inner & terry cotton outer 7150-7342-000	With breathable inner & terry cotton outer 7150-7343-000	With breathable inner & terry cotton outer 7150-7344-000
With wipeable incontinence inner & terry cotton outer 7150-7345-000	With wipeable incontinence inner & terry cotton outer 7150-7346-000	With wipeable incontinence inner & terry cotton outer 7150-7347-000	With wipeable incontinence inner & terry cotton outer 7150-7348-000
			
Spare terry cotton cover 7150-7221-000	Spare terry cotton cover 7150-7222-000	Spare terry cotton cover 7150-7223-000	Spare terry cotton cover 7150-7224-000
Spare wipeable inner cover 7150-7225-000	Spare wipeable inner cover 7150-7226-000	Spare wipeable inner cover 7150-7227-000	Spare wipeable inner cover 7150-7228-000

Side Lying Leg Support

Assists with pelvic and thigh positioning in side-lying

Size 1	Size 2	Size 3
A: 270mm B: 380mm C: 280mm D: 80mm E: 130mm F: 80mm	A: 400mm B: 560mm C: 340mm D: 100mm E: 200mm F: 90mm	A: 400mm B: 610mm C: 400mm D: 110mm E: 150mm F: 70mm
With breathable inner & terry cotton outer (standard) 7150-7971-000	With breathable inner & terry cotton outer (standard) 7150-7972-000	With breathable inner & terry cotton outer (standard) 7150-7973-000
With wipeable incontinence inner & terry cotton outer 7150-7981-000	With wipeable incontinence inner & terry cotton outer 7150-7982-000	With wipeable incontinence inner & terry cotton outer 7150-7983-000
Spare terry cotton cover 7150-7991-000	Spare terry cotton cover 7150-7992-000	Spare terry cotton cover 7150-7993-000

Foot Support

Provides positioning & skin protection for the feet (pair)

Size 1	Size 2	Size 3	
A: 210mm (max) 7150-5351-000	A: 250mm (max) 7150-5352-000	A: 300mm (max) 7150-5353-000	

Pillows

Roll Pillow

Small A: 150mm B: 400mm	Large A: 180mm B: 500mm	
With breathable inner & terry cotton outer 7150-7851-000	With breathable inner & terry cotton outer 7150-7852-000	
With wipeable incontinence inner & terry cotton outer 7150-7855-000	With wipeable incontinence inner & terry cotton outer 7150-7856-000	
Spare terry cotton cover 7150-7857-000	Spare terry cotton cover 7150-7858-000	

Neck Support Pillow

(Suitable for air & foam mattresses). 300mm x 350mm x 100mm 7150-4050-000	
Spare terry cotton cover 7150-4055-000	

Horseshoe Pillow

Small A: 190mm B: 460mm	Large A: 200mm B: 500mm	
With breathable inner & terry cotton outer 7150-7711-000	With breathable inner & terry cotton outer 7150-7712-000	
With wipeable incontinence inner & terry cotton outer 7150-7721-000	With wipeable incontinence inner & terry cotton outer 7150-7722-000	
Spare terry cotton cover 7150-7731-000	Spare terry cotton cover 7150-7732-000	

Terry Towling Pillow

Small L: 480mm W: 350mm 7150-7501-000	Large L: 680mm W: 450mm 7150-7502-000	
Spare terry cotton cover 7150-7511-000	Spare terry cotton cover 7150-7512-000	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: SimpleStuffWorks_20241031