

Simple Stuff Works



SLEEPING & POSITIONING SYSTEM

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

Client details

Client/User name:			
DOB:		Date:	
Therapist:			

Notes

Further Information

CHOOSE A **PACKAGE** TO USE ON **FOAM MATTRESS**


Packages

Infant Package	Child Package	Adult Package
<p>For infants from 0 to 5 years of age</p> <ul style="list-style-type: none"> • 2x cot Bed Mesh • 2x cot Toppers • 2x cot cotton sheets • 2x small Lateral Supports • 2x small Wedges • Small Horseshoe pillow • Simple Stabiliser Size 1 • Travel bag <p>7150-7001-C00</p>	<p>For children from 5 to 12 years of age</p> <ul style="list-style-type: none"> • 2x single Bed Mesh • 2x single Toppers • 2x single cotton sheets • 2x medium Lateral Supports • 2x small Wedges • 2x small Lateral Supports • Simple Stabiliser Size 2 • Small Terry Towelling Pillow • Travel bag <p>7150-7002-S00</p>	<p>For adults from 12+ years of age</p> <ul style="list-style-type: none"> • 2x single Bed Mesh • 2x single Toppers • 2x single cotton sheets • 2x medium Lateral Supports • 2x medium Wedges • 2x large Lateral Supports • Large Terry Towelling Pillow • Travel bag <p>7150-7003-S00</p>




OR REFER TO **PAGE 3** TO **CONFIGURE** YOUR OWN


Sheet Sets

Cot Bed	Single Bed	
<p>Set Includes:</p> <ul style="list-style-type: none"> Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top) 2 x cot temperate toppers 2 x cot fitted sheets <p>7150-7021-C00</p>	<p>Set Includes:</p> <ul style="list-style-type: none"> Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top) 2 x single temperate toppers 2 x single fitted sheets <p>7150-7022-S00</p>	

Sheets


Stabilising Mesh (Set of 2)		
<p>Cot Bed</p> <p>Base layer: 1500 mm (L) x 700 mm (W) Top layer: 1800 mm (L) x 1300 mm (W)</p> <p>7150-7011-C00</p>	<p>Single Bed</p> <p>Base layer: 1800 mm (L) x 1000 mm (W) Top layer: 1800 mm (L) x 1500 mm (W)</p> <p>7150-7012-S00</p>	

Temperate Topper (Set of 2)		
<p>Cot Bed</p> <p>1000 mm (L) x 590 mm (W)</p> <p>7150-7051-C00</p>	<p>Single Bed</p> <p>1000 mm (L) x 1000 mm (W)</p> <p>7150-7052-S00</p>	

Fitted Sheets			
<p>Cot Bed</p> <p>1400 mm (L) x 700 mm (W)</p> <p>7150-7121-C00</p>	<p>Single Bed/King Single</p> <p>1900 mm (L) x 900 mm (W)</p> <p>7150-7122-S00</p>	<p>Double/Queen Bed</p> <p>2030 mm (L) x 1530 mm (W)</p> <p>7150-7123-S00</p>	

Flat Sheet
<p>Single</p> <p>1620 mm (L) x 2600 mm (W)</p> <p>7150-7126-S00</p>

Wedges


Wedges (pair)		
With padded covers.		
<p>Small</p> <p>A: 220 mm B: 110 mm C: 110 mm</p> <p>7150-7351-000</p>	<p>Medium</p> <p>A: 360 mm B: 80 mm C: 90 mm</p> <p>7150-7352-000</p>	

Lateral Supports

Padded Lateral Supports




Also available in an option for air mattresses, designed to grip directly to the mattress surface.

Small (Foam) 7150-7201-000	Medium (Foam) 7150-7202-000	Large (Foam) 7150-7203-000
Small (Air) 7150-7211-000	Medium (Air) 7150-7212-000	Large (Air) 7150-7213-000
A: 180 mm B: 120 mm C: 180 mm	A: 220 mm B: 140 mm C: 220 mm	A: 280 mm B: 220 mm C: 280 mm



Simple Stabiliser

Supports hips and thighs in neutral supine alignment

Size 1	Size 2	Size 3	Size 4
H: 110mm D: 110mm W: 420mm	H: 190mm D: 350mm W: 460mm	H: 250mm D: 400mm W: 650mm	H: 280mm D: 400mm W: 650mm
To fit approx leg circumference of 250mm. Suggested age range 2 - 4yrs.	To fit approx leg circumference of 350mm. Suggested age range 2 - 10yrs.	To fit approx leg circumference of 550mm. Suggested age range 5 - 13yrs.	To fit approx leg circumference of 600mm. Suggested age range 13 - adult.
With breathable inner & terry cotton outer 7150-7341-000	With breathable inner & terry cotton outer 7150-7342-000	With breathable inner & terry cotton outer 7150-7343-000	With breathable inner & terry cotton outer 7150-7344-000
With wipeable incontinence inner & terry cotton outer 7150-7345-000	With wipeable incontinence inner & terry cotton outer 7150-7346-000	With wipeable incontinence inner & terry cotton outer 7150-7347-000	With wipeable incontinence inner & terry cotton outer 7150-7348-000
			
Spare terry cotton cover 7150-7221-000	Spare terry cotton cover 7150-7222-000	Spare terry cotton cover 7150-7223-000	Spare terry cotton cover 7150-7224-000
Spare wipeable inner cover 7150-7225-000	Spare wipeable inner cover 7150-7226-000	Spare wipeable inner cover 7150-7227-000	Spare wipeable inner cover 7150-7228-000

Side Lying Leg Support

Assists with pelvic and thigh positioning in side-lying

Size 1	Size 2	Size 3
A: 270mm B: 380mm C: 280mm D: 80mm E: 130mm F: 80mm	A: 400mm B: 560mm C: 340mm D: 100mm E: 200mm F: 90mm	A: 400mm B: 610mm C: 400mm D: 110mm E: 150mm F: 70mm
With breathable inner & terry cotton outer (standard) 7150-7971-000	With breathable inner & terry cotton outer (standard) 7150-7972-000	With breathable inner & terry cotton outer (standard) 7150-7973-000
With wipeable incontinence inner & terry cotton outer 7150-7981-000	With wipeable incontinence inner & terry cotton outer 7150-7982-000	With wipeable incontinence inner & terry cotton outer 7150-7983-000
Spare terry cotton cover 7150-7991-000	Spare terry cotton cover 7150-7992-000	Spare terry cotton cover 7150-7993-000

Foot Support

Provides positioning & skin protection for the feet (pair)

<p>Size 1 A: 210mm (max) B: 200mm 7150-5351-000</p>	<p>Size 2 A: 250mm (max) B: 220mm 7150-5352-000</p>	<p>Size 3 A: 300mm (max) B: 300mm 7150-5353-000</p>	
--	--	--	--

Pillows

Roll Pillow

Small A: 150mm B: 400mm	Large A: 180mm B: 500mm	
With breathable inner & terry cotton outer 7150-7851-000	With breathable inner & terry cotton outer 7150-7852-000	
With wipeable incontinence inner & terry cotton outer 7150-7855-000	With wipeable incontinence inner & terry cotton outer 7150-7856-000	
Spare terry cotton cover 7150-7857-000	Spare terry cotton cover 7150-7858-000	

Neck Support Pillow

(Suitable for air & foam mattresses). 300mm x 350mm x 100mm 7150-4050-000	
Spare terry cotton cover 7150-4055-000	

Horseshoe Pillow

Small A: 190mm B: 460mm	Large A: 200mm B: 500mm	
With breathable inner & terry cotton outer 7150-7711-000	With breathable inner & terry cotton outer 7150-7712-000	
With wipeable incontinence inner & terry cotton outer 7150-7721-000	With wipeable incontinence inner & terry cotton outer 7150-7722-000	
Spare terry cotton cover 7150-7731-000	Spare terry cotton cover 7150-7732-000	

Terry Towling Pillow

Small L: 480mm W: 350mm 7150-7501-000	Large L: 680mm W: 450mm 7150-7502-000	
Spare terry cotton cover 7150-7511-000	Spare terry cotton cover 7150-7512-000	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: SimpleStuffWorks_20241121