

Juni Chair

Whether in the classroom, nursery or at home, better posture while seated creates a more positive experience for your child. It can help develop fine motor skills, increase confidence and the ability to engage and learn.

Confidence, engagement, learning

Sturdy frame helps to weight and prevent chair movements

Gentle contouring provides support and gives stability for sitting, promoting engagement and focus on the task at hand

Juni's low seat to floor height allows children to sit with their feet on the floor

Tailor Juni to your child's needs:

With exceptional adjustment of Juni, along with a great range of optional accessories, your child's Juni can be tailored to them, supporting them to achieve their goals and more!

Colour Options:

Base Package includes seat unit and hip belt. Available in four colours.



Black (-099)



Blue (-060)



Red (-022)




Lime (-054)

Size	Code	Seat Height	Seat Width	Seat Depth	Back Support Height	Weight Capacity	Product Weight
1	3209-1001-0XX	16-22cm	24cm	18-24cm	26cm	30kg	4kg
2	3209-1002-0XX	20-28cm	28cm	22-30cm	29cm	40kg	6kg
3	3209-1003-0XX	24-34cm	32.5cm	26-36cm	32.5cm	50kg	8kg
4	3209-1004-0XX	30-40cm	36.5cm	32-42cm	40cm	60kg	10kg

ACCESSORIES

Side Cushions for Size 1 and 2 Juni
 25mm 3209-0001-025
 50mm 3209-0001-050
 75mm 3209-0011-050



Side Cushions for Size 3 and 4 Juni
 25mm 3209-0003-025
 50mm 3209-0003-050
 75mm 3209-0013-050

Standard Skis
 1 3209-5021-000
 2 3209-5022-000
 3 3209-5023-000
 4 3209-5024-000



Mobile Skis
 1 3209-5051-000
 2 3209-5052-000
 3 3209-5053-000
 4 3209-5054-000



Sliding Footboard
 1 3209-5071-000
 2 3209-5072-000
 3 3209-5073-000
 4 3209-5074-000



Adjustable and Foldable Table
 1 3209-7191-000 (Height 31-45cm)
 2 3209-7192-000 (Height 40-60cm)



Full range of accessories available in the Junii Chair Scripting and Ordering Guide.