

Zoomi+ High Chair

Zoomi+ is the ideal feeding chair for children and adults who require a little extra postural support. With its mobile skis, quick release tray and padded seat, the Zoomi+ comes with some great benefits. It can be easily maneuvered from room-to-room, enable mealtime inclusion and provide the user with greater comfort.

The one-size-fits all feeding chair for all ages!



Wipe-clean surfaces as standard make cleaning easy.

Sturdy frame design provides stability and prevents chair movements.

Gentle contouring provides support and gives stability.

Quick-release tray supports upper limbs while providing a space for activities

Mobile skis and an integrated handle mean the Zoomi+ can be moved around home with ease.

ACCESSORIES

Centrepoint Harness
0 1409-5660-017 (Pad height 22cm)
1 1409-5661-017 (Pad height 28cm)



2-Point Hip Belt X-Small
1405-6220-027 (Pad width 30cm)



Transfer Step
Black 3212-5901-099
Blue 3212-5901-060
Red 3212-5901-022
Lime 3212-5901-054



Extended Slotted Footboard
Black 3212-5021-099
Blue 3212-5021-060
Red 3212-5021-022
Lime 3212-5021-054



Extended Contoured Padded Seat
Black 3212-8991-099
Blue 3212-8991-060
Red 3212-8991-022
Lime 3212-8991-054



Zoomi+ Package Contains:

- Seat Unit
- High Back Support
- Contoured Padded Seat
- Arm Supports
- Foot Support
- Mobile Skis
- 5-Point Harness
- Quick-Release Tray

Colour Options:



Black (-099)

Blue (-060)

Red (-022)

Lime (-054)

ANTI-MICROBIAL PROPERTIES

The Zoomi+'s wooden frame has anti-microbial properties that ensure surfaces are kept hygienic, which helps to reduce any risk of cross-infection and is especially useful in multi-user environments such as nurseries and schools.



Code	Frame Height	Seat Width	Frame Depth	Seat Depth	Back Pad Height	Foot Support Height	Weight Capacity	Product Weight
3212-1001-0XX	81 cm	43 cm	55 cm	15-26 cm	27 cm	26-60 cm	85 kg	16 kg