

Symmetrisleep



SLEEPING & POSITIONING SYSTEM

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:		Notes
Date:		
DOB:		

- Contents:**
- Symmetrisleep Positioning System
 - Symmetrisleep Thermal Regulating Options
 - Lowzone Posture Cushions

Velco Receptor Sheet

Cot	Single	Single King	Double
Stretchy velcro-like fitted sheet for securing of Positioning Brackets and other supports.			
Length 1000mm Width 750mm	Length 1900mm Width 1000mm	Length 2050mm Width 1060mm	Length 2050mm Width 1500mm
7010-8005-CTO	7010-8005-SNO	7010-8005-KSO	7010-8005-D00



Positioning Brackets (Pair)

Small Qty:	Medium Qty:	Large Qty:	Extra Large Qty:
Brackets that attach to Velco Receptor Sheet.			
Colour: Red	Colour: Yellow	Colour: Blue	Colour: Green
Width 100mm Height 80mm Length 130mm	Width 160mm Height 130mm Length 130mm	Width 150mm Height 250mm Length 250mm	Width 150mm Height 280mm Length 250mm
7011-7200-001	7011-7220-001	7011-7225-001	7011-7230-001



Coolover Bracket Pads - Padding to be installed on Positioning Brackets			
Small Qty:	Medium Qty:	Large Qty:	Extra Large Qty:
7012-7165-001	7012-7175-001	7012-7185-001	7012-7195-001



Monoblok Knee Supports

Monoblok Supine Knee Support

Small 7013-7552-000		Medium 7013-7553-000		Large 7013-7554-000	
Length	100mm	Length	200mm	Length	250mm
Width	360mm	Width	400mm	Width	600mm



Monoblok Rolling Knee Support

Small 7013-7562-000		Medium 7013-7563-000		Large 7013-7564-000	
Length	180mm	Length	200mm	Length	250mm
Width	100mm	Width	210mm	Width	270mm



Monoblok Knee Support Accessories

Coolover Legwrap Liners - 10mm (suits either of Monoblok supports)

Small (80mm Diameter) 7013-7411-CVR	Medium (80mm Diameter) 7013-7412-CVR	Large (160mm Diameter) 7013-7413-CVR
---	--	--



Monoblok Ramp Wedges

Medium 7013-7573-000	Large 7013-7575-000
Length 400mm Width 600mm	Length 500mm Width 700mm



Monoblok Safety Straps to secure to mattress

Single Strap 7013-7585-000	Double Strap 7013-7586-000
-------------------------------	-------------------------------

Airmantle

Mattress topper for draping over brackets and supports for air circulation/thermal comfort

Single	
7014-8010-S00	
Length	1900mm
Width	1000mm
Depth	20mm



Coolover Sheet

Temperature regulation sheets.

Cot		Single	
Length	1300mm	Length	1900mm
Width	600mm	Width	1060mm
7015-7120-C00		7015-7120-S00	

Single King		Double/Queen		Flat Sheet	
Length	2100mm	Length	2150mm	Length	3200mm
Width	1050mm	Width	1500mm	Width	2100mm
7015-7120-SK0		7015-7120-D00		7015-7120-SFS	



Symmetrisleep Stickers











Provides the carers with customised user instruction sheets

- Sheets allow you to choose a lying position and then select the parts you need, to assist with planning and recording sleep system set-up.
 - Add the selected model to the chart, and accessories can be placed underneath as well as on top, to mirror the sleep system set-up of the client.
 - Space to add name and instructions.
- 7010-0061-000



Lowzone Posture Cushions

Wipe-down PU fabric pillows with chipped memory foam filling accessible via zip opening

Square 7016-7500-001		Med Oblong 7016-7580-001		Lrg Oblong 7016-7610-001		Sml Horseshoe 7016-7720-001	
Length	300mm	Length	450mm	Length	750mm	Length	370mm
Width	300mm	Width	300mm	Width	450mm	Width	250mm
							
Coolover Cover 7016-7500-CVR		Coolover Cover 7016-7580-CVR		Coolover Cover 7016-7610-CVR		Coolover Cover 7016-7720-CVR	
Med Horseshoe 7016-7710-001		Lrg Horseshoe 7016-7700-001		Horseshoe Extended 7016-7702-001		Sml Roll 7016-7800-001	
Length	650mm	Length	800mm	Length	650mm	Length	300mm
Width	305mm	Width	350mm	Width	300mm	Width	80mm
							
Coolover Cover 7016-7710-CVR		Coolover Cover 7016-7700-CVR		Coolover Cover 7016-7702-CVR		Coolover Cover 7016-7800-CVR	
Med Roll 7016-7850-001		Med Roll Extended 7016-7860-001					
Length	750mm	Length	370mm				
Width	450mm	Width	250mm				
							
Coolover Cover 7016-7850-CVR		Coolover Cover 7016-7860-CVR					

Notes

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Symmetrisleep 20241030_NZ