



Helio KIDS
ORDER FORM - AUSTRALIA(UNPRICED)

1- Purchasing information

P.O. number:* _____ Date: _____

Company Name:* _____

Purchasing contact:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone Email

Marked for: _____

3- Shipping information

Ship to:* _____

Name: _____

Address: _____

City: _____ State: _____

Zip Code: _____

Phone: _____

2- Clinician information

Clinician:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone Email

MA083 MA084 MA085 MA086 MA087 MA088 MA089

English manual (STD) German manual Italian manual Dutch manual Spanish manual French manual Portuguese manual

MODEL		AUD (taxes not included)					
WC06	<input type="checkbox"/>	Ultralight Folding Carbon Fiber Helio KIDS Wheelchair Including fixed seat depth OR depth adjustable back (79,4 kg - 175 lb weight capacity)				BASE FRAME	
WC10	<input type="checkbox"/>	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber				Upcharge	
FREE Growth Kit within the first 5 years of the Helio KIDS wheelchair ownership The growth kit includes carbon crossbraces, crossbrace links, seat sling, back upholstery, front riggings and footplates.							
SEAT WIDTH ¹							
SW	<input type="checkbox"/>	30.48 (12") <input type="checkbox"/>	33.02 (13") <input type="checkbox"/>	35.56 (14") <input type="checkbox"/>	38.10 (15") <input type="checkbox"/>	40.64 (16") <input type="checkbox"/>	NC

[1. Width Calculator](#)

↓ **Caution: complete only one of the following two sections** ↓

SEAT DEPTH (Available with WC-19)			
SD12	<input type="checkbox"/>	30.48cm/12" Seat sling - frame of 35.56 cm/14"	NC
SD13	<input type="checkbox"/>	33.02.02cm/13" Seat sling - frame of 35.56 cm/14"	NC
SD14	<input type="checkbox"/>	35.56cm/14" Seat sling - frame of 35.56 cm/14"	NC
SD15	<input type="checkbox"/>	38.10cm/15" Seat sling - frame of 38.10cm/15"	NC
SD16	<input type="checkbox"/>	40.64cm/16" Seat sling - frame of 40.64cm/16"	NC
FIXED BACK CANE (MUST select ONE push handle options)			
BC11-1215	<input type="checkbox"/>	Straight back cane 30.48-38.10cm (12-15")	NC
BC11-1518	<input type="checkbox"/>	Straight back cane 38.10-45.72cm (15-18")	NC
BC11-1821	<input type="checkbox"/>	Straight back cane 45.72-53.34cm (18-21") STD	NC
BC12-1315	<input type="checkbox"/>	Back cane 8° bend 34.29-38.10cm (13.5-15")	NC
BC12-1518	<input type="checkbox"/>	Back cane 8° bend 38.10-45.72cm (15-18")	NC
BC12-1821	<input type="checkbox"/>	Back cane 8° bend 45.72-53.34cm (18-21") STD	NC
Specify height: ⁵		measured from seat sling to top of back upholstery	
ANGLE ADJUSTABLE BACK CANE (MUST select ONE push handle options)			
BC13-1518	<input type="checkbox"/>	Straight back cane 38.10-45.72cm (15-18")	NC
BC13-1821	<input type="checkbox"/>	Straight back cane 45.72-53.34cm (18-21") STD	NC
BC14-1518	<input type="checkbox"/>	Back cane 8° bend 38.10-45.72cm (15-18")	NC
BC14-1821	<input type="checkbox"/>	Back cane 8° bend 45.72-53.34cm (18-21") STD	NC
Select angle for the back cane:		85° <input type="checkbox"/> 90° <input type="checkbox"/> 95° STD <input type="checkbox"/> 100° <input type="checkbox"/> 105° <input type="checkbox"/> 110° <input type="checkbox"/>	
Specify height : ⁵			

⁵- If not indicated, default height setting is at the lowest position

SEAT DEPTH- Depth and angle adjustable (Not available with WC-19) ³			
SD14	<input type="checkbox"/>	30.48 - 38.10 cm (12-15") (35.56/14" seat sling/frame)	NC
SD15	<input type="checkbox"/>	33.02 - 40.64 cm (13-16") (38.10/15" seat sling/frame)	NC
SD16	<input type="checkbox"/>	35.56 - 43.18 cm (14-17") (40.64/16" seat sling/frame)	NC
SD17	<input type="checkbox"/>	38.10 - 45.72 cm (15-18") (43.18/17" seat sling/frame)	NC
SD18	<input type="checkbox"/>	40.64 - 48.26 cm (16-19") (45.72/18" seat sling/frame)	NC
Specify depth ⁴ :			
³ - Maximum weight of 175 lb with a depth adjustable backrest			
⁴ - If not indicated, the depth will be equal to the seat sling			

DEPTH AND ANGLE ADJUSTABLE BACK CANE (MUST select ONE push handle options)			
BC15-912	<input type="checkbox"/>	Straight back cane 22.86-30.48cm (9-12")	NC
BC15-1215	<input type="checkbox"/>	Straight back cane 30.48-38.10cm (12-15")	NC
BC15-1518	<input type="checkbox"/>	Straight back cane 38.10-45.72cm (15-18")	NC
BC15-1821	<input type="checkbox"/>	Straight back cane 45.72-53.34cm (18-21") STD	NC
BC16-1215	<input type="checkbox"/>	Back cane 8° bend 30.48-38.10cm (12-15")	NC
BC16-1518	<input type="checkbox"/>	Back cane 8° bend 38.10-45.72cm (15-18")	NC
BC16-1821	<input type="checkbox"/>	Back cane 8° bend 45.72-53.34cm (18-21") STD	NC
Select angle for the back cane:		90° <input type="checkbox"/> 95° STD <input type="checkbox"/> 100° <input type="checkbox"/> 105° <input type="checkbox"/>	
Specify height : ⁵			

TRANSIT TIE DOWN			
AC04	<input type="checkbox"/>	Transit Tie Down (Unoccupied)	NC
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt ⁶ Short belt 50-104 cm (20" - 41")	Upcharge
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt ⁶ Long belt 50-116 cm (20" - 46")	Upcharge
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt ⁶	Upcharge

⁶- Not available with depth adjustable back

Helio Kids FRONT Seat-to-floor Height - Please select your choice

Forks / Position		7,62cm (3") caster			10,16cm (4") caster			12,7cm (5") caster			15,24cm (6") caster			17,78cm (7") caster			20,32cm (8") caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"
4" 10,2 cm	P1	33.02cm 13"	35.56cm 14"	38.10cm 15"	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"												
	P2	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"									
5" 12,7 cm	P1																		
	P2				35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"									
	P3	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"						
	P4	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"	40.64cm 16"	43.18cm 17"	45.72cm 18"	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"			
7" 17,8 cm	P1									41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"				
	P2							41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"			
	P3							43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	
	P4	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"
	P5	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"	49.53cm 19.5"	52.07cm 20.5"	54.6cm 21.5"

Height available +/- 0,65cm (1/4") (see stem bolt and fork combination at the end of the order form)

CASTERS							
CA01-03	<input type="checkbox"/>	Composite wheel w/ PU tire	3" x 1"			NC	
CA01-04	<input type="checkbox"/>		4" x 1"			NC	
CA01-05	<input type="checkbox"/>		5" x 1"			NC	
CA01-06	<input type="checkbox"/>		6" x 1"			NC	
CA01-07	<input type="checkbox"/>		7" x 1"			NC	
CA01-08	<input type="checkbox"/>		8" x 1"			NC	
CA10	<input type="checkbox"/>		Pneumatic casters	6" x 1 1/4"			Upcharge
CA11	<input type="checkbox"/>			8" x 1 1/4"			Upcharge
CA14	<input type="checkbox"/>	Newton UltraCaster Composite	4" x 1 1/2"			Upcharge	
CA15	<input type="checkbox"/>		5" x 1 1/2"			Upcharge	
CA16	<input type="checkbox"/>		6" x 1 1/2"			Upcharge	
CA17	<input type="checkbox"/>		4" x 1 1/2"			Upcharge	
CA18	<input type="checkbox"/>	Newton UltraCaster Aluminum	5" x 1 1/2"			Upcharge	
CA19	<input type="checkbox"/>		6" x 1 1/2"			Upcharge	
CASTER OPTIONS							
CA12	<input type="checkbox"/>	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹				Upcharge	
CA13	<input type="checkbox"/>	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)				Upcharge	

1- Not available with WC-19

KIDS REAR Seat-to-floor Height

WHEEL SIZES	Height available	Specify desired height:
20" (451) wheel <input type="checkbox"/>	31.75cm / 12.5" to 43.81cm / 17.25"	<input style="width: 100%; height: 20px;" type="text"/>
22" (501) wheel <input type="checkbox"/>	34.29cm / 13.5" to 46.35cm / 18.25"	
24" (540) wheel <input type="checkbox"/>	36.83cm / 14.5" to 48.89cm / 19.25"	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1.27cm (1/2") of selected seat to floor height

Footrest							
SWING-IN/SWING-OUT FOOTREST							
FR03-01	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 34.29 - 43.18cm (13.5" - 17")			Set at: *		NC
FR03-02	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5")			Set at: *		NC
FR03-03	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 39.29 - 52.07cm (15.5" - 20.5")			Set at: *		NC
FR04-01	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 34.29 - 43.18cm (13.5" - 17")			Set at: *		NC
FR04-02	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5")			Set at: *		NC
FR04-03	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 39.37 - 52.07cm (15.5" - 20.5")			Set at: *		NC
FR15-01	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 34.29 - 40.64cm (13.5" - 16")			Set at: *		NC
FR15-02	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5")			Set at: *		NC
FR15-03	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 39.37 - 52.07cm (15.5" - 20.5")			Set at: *		NC
FR06	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Standard (no extension) ¹		19.05 - 27.94cm (7.5" - 11")	Set at: *		NC
FR07	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Short extension 2" ^{1,2}		24.13 - 33.02cm (9.5" - 13")	Set at: *		NC
FR08	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Medium extension 4" ^{1,2}		24.13 - 36.83cm (9.5" - 14.5")	Set at: *		NC
FOOTREST OPTIONS							
FR09	<input type="checkbox"/>	High mount (footplate attached to the hanger with clamps) ³		12.70 - 34.29cm (5" - 13.5")	Set at: *		Upcharge
FR10	<input type="checkbox"/>	Elevating legrest with calf pad, length adjustable ⁴ Left <input type="checkbox"/> Right <input type="checkbox"/>		35.56 - 55.88 cm (14"-22")	Set at:		Upcharge
FR99	<input type="checkbox"/>	Omit Footrest					NC

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount

* +/- 1/2" (1.27 cm)

4- The elevating legrest replace the footrest side selected; Not compatible with seat width of 12" and 13" and one piece footplate

FOOTPLATES			
FP01	<input type="checkbox"/>	Standard - Composite ¹	NC
FP02	<input type="checkbox"/>	Newton adjustable angle - Composite	Upcharge
FP03	<input type="checkbox"/>	Newton adjustable angle - Aluminum ²	Upcharge
FP04	<input type="checkbox"/>	Newton adjustable angle - Carbon fiber with rubber inserts	Upcharge
FP23	<input type="checkbox"/>	Newton adjustable angle extra deep - Aluminum	Upcharge
FP06	<input type="checkbox"/>	One piece adjustable angle flip-up - Aluminum	Upcharge
FP07	<input type="checkbox"/>	One piece adjustable angle flip-up - Carbon fiber with rubber inserts	Upcharge
FOOTPLATE ACCESSORIES			
FP08	<input type="checkbox"/>	Length adjustable heel loop (not available with one piece footplate)	Upcharge
FP09	<input type="checkbox"/>	Calf strap	Upcharge
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap	Upcharge

1- Not available with front mount 2 - The Newton adjustable angle (Aluminum) is the only available option with a 12 and 13" wide frame

SWING AWAY RESIDUAL LIMB SUPPORT			
FR27	<input type="checkbox"/>	Residual limb support - flat 15.24x20.32cm (6"x8")	Left <input type="checkbox"/> Right <input type="checkbox"/> Upcharge
FR28	<input type="checkbox"/>	Residual limb support - flat 20.32x24.4cm (8"x10")	Left <input type="checkbox"/> Right <input type="checkbox"/> Upcharge
FR29	<input type="checkbox"/>	Residual limb support - flat 25.4x25.4cm (10"x10")	Left <input type="checkbox"/> Right <input type="checkbox"/> Upcharge
FR30	<input type="checkbox"/>	Residual limb support - flat 25.4x35.56cm (10"x14")	Left <input type="checkbox"/> Right <input type="checkbox"/> Upcharge
FR32	<input type="checkbox"/>	Residual limb support - contour 15.24x20.32cm (6"x8")	Left <input type="checkbox"/> Right <input type="checkbox"/> Upcharge
FR33	<input type="checkbox"/>	Residual limb support - contour 20.32x24.4cm (8"x10")	Left <input type="checkbox"/> Right <input type="checkbox"/> Upcharge

PUSH HANDLES OPTIONS (ONE option at a time)			
PH01	<input type="checkbox"/>	Push handles	STD
PH05	<input type="checkbox"/>	Fold down push handles ¹	Upcharge
PH06-920	<input type="checkbox"/>	clamp on height adjustable push handles 22.86-50.80cm (9-20") ^{1,2}	Upcharge
PH06-923	<input type="checkbox"/>	clamp on height adjustable push handles 22.86-58.40cm (9-23") ^{1,2}	Upcharge
PH02	<input type="checkbox"/>	Omit push handles	NC

STROLLER HANDLE AND STABILIZER BAR (only available with push handle)			
PH04	<input type="checkbox"/>	Stroller handle	Upcharge
PH03	<input type="checkbox"/>	Newton folding stabilizer bar (not available with stroller handle)	Upcharge

1- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock; 2- Not compatible with depth adjustable back, 8" bend back canes and amputee axle plates

BACK UPHOLSTERY			
BH07	<input type="checkbox"/>	Slip-on Active-Mesh back upholstery	NC
BH08	<input type="checkbox"/>	Adjustable tension Active-Mesh back upholstery	Upcharge
BH99	<input type="checkbox"/>	Omit Back upholstery	STD

SEAT CUSHION			
SC01	<input type="checkbox"/>	5.08 cm (2") Cushion	Upcharge
SC02	<input type="checkbox"/>	7.62 cm (3") Cushion	Upcharge
SC99	<input type="checkbox"/>	Omit seat cushion	STD

Rear wheels

ONE ARM DRIVE ¹ Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims, permanent axles, 0° camber			
HR22	<input type="checkbox"/>	One arm drive - aluminum inner handrim ² 24"	Left <input type="checkbox"/> Right <input type="checkbox"/> Upcharge
HR60	<input type="checkbox"/>	One arm drive - plastic coated inner handrim ² 24"	Left <input type="checkbox"/> Right <input type="checkbox"/> Upcharge
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)	NC

REAR WHEELS Available			
RW01	<input type="checkbox"/>	Mag ²	20", 22" 24" Upcharge
RW02	<input type="checkbox"/>	Newton One - spoke wheel ³	20", 22" 24" NC
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel	20", 22", 24" Upcharge
RW04	<input type="checkbox"/>	Spinergy Spox Bik STD <input type="checkbox"/> Wht ⁴ <input type="checkbox"/> Blue ⁴ <input type="checkbox"/> Red ⁴ <input type="checkbox"/> Yel ⁴ <input type="checkbox"/>	22", 24" Upcharge
RW05	<input type="checkbox"/>	Spinergy LX Bik STD <input type="checkbox"/> Pur ⁴ <input type="checkbox"/> Wht ⁴ <input type="checkbox"/> Blue ⁴ <input type="checkbox"/> Red ⁴ <input type="checkbox"/> Yel ⁴ <input type="checkbox"/> Pink ⁴ <input type="checkbox"/> Grn ⁴ <input type="checkbox"/> Org ⁴ <input type="checkbox"/>	22", 24" Upcharge
RW07	<input type="checkbox"/>	Spinergy CLX ⁵ Bik STD <input type="checkbox"/> Pur ⁴ <input type="checkbox"/> Wht ⁴ <input type="checkbox"/> Blue ⁴ <input type="checkbox"/> Red ⁴ <input type="checkbox"/> Yel ⁴ <input type="checkbox"/> Pink ⁴ <input type="checkbox"/> Grn ⁴ <input type="checkbox"/> Org ⁴ <input type="checkbox"/>	24" Upcharge

TIRES Available			
TI01	<input type="checkbox"/>	Soft Urethane - 1 3/8"	20", 22" 24" STD
TI02	<input type="checkbox"/>	Soft Urethane - knobby 1 3/8 "	22", 24" Upcharge
TI04	<input type="checkbox"/>	Soft Urethane - Shox 1"	22", 24" Upcharge
TI05	<input type="checkbox"/>	Pneumatic - 1 3/8"	20", 22", 24" NC
TI08	<input type="checkbox"/>	Pneumatic with Airless insert - 1 3/8" ⁶	22", 24" Upcharge
TI06	<input type="checkbox"/>	Pneumatic HP - 1"	20", 22", 24" Upcharge
TI07	<input type="checkbox"/>	Pneumatic HP - Puncture resistant, black 1"	24" Upcharge

1- Offered with handrims on each side of the wheelchair; Select outer handrim in the handrims section; 2- Not compatible with HP tires; 3- Not compatible with Shox tires;

4- Colour spokes require an additional 10 days lead-time; 5- Only available with 1" tires; 6- Only available with Mag, Newton One and Newton Gravity rear wheels

Handrims

HANDRIM			Available		
HR01	<input type="checkbox"/>	Aluminum anodized	20", 22", 24"		STD
HR61	<input type="checkbox"/>	Aluminum black hard anodized ²	20", 22", 24"		Upcharge
HR02	<input type="checkbox"/>	Plastic coated ¹	20", 22", 24"		Upcharge
HR03	<input type="checkbox"/>	High friction coated ^{1,2}	22", 24"		Upcharge
HR07 to HR18	<input type="checkbox"/>	Newton Air Grip ² See on page 6 for color choice and price	20", 22", 24"		
HR28	<input type="checkbox"/>	Surge ^{1,2} (1 1/2" Oval handrim w/ Gripton high friction strip)	**Limited availability	22", 24"	Upcharge
HR05	<input type="checkbox"/>	Surge LT ^{1,2} (1 1/8" Oval handrim w/ Gripton high friction strip)	**Limited availability	22", 24"	Upcharge
HR04	<input type="checkbox"/>	Natural Fit ^{1,2}	**Limited availability	22", 24"	Upcharge
HR62	<input type="checkbox"/>	Tour H ^{1,2} (1 5/16" Oval handrim with high friction strip)		22", 24"	Upcharge
HR63	<input type="checkbox"/>	Nova H ^{1,2} (1 1/4" Oval handrim with high friction strip)		22", 24"	Upcharge
HR64	<input type="checkbox"/>	Optimum H ^{1,2} (1 1/8" heart-shaped handrim)		22", 24"	Upcharge
HR99	<input type="checkbox"/>	Omit handrim			NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels

Handrims option

ASSEMBLY POSITION					
HR29	<input type="checkbox"/>	Standard Handrim position Full tab			STD
HR030	<input type="checkbox"/>	Narrow Handrim position no cut Full tab			NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut ³ Short tab 2 positions			Upcharge
HR31	<input type="checkbox"/>	Super Narrow Handrim ³ Short tab 2 positions STD <input type="checkbox"/> Short tab 1 position <input type="checkbox"/>			Upcharge
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS					
HR23	<input type="checkbox"/>	Standard grip	**Limited availability		Upcharge
HR27	<input type="checkbox"/>	Super grip	**Limited availability		Upcharge
HR099	<input type="checkbox"/>	Omit thumb piece			NC


3- Only available with Aluminum and Air Grip handrims;

AXLE					
AX01	<input type="checkbox"/>	Quick release axle ¹			STD
AX02	<input type="checkbox"/>	Permanent axle			NC
AX03	<input type="checkbox"/>	Quad release axle ¹			Upcharge
AXLE PLATE					
AX04	<input type="checkbox"/>	Amputee adjustable axle plate ²			Upcharge

1- Not required with one arm drive 2 - May interfere with swing away armrest, depth adjustable back and limit seat to floor height available

REAR WHEEL CAMBER					
CM00	<input type="checkbox"/>	0° camber			STD
CM03	<input type="checkbox"/>	3° camber			NC
CM06	<input type="checkbox"/>	6° camber (not available with amputee axle plate)			Upcharge

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

X =	Most stable					Most responsive						
	2,5 1"	3,2 1/4"	3,8 1/2"	4,4 3/4"	5,1 2"	5,7 1/4"	6,3 1/2"	7 3/4"	7,6 3"	8,3 1/4"	8,9 1/2"	9,5 3/4"
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

If not indicated, default setting is at 3,8 cm (1.5")

Possible conflict between a seat depth under 40,6 cm (16") and a center of gravity over 6,3 cm (2.5").

Possible conflict between a 66 cm (26") rear wheel and a center of gravity over 6,3 cm (2.5").

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

WHEEL LOCK					
WL08	<input type="checkbox"/>	Newton Wheel lock - push to lock			STD
WL16		Aluminum push to lock with extension			Upcharge
WL02	<input type="checkbox"/>	Pull to lock			NC
WL04	<input type="checkbox"/>	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)			Upcharge
WL13 WL14	<input type="checkbox"/>	6" removable extension handle			Upcharge
WL05	<input type="checkbox"/>	Scissors Lock			Upcharge
WL06	<input type="checkbox"/>	Attendant Lock (additional set of wheel lock mounted to the back)			Upcharge
WL15	<input type="checkbox"/>	Unilateral wheel lock - pull to lock ¹ Left <input type="checkbox"/> Right <input type="checkbox"/>			Upcharge
WL99	<input type="checkbox"/>	Omit wheel locks			NC

1- For more information about incompatibilities and restriction, refer to the table at page 6 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_intl/support-and-education/how-to-documents

Armrest

HEIGHT ADJUSTABLE "T" ARMREST ¹				
AR13-610	<input type="checkbox"/>	25.4cm (10") Desk length armpad ⁵	15.25 - 25.40cm (6" - 10")	Upcharge
AR03-812	<input type="checkbox"/>	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")	Upcharge
AR03-1014	<input type="checkbox"/>	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")	Upcharge
AR14-610	<input type="checkbox"/>	35.56cm (14") Full length armpad ⁵	15.25 - 25.40cm (6" - 10")	Upcharge
AR04-812	<input type="checkbox"/>	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")	Upcharge
AR04-1014	<input type="checkbox"/>	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")	Upcharge
LOCKING FLIP-UP ARMREST ²				
AR07	<input type="checkbox"/>	Locking Flip-Up Armrest		Upcharge
ARMREST OPTIONS				
AR12	<input type="checkbox"/>	35.56cm (14") Full Length Armpad - Gel Ovation (pair) ⁴		Upcharge
AR99	<input type="checkbox"/>	Omit armrest		NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- May interfere with rigid back hardware; Only available with 18"-21" straight aluminum back canes; Not available with back upholstery

3- Not available with locking flip-up armrest; 4- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (not available with "T" Armrest)				
SG01	<input type="checkbox"/>	Plastic side guards (removable)		Upcharge
SG03	<input type="checkbox"/>	Carbon fiber side guards (removable)		Upcharge
SG06	<input type="checkbox"/>	Plastic fender side guards - straight bracket (removable)		Upcharge
SG08	<input type="checkbox"/>	Plastic fender side guards - offset bracket (removable) ¹		Upcharge
SG04	<input type="checkbox"/>	Carbon fiber fender side guards (removable)		Upcharge
SG99	<input type="checkbox"/>	Omit side guard		NC
SIDE GUARDS OPTION				
SG07	<input type="checkbox"/>	Side guards clip (pair)		NC

1- Only compatible with 0° camber; Increase distance between fender side guards by 1 1/8";

POSITIONING BELT				
SB01-S	<input type="checkbox"/>	Velcro adjustable - Standard 88.9cm (35")		Upcharge
SB01-L	<input type="checkbox"/>	Velcro adjustable - Long 152.40cm (60")		Upcharge
SB02-S	<input type="checkbox"/>	Auto buckle - Standard 88.9cm (35")		Upcharge
SB02-L	<input type="checkbox"/>	Auto buckle - Long 152.40cm (60")		Upcharge
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Medium 142.24cm (56")		Upcharge
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Long 154.94cm (61")		Upcharge
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Small 129.54cm (51")		Upcharge
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Medium 142.24cm (56")		Upcharge
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Long 154.94cm (61")		Upcharge
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Small 129.54cm (51")		Upcharge
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Medium 142.24cm (56")		Upcharge
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Long 154.94cm (61")		Upcharge
SB03	<input type="checkbox"/>	Hardware attachment for 2 extra points (1 pair)		Upcharge
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 73.66 (18")		Upcharge
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 93.98cm (24")		Upcharge
SB99	<input type="checkbox"/>	Omit positioning belt		STD

HEADREST ¹				
AC01-S	<input type="checkbox"/>	Profiled Headrest - Small 8.8x15.2cm (3.5"X6")		Upcharge
AC01-M	<input type="checkbox"/>	Profiled Headrest - Medium 11.4x17.8cm (4.5"X7")		Upcharge
AC01-L	<input type="checkbox"/>	Profiled Headrest - Large 14x2.3cm(5.5"X8")		Upcharge
AC28	<input type="checkbox"/>	Headrest mounting hardware with removable non folding stabilizer bar ¹		Upcharge

1 - Only available with push handles - Not available with a stroller handle or a folding stabilizer bar


ACCESSORIES													
AC03	<input type="checkbox"/>	Anti-tippers								Upcharge			
AC29	<input type="checkbox"/>	Swing-away anti-tipper ¹	Left <input type="checkbox"/>	Right <input type="checkbox"/>						Upcharge			
AC27	<input type="checkbox"/>	Tip assist ²								Upcharge			
AC06	<input type="checkbox"/>	Cane and crutch holder	Left <input type="checkbox"/>	Right <input type="checkbox"/>						Upcharge			
AC07	<input type="checkbox"/>	Oxygen tank holder ⁴ - D and E Type - installed on the left								Upcharge			
AC08	<input type="checkbox"/>	Spoke guards ³	Available: 20", 22", 24", 26"							Upcharge			
AC09	<input type="checkbox"/>	Footrest neoprene impact guard (pair) - Black								Upcharge			
AC10	<input type="checkbox"/>	Seat pouch								Upcharge			
AC30	<input type="checkbox"/>	Newton multi-purpose pouch								Upcharge			
AC23	<input type="checkbox"/>	BodyPoint Mobility Bag								Upcharge			
AC26	<input type="checkbox"/>	Newton Backpack								Upcharge			
AC12	<input type="checkbox"/>	Velcro Closing Strap								Upcharge			
AC13	<input type="checkbox"/>	Tool kit								Upcharge			
AC16	<input type="checkbox"/>	Motion Composites Polo Shirt	W/S <input type="checkbox"/>	W/M <input type="checkbox"/>	W/L <input type="checkbox"/>	M/S <input type="checkbox"/>	M/M <input type="checkbox"/>	M/L <input type="checkbox"/>	M/XL <input type="checkbox"/>	M/XXL <input type="checkbox"/>	M/XXXL <input type="checkbox"/>		Upcharge

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 34,3cm (13 1/2") and 50,2cm (19 3/4")
 2 - Cannot be combined with Anti Tippers
 3 - Not available with Mag wheels and Spinergy; 4- Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down

Colors																					
Options ¹	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Upcharge
Frame ²	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	NC
FC99	<input type="checkbox"/>	Custom Color	Specify PPG Paintit color:																additional 10 days lead-time	Upcharge	

1- Finish between Airgrip and frame might be different 2- With carbon fiber crossbrace;

Fork and Stem bolt combinations



Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
22 in	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Notes :