

# Jenx Junior+



SEATING SYSTEM FOR MILD,  
MODERATE AND COMPLEX NEEDS

Distributed by  
**medifab**

SCRIPTING AND ORDERING GUIDE


October 2023

# Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

## User Measurements

A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	






Product Specifications	Size 1	Size 2	Size 3
Average age	18 Months – 5 Years	3 – 8 Years	5 – 12 Years
Weight capacity (kg)	50	50	50
Seat width (cm)	15* – 37.5	15* – 37.5	15* – 37.5
Seat depth (cm)	21 – 30.5	24 – 35.5	26 – 44
Chest support width (cm)	18.5 – 28	18.5 – 28	21 – 30
Seat-top of backpad height (cm)	34.5 – 51**	34.5 – 51**	41 – 57.5**
Footplate angle (°)	-35 – 5	-35 – 5	-35 – 5
Footplate height (cm)	19.5 – 27	23.5 – 35	25 – 40.5
Backrest recline (°)	-15 – 15	-15 – 15	-15 – 15
Tilt range (°)	30 (Posterior) – 10 (Anterior)	30 (Posterior) – 10 (Anterior)	30 (Posterior) – 10 (Anterior)


\* 15cm achieved only when using the hip reducing pads on Page 3.

\*\* 51cm back height in Size 1 + 2 and 57.5cm in Size 3 only achieved with back post extension on Page 4.

# Colour Options


Red	Lime	Blue
		
Fabric Cover, Red -022	Fabric Cover, Lime -054	Fabric Cover, Blue -066
Wipeable Cover, Red -W22		Wipeable Cover, Blue -W66

Medifab recommends that items outlined in blue are selected as standard package.






Therapists Normally Choose


# Base

Y-Base
<p>Jenx Y-Base with height adjustment of 34 - 76cm, lockable swivel castors, tilt-in-space, 30° posterior, 10° anterior 3615-8740-000</p> 



# Seat Unit


Size 1	Size 2	Size 3
<p>Contoured PU seat with removeable abduction pad and memory foam seat pad, adjustable and flip-away PU hip pads, angle-adjustable footplate, 4-point hip belt, height-adjustable PU back pad, an independently adjustable PU lumbar/sacral pad, arm supports, tray table and push handle.</p>		
<p>Age range: 18 months–4 years. Max user weight: 50kg.</p>	<p>Age range: 3–8 years. Max user weight: 50kg.</p>	<p>Age range: 6–12 years. Max user weight: 50kg.</p>
3610-1001-000	3610-1002-000	3610-1003-000
		
Covers Required:	Covers Required:	Covers Required:
Seat Cover Fabric 3612-0421-XXX	Seat Cover Fabric 3612-0423-XXX	Seat Cover Fabric 3612-0423-XXX
Back Cover Fabric 3612-2421-XXX	Back Cover Fabric 3612-2421-XXX	Back Cover Fabric 3612-2423-XXX
Hip Cover Fabric 3612-1072-XXX	Hip Cover Fabric 3612-1072-XXX	Hip Cover Fabric 3612-1072-XXX
Cover set (wipe clean only) 3612-0521-WXX	Cover set (wipe clean only) 3612-0522-WXX	Cover set (wipe clean only) 3612-0523-WXX

# Hip Width Reducing Pads


<b>Jenx Junior+ Hip width reducing pads - pair.</b>	
Reduces Hip Width by 2.5cm on each side. 3612-1040-000	

# Head Supports



Multigrip - 1 Link	Multigrip - 2 Link	Multigrip - 3 Link	
Single link wings for some degree of lateral support. 33cm W x 5cm D x 7.5cm H 5356-4601-000	2 Links in each wing for extra support and positioning. 47.5cm W x 5cm D x 12cm H 5356-4602-000 	3 adjustable links in each wing for ultimate lateral headrest control and comfort. 62cm W x 5cm D x 12cm H 5356-4603-000	
Wipeable Cover (Optional) 5356-4671-000	Wipeable Cover (Optional) 5356-4672-000	Wipeable Cover (Optional) 5356-4673-000	

<b>Oval Head Support</b>	
19cm W x 10cm H 5356-4611-XXX	

# Shoulder Protractors

Size 0	Size 1	Size 2	
Can be used on all three size seat units			
Pad size 13cm high x 13.5cm deep 3612-2750-000	Pad size 21cm high x 14.5cm deep 3612-2751-000	Pad size 23.5cm high x 17.5cm deep 3612-2752-000	
Covers Required 3612-2760-XXX	Covers Required 3612-2761-XXX	Covers Required 3612-2762-XXX	

# Trunk Supports

Flexi Thoracics (Pair)			
<b>Size 1</b> 18cm long x 8cm high 3615-3251-000	<b>Size 1 Long</b> 22cm (L) x 8cm (H) 3615-3252-000 	<b>Size 2</b> 22cm long x 11cm high 3615-3253-000	
Covers Required 3615-3271-XXX	Covers Required 3615-3272-XXX	Covers Required 3615-3273-XXX	

# Waistcoat Harness

Size 1	Size 2	Size 3
<b>A</b> 5cm <b>B</b> 28cm <b>C</b> 12cm <b>D</b> 13.5cm 3615-6701-000	<b>A</b> 5cm <b>B</b> 32cm <b>C</b> 14cm <b>D</b> 15.5cm 3615-6702-000	<b>A</b> 6.5cm <b>B</b> 37cm <b>C</b> 17.5cm <b>D</b> 18cm 3615-6703-000



# Back Post Extension

7.5cm Back Post Extension for extra backrest height  
 3610-2480-000



# Abduction Block

Size 0	Size 1
Length: 9cm Height above Seat: 8.5cm Overall height: 14.5cm Width Front: 9cm Width Rear: 7cm 3615-1650-000	Length: 18cm Height above Seat: 7.5cm Overall height: 15.5cm Width Front: 10cm Width Rear: 8.5cm 3615-1651-000



# Leg Positioning


**Leg Control Pads**

Control leg position independently of pelvic position, flip completely out of the way to facilitate smooth and easy transfers.  
 12.5cm H x 7.5cm D  
 3615-1861-000


**Covers Required**  
 3615-3171-XXX




# Sandals

Size 1	Size 2	Size 3	
17cm long x 8cm wide 5356-6081-000	22cm long x 9.5cm wide 5356-6082-000	27.5cm long x 10.5cm wide 5356-6083-000	


# Sandal Raiser Blocks

Size 1	Size 2	Size 3	
Kit of 4x1.5cm blocks 5356-6091-000	Kit of 4x1.5cm blocks 5356-6092-000	Kit of 4x1.5cm blocks 5356-6093-000	

# Ankle Snugs


Size 1	Size 2 & 3	
Inner Circumference: 12cm - 22.5cm 5356-6131-000	Inner Circumference: 16.5cm - 29.5cm 5356-6132-000	

# Foot Plate Pad


Foot Plate Pad	
Footplate pad not compatible with sandals. 4cm Thick 3615-5301-000	

# Accessories for Arm Support & Tray

## Tray Infill

Great for attaching suction toys.	
<b>Size 1</b> 3610-7180-000	
<b>Size 2</b> 3610-7182-000	

## Half Tray Padding

Comfort for forearms and elbows, whilst still enjoying some hard tray surface	
<b>Size 1</b> 3612-7441-000	
<b>Size 2</b> 3612-7442-000	


## Full Tray Padding

Added comfort for hands also	
<b>Size 1</b> 3612-7451-000	
<b>Size 2</b> 3612-7452-000	

## Grab Bar

Used for children who need to fix one hand to increase functionality. 54cm W x 9.5cm D x 7 cm H 3615-7200-000	
---	---

## Activity Frame

Great for encouraging head control and hand to eye co-ordination. 54cm W x 16cm D x 35 cm H 3615-7475-000	
---	---

## Elbow Blocks

Used to encourage arms into a midline position and/or to block excessive retraction. 10cm W x 14cm H 3615-7350-000	
--	---

## Arm Support Covers

3615-5671-XXX	
---------------	---

# medifab

Distributed by Medifab Ltd

22 Detroit Drive  
Rolleston 7675  
New Zealand

0800 543 343

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: JenxJunior+\_20231025