

Flux Walker

The Flux is a lightweight and sturdy posterior walking aid designed for children and young teens on the move. The walker is offered with a range of accessories for function and promotion of upright walking patterns. The perfect companion to children on the move, both indoors and outside!

Dependable support for children & adults



Easy to adjust hand supports allow Flux to grow with the user through tool free adjustment for height.

Opens out in walking direction and optically and physically free from interfering structural parts

Anti roll back mechanism included as standard with added benefit of being turned on and off depending on user.

Friction brakes and bumper wheels included as standard

Folding mechanism with **safety lock** eases transport and storage

Which size?



Standard Package Includes:

- Standard Handlebar
- Front Swivel Wheels
- Polyurethan Tyres
- Bumper Wheels
- Friction Brakes
- Anti-Reverse Mechanism



ACCESSORIES

Fold-Up Seat

Size 1-2 5640-0201-000
Size 3 5640-0204-000



Multi Adjustable Hip Pad

5640-1900-000



Dynamic Hip Control

Soft Suspension (Yellow) 5640-1952-000
Hard Suspension (Red) 5640-1953-000
For Orthoses 5640-1951-000



Add Hip Belt (Required)

Small, Hip 53-83 cm 5640-6201-000
Large, Hip 65-115 cm 5640-6202-000



Forearm Supports, 2 Axial Adjust

Small 5640-5500-000
Standard 5640-5501-000



Universal Handlebars

Size 1 5641-5751-000
Size 2 5642-5752-000
Size 3 5643-5753-000



Not compatible with forearm supports

Driving Brake & Wheel Lock

5640-8451-000



Anti-Tipper

5640-8421-000



Anti-Tipper & Castor

5640-8425-000



Hand Disk

5640-5550-000



Swivel Lock

5640-8440-000



Size	Code	Colour	Standard Grip Height	Universal Grip Height	Frame Width	Frame Length	Frame Weight	Max User Weight
1	5641-0001-000	Luminous Red	48-61 cm	43-68 cm	64 cm	65 cm	6.5 kg	60 kg
2	5642-0001-000	Sky Blue	57-81 cm	51-88 cm	69 cm	80 cm	6.8 kg	75 kg
3	5643-0001-000	Gentian Blue	69-87 cm	62-94 cm	74 cm	88 cm	7.1 kg	75 kg