

# FLUX Walker



## VERSATILE POSTERIOR WALKING AID

Distributed by  
**medifab**

NEW ZEALAND  
SCRIPTING AND ORDERING GUIDE

# Client Details

Client/User Name		Notes
Date		
NDIS/Funding		

# Measurements

	Size 1	Size 2	Size 3	Size 4
Standard grab rail grip height	48-61cm	57-81cm	69-87cm	82-102cm
Standard grab rail grip width	41cm	46cm	51cm	51cm
Universal grip grab rail grip height	43-68cm	51-88cm	62-94cm	75-102cm
Universal grip grab rail grip width	29-41cm	34-46cm	39-51cm	39-51cm
Total frame width	64cm	69cm	74cm	78cm
Total frame length	65cm	80cm	88cm	91cm
Total load capacity	60kg	75kg	75kg	100kg
Front / rear wheel diameter	15cm / 20cm	15cm / 20cm	15cm / 20cm	20cm / 30cm

# Package & Size Options

Package Includes	Size 1	Size 2	Size 3	Size 4
<ul style="list-style-type: none"> <li>Standard Handlebar</li> <li>Front Swivel Wheels</li> <li>Polyurethan Tyres</li> <li>Bumper Wheels</li> <li>Friction Brakes</li> <li>Anti-Reverse Mechanism</li> </ul>	5641-0001-000	5642-0001-000	5643-0001-000	5644-0001-000
				
	Luminous Red	Sky Blue	Gentian Blue	Grey Aluminium

# Accessories

## Anti-Tippers without Castors

**Size 1-3**  
5640-8421-000



**Size 4**  
5640-8424-000



## Anti-Tippers with Castors

**Size 1-3**  
5640-8425-000




**Size 4**  
5640-8428-000



## Fold Up Seat

**Size 1-2**  
5640-0201-000

**Size 3-4**  
5640-0204-000



## Carrier

**Size 1**  
5640-8061-000

**Size 2**  
5640-8063-000



## Sling Seat

**Small**  
Hip Size: 57-77cm  
5640-1921-000

**Large**  
Hip Size: 73-103cm  
5640-1922-000



## Sling Seat Rail

**Size 1**  
5640-1931-000

**Size 2**  
5640-1932-000

**Size 3**  
5640-1933-000



## Flux Basket

5640-8051-000



## Hand Disc

5640-5550-000



## Driving Brake

**Size 1-3**  
Includes wheel lock.  
5640-8451-000

**Size 4**  
Includes wheel lock.  
5640-8454-000




## Swivel Lock



To hold castors straight  
when brisk walking.  
2880-8440-000



# Accessories Continued

Dynamic Hip Control (Hip Belt Required)		Hip Belt (Required)	
<b>Suspension, soft</b> 5640-1952-000		<b>Small</b> Hip Size: 53-83cm 5640-6201-000	
<b>Suspension, hard</b> 5640-1953-000		<b>Large</b> Hip Size: 65-115cm 5640-6202-000	
<b>Kit for orthoses</b> 5640-1951-000			

Universal Handlebar			
<b>Size 1</b> 5641-5751-000	<b>Size 2</b> 5642-5752-000		
<b>Size 3</b> 5643-5753-000	<b>Size 4</b> 5644-5754-000		

Forearm Supports - 2 Axial Adjust		Hip Pad	
22 x 10.5cm 5640-5501-000		Multi-Adjustable 5640-1900-000	

## Notes

---



---



---



---



---



---



---



---



---



---

# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Flux\_20231215\_NZ