

# Duro Stroller



## PAEDIATRIC ADAPTIVE STROLLER

AUSTRALIAN  
SCRIPTING AND ORDERING GUIDE

Distributed by  
**medifab**

January 2024

# Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

# Product Details



Specifications	Size 1	Size 2
<b>Order Code</b>	2780-0001-052	2780-0002-052
<b>Seat Specifications</b>		
Seat width	23-34cm	28-39cm
Seat depth	26-40cm	31-45cm
Seat angle	16 or 20°	16 or 20°
Back height	45-63cm	61-83cm
Back angle	78-120°	78-120°
Knee angle	85-160°	85-160°
Lower leg length	20-44cm	20-44cm
Folded size (L x W x H)	84 x 63 x 32cm	84 x 68 x 38cm
<b>Base Specifications</b>		
Max weight	75kg	75kg
Front wheel diameter	20cm	20cm
Rear wheel diameter	30cm	30cm

### Stroller Package Includes:


- Swivelling Front Wheels
- Foot Support Bracket, Angle Adjustable, Flip Up
- Foot Support With Heel Stop
- Accessories Bag (Max Weight: 3kg)
- Tie Down Kit
- Retainer for Accessories
- Hip Guides (Low Version)
- **Note: Upholstery Not Included In Package**






# Upholstery Colour Options

Turquoise		Black	
Seat Upholstery (required)			
<b>Size 1</b> 2780-8801-349		<b>Size 1</b> 2780-8801-346	
<b>Size 2</b> 2780-8802-349		<b>Size 2</b> 2780-8802-346	

# Head Supports

Head Support Pads	
<b>Size 1</b> 11cm Deep x 16cm High 2880-4471-346	
<b>Size 2</b> 13cm Deep x 18cm High 2880-4472-346	

Occipital Head Support	Standard Lateral Head Support	Extended Lateral Head Support	
<b>Merino/Spacer Fabric</b> 2010-4110-099	<b>Small (17cm Internal width)</b>		
	<b>Merino/Spacer Fabric</b> 2010-4210-099		<b>Merino/Spacer Fabric</b> 2010-4310-099
<b>Wipe Down Cover</b> 2010-4110-209	<b>Wipe Down Cover</b> 2010-4210-209		<b>Wipe Down Cover</b> 2010-4310-209
	<b>Large (21cm Internal width)</b>		
<b>Wipe Down Cover</b> 2010-4110-209	<b>Merino/Spacer Fabric</b> 2010-4211-099		<b>Merino/Spacer Fabric</b> 2010-4311-099
	<b>Wipe Down Cover</b> 2010-4211-209		<b>Wipe Down Cover</b> 2010-4311-209

# Lateral Supports



Standard Lateral Supports		
Height and angle adjustable. All sizes are compatible with both Duro sizes.		
<b>Size 1</b> 8cm(H) x 10cm(L) 2880-3010-346	<b>Size 2</b> 10cm(H) x 12cm(L) 2880-3020-346	<b>Size 3</b> 12cm(H) x 13cm(L) 2880-3030-346





# Lateral Supports (continued)

Spex Teardrop Lateral Supports		Spex Broad Lateral Supports	
13.5cm(H) x 7.5cm(L) 2880-3008-500		10.5cm(H) x 13cm(L) 2880-3008-310	



# Harness Options


Spex Centrepoint Harness		5-Point Harness	
<b>X-Small</b> Recommended for size 1 Duro Total pad dimensions: 26cm(H) x 22cm(W) 1409-6660-017		With front pull adjusters and neck padding 2880-6511-346	
<b>Small</b> Recommended for size 1 & 2 Duro Total pad dimensions: 36cm(H) x 32cm(W) 1409-6661-017			
<b>Medium</b> Recommended for size 2 Duro Total pad dimensions: 36cm(H) x 32cm(W) 1409-6662-017			

# Hip Belt Options


Spex Push Button		Pen Popper	
<b>Small</b> Pad width 36cm 1405-6211-028		<b>Small</b> Pad width 36cm 1405-6211-025	
<b>Medium</b> Pad width 46cm 1405-6212-028		<b>Medium</b> Pad width 46cm 1405-6212-025	

# Tray & Grab Rail Options

Therapy Tray		Therapy Tray Infill	
Cannot be used with Grab Rail			
<b>To fit Size 1</b> 2880-7111-000	<b>To fit Size 2</b> 2880-7112-000		
		For smaller trunk cut out. Fits both tray sizes. 2880-7180-000	

Grab Rail	
Cannot be used with therapy tray	
<b>To fit Size 1</b> 2880-7251-346	<b>To fit Size 2</b> 2880-7252-346
	


# Canopy

Includes mounts and rain cover	
<b>Size 1</b> 2780-7501-346	
<b>Size 2</b> 2780-7502-346	

# Hand Strap

Wrist safety tether strap	
2005-6930-000	

# Swivel Lock

To hold castors straight when brisk walking	
2880-8440-000	

# medifab

Distributed by Medifab Ltd

26 Pardoe Street  
East Devonport  
Tasmania 7310

+61 3 6427 0105

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY: