

Duro Stroller

Whether navigating city footpaths or the outdoors terrain, Duro is a stable stroller solution with a great range of support and positioning accessories that will enable you to embrace the outdoors.

Configurable rehab stroller for everyday challenges



Extremely agile with great manoeuvring due to compact mobility base

Adjustable back angle

2 options for seat angle

Generous growth capacity provides long-term use for growing users

Individually customisable to suit requirements, with proven accessories from the BINGO Evo Stroller

Very stable stroller option

Adjustable knee angle



Super compact folding makes the Duro easy to transport and stow in the boot of your car!



Reclining back support means you can relax back in the Duro when being out and about gets tiring!

Package Includes:

- Hoggi DURO Rehab-Buggy,
- Hip support, aluminium
- Seat and backrest upholstery
- Spex hip belt
- Canopy
- PU wheels with black rim, swivelling front wheels
- Supporting brackets
- Foot rest bracket angle adjustable, flip up type
- Foot rest with heel stop
- "Soft" accessories bag
- Tie down kit - according to ISO 7176-19
- Retainer for accessories grab rail or therapy tray



ACCESSORIES

Head Support Pads
1 2880-4471-346
2 2880-4472-346



Grab Rail
1 2880-7251-346
2 2880-7252-346



Standard Lateral Head Support
Sml 2010-4210-XXX
Lrg 2010-4211-XXX



Extended Lateral Head Support
Sml 2010-4310-XXX
Lrg 2010-4311-XXX



Occipital Head Support
2010-4110-XXX



Head Support Fabric Options:

- Merino/Spacer -099
- Wipe Down -209




Std Lateral Trunk Support
1 2880-3010-346
2 2880-3020-346
3 2880-3030-346

Contoured Lateral Trunk Support
2006-3030-099
Teardrop Lateral Trunk Support
1 2006-3110-099
2 2006-3120-099

Wonderseat Padded Centrepoint Harness
1 1409-6661-017
2 1409-6662-017



5-Point Anterior Chest Support
2880-6511-000



Swivel Lock
2880-8440-000



Therapy Tray
1: 2880-7111-000
2-3: 2880-7112-000
Infil: 2880-7180-000



Full range of accessories available in the Duro Stroller Scripting and Ordering Guide.

Size	Code	Weight Limit	Seat Width	Seat Depth	Seat Angle	Back Support Height	Back Support Angle	Lower Leg Length	Wheel Diameter	Folded Size (L x W x H)
1	Turquoise - 2780-8801-349	75 kg	23 cm	26 cm	16° or 20°	45 cm	78	20 cm	20 cm front, 30 cm rear	84x63x32 cm
	Black - 2780-8801-346		34 cm	40 cm		63 cm	120°	44 cm		
2	Turquoise - 2780-8802-349	75 kg	28 cm	31 cm	16° or 20°	61 cm	78	20 cm	20 cm front, 30 cm rear	84x68x38 cm
	Black - 2780-8802-346		39 cm	45 cm		83 cm	120°	44 cm		