


Jenx Junior+

Comfortable, contoured and stable modular seating system, designed for children with mild to complex posture support needs. Outstanding levels of adjustment with a wide variety of support options.

Functionality, flexibility and style




Carefully designed to draw attention to the child rather than the equipment itself

Flip-up angle adjustable foot support to accommodate dorsi or plantar flexion


Pelvic control system with independently adjustable contoured seating

All support pads and accessories are designed to **quickly & easily flip away** without being removed



Standard Package Includes

Contoured PU seat with removable abduction pad and memory foam seat pad, adjustable and flip-away PU hip pads, angle-adjustable footplate, 4-point hip belt, height-adjustable PU back pad, an independently adjustable PU lumbar/sacral pad, arm supports, tray table and push handle.



Y-Base

- Versatile and flexible, with lockable swivel castors
- High-low range of 25.8-75 cm seat height from floor
- Tilt from 30° posterior - 10° anterior
- Base footprint of 65 x 78 cm

3615-8740-000



Low seat-to-floor height for peer interaction.




Ease of transfer with tilt and swing-away supports.




Open design for comfort and support.

ACCESSORIES


Multigrip Head Support
 1 Link 5356-4601-000
 2 Link 5356-4602-000
 3 Link 5356-4603-000




Waistcoat Harness
 1 3615-6701-000
 2 3615-6702-000
 3 3615-6703-000




Abduction Block
 0 3615-1650-000
 1 3615-1651-000




Shoulder Protractors
 0 3612-2750-000
 1 3612-2751-000
 2 3612-2752-000




Flexi Lateral Supports
 1 3615-3251-000
 1 (Long) 3615-3252-000
 2 3615-3253-000




Activity Frame
 3615-7475-000




Tray Infill
 1 3610-7180-000
 2 3610-7182-000




Tray Elbow Blocks
 3615-7350-000




Back Post Extension
 3610-2480-000




Hip Width Reducing Pads
 3612-1040-000



Wipeable Cover, Blue
 1 3612-0521-W66
 2 3612-0522-W66
 3 3612-0523-W66



Wipeable Cover, Red
 1 3612-0521-W22
 2 3612-0522-W22
 3 3612-0523-W22



Size	Code	Approx Age Range	Max User Weight	Back Angle	Seat Width	Seat Depth	Back Support Height	Lower Leg Length
1	3612-0001-000	18 mths-5 yrs	50 kg	15° recline	15*-37.5 cm	21-30.5 cm	34.5-51 cm**	19.5-27 cm
2	3612-0002-000	3-8 yrs				24-35.5 cm	34.5-51 cm**	25-35 cm
3	3612-0003-000	5-12 yrs				26-44 cm	41-57.5 cm**	25-40.5 cm

* 15cm achieved only when using the hip reducing pads.

** 51cm back height in Size 1 + 2 and 57.5cm in Size 3 only achieved with back post extension.