

Ella Stander

A multiplicity of standing options and huge growth range maximises your investment in this versatile aid for individuals with moderate to severe disabilities.

Versatility for Prolonged Standing



Highly durable frame:
easy to assemble and adjust

Easy care upholstery
can be removed for machine washing

Stepless support
adjustments ensure accurate and comfortable positioning

Individual foot and knee adjustments

Can be set up for **supine or prone standing**: change over set-up with only 3 steps

Knee and ankle flexion adjustment for:

- Angles in all directions (including abduction)
- Height
- Width
- Depth



Good hoist access and swing-away side supports to facilitate easy transfers for supine standing



Prone set-up

Supine set-up

ACCESSORIES

Shoulder Protractors
5220-8460-000



Harness Fitting Kit
5220-8460-000



Swing Away Laterals
5220-XXXX-000



Tray with Bowl
5220-8460-000



Knee Support Rotation Module
5220-7200-000




Rotational Lumbar Support
S 5220-2601-000
M 5220-2602-000
L 5220-2603-000



ELECTRIC POWER ASSIST

Assisted 90° electric tilt function makes ease of adjustment, reduces manual handling and improves client experience.



EASY TO ORDER

Comes standard as fully complete package. Includes: headrest, tray, body supports and knee cuffs and offset brackets.

Size	Code	Max User Weight	Approx User Height	Footplate to Axilla	Chest / Hip Width	Foot to Patella	Foot to Pelvis	Abduction (Max Angle)	Kneepads to Footrest
1	5221-0020-000	30 kg	60-95 cm	43-68 cm	12-20 cm	14-45 cm	30-54cm	50°: 60 cm* 42°: 95 cm*	0-25°
2	5222-0020-000	40 kg	90-130 cm	57-97.5 cm	17-25 cm	18-55 cm	47.5-78.5cm	44°	0-25°
3	5223-0020-000	55 kg	120-155 cm	80-128 cm	19-30 cm	20-65 cm	50-104 cm	40°	0-25°
4	5224-0020-000	80 kg	130-180 cm	90-150 cm	25-41 cm	25-85 cm	62.5-122.5cm	40°	0-25°

*User Height