

# Q500M

POWERCHAIR SCRIPTING AND ORDERING GUIDE



Distributed by **medfab**

## CLIENT DETAILS

Client/User Name	
Assessor	
Funding Body	
User Weight	
Date	
Notes	

### Q500M BASE PACKAGE

- R-Net Expandable Electronics
- Lights
- 13" Drive Wheels
- 7" Castors
- Transit Tie Downs
- 4 Pole Motors
- Q500M Mini - 50Ah Gel Batteries
- Q500M - 60Ah Gel Batteries
- Battery Charger
- Standard Lap Belt

## BASE SIZE

Q500M	Q500M Mini
-------	------------

## SEAT SELECTION

Seat Type	
Sedeo Lite Seat Top (N/A with recline)	(Width: 15" - 20" Depth: 16" - 22")
Sedeo Pro Short Seat Top	(Width: 16" - 22" Depth: 14" - 18")
Sedeo Pro Long Seat Top	(Width: 16" - 22" Depth: 17" - 22")

## SEAT SIZE SELECTION

Seat Width:	Seat Depth:
-------------	-------------

## SEAT TO FLOOR HEIGHT SELECTION

(If you've selected the Q500M Mini please choose size 17 - 19)

16	16.5	17	17.5	18
18.5	19	Other:		

## POWER SEAT FUNCTIONS

Elevate 12"	Tilt 50°
Power Recline 130°	

## LEGREST OPTIONS

Legrest Type	
Swing Away 70°	Swing Away Manual Elevating
Swing Away Power Elevating	Fixed Centre Mount
Powered Center Mount	

### Stump Support

Left	Right
------	-------

## ARM SUPPORT OPTIONS

### Arm Support Assembly

Sedeo Lite Cantilever Flip Back Armrests (N/A with PWR recline)	Single Post Flip Back Armrests
Sedeo Pro Cantilever Flip Back Armrests	Dual Locking Armrest Post for Cantilever Flip Back

### Armrest Pads

Sedeo Lite Cantilever Armrest Pads 380mm (N/A with PWR recline)	Sedeo Pro Single Post Contour Arm Pads 300mm
Sedeo Pro Single Post Contour Arm Pads 400mm	Comfort Gel 2" Armrest Pads
Comfort Foam 2" Armrest Pads	Gel Arm Troughs Medium
Gel Arm Troughs Large	Spex Comfi Wide Armrests

## JOYSTICK OPTIONS

### Joystick

Rnet Joystick Controller Standard LED	Rnet Joystick Controller LCD Screen CJSM1
Rnet Joystick Controller Advanced LED Screen CJSM2 with Assignable Buttons, Bluetooth & IR	Attendant Control

### Joystick Mount

Swing Away Interchangeable	Round Style Swing Away Joystick Mount - Left
Midline Joystick Mount	Round Style Swing Away Joystick Mount - Right

### Joystick Knob Option

Goal Post	T-Bar	Mushroom
Solid Ball	Stick	Large Softball

## ACCESSORIES

Armrest Bag	Ram Tablet Mount
Side Guards	Ram Cup Holder
Tray Table	Rear View Mirror
Transfer Bar - Pair	Control 5 Switch Box - Button Controls
Ram Phone Mount	Control 5 Switch Box - Toggle Controls

### Colour Accents

Black	Blue	Green
Orange	Red	White