

Z500 Mini

POWERCHAIR SCRIPTING AND ORDERING GUIDE

Distributed by
medfab



CLIENT DETAILS

Client/User Name	
Assessor	
Funding Body	
User Weight	
Date	
Notes	

Z500 MINI BASE PACKAGE

- R-Net Expandable Electronics
- Lights
- 13" Drive Wheels
- 7" Castors
- Transit Tie Downs
- 4 Pole Motors
- 50Ah Gel Batteries
- Battery Charger
- Standard Lap Belt

SEAT SELECTION

Seat Type

ASAP Adjustable Seat Top (Width: 10" - 13" Depth: 10" - 16")

ASAP Adjustable Seat Top (Width: 12" - 16" Depth: 12" - 19")

SEAT SIZE SELECTION

Seat Width:

Seat Depth:

SEAT TO FLOOR HEIGHT SELECTION

17 17.5 18 18.5 19

Other:

POWERED SEATING OPTIONS

Tilt 50°

Elevate 12"

LEGREST OPTIONS

Hanger Options

Swing Away 90° Hangers

Swing Away 80° Hangers

Swing Away 70° Hangers

Manual Elevating Hangers

Fixed Centre Mount Footplate

ARM SUPPORT OPTIONS

Arm Support Assembly

Dual Post Flip
Back Armrests

Single Post
Armrests

Cantilever
Armrests

Armrest Pads

Desk Length
Arm Pad

Full Length
Arm Pads

Spex Comfi
Wide Armrests

JOYSTICK OPTIONS

Joystick

Rnet Joystick Controller
Standard LED

Rnet Joystick Controller
LCD Screen CJSM1

Rnet Joystick Controller
Advanced LED Screen
CJSM2 with Assignable
Buttons, Bluetooth & IR

Attendant Controller

Joystick Mount

Fixed Joystick Controller

Swing Away Joystick
Mount - Left

Midline Joystick Mount

Swing Away Joystick
Mount - Right

Joystick Knob Option

Goal Post

T-Bar

Mushroom

Solid Ball

Stick

Large Softball

ACCESSORIES

Armrest Bag

Ram Tablet Mount

Tray Table

Ram Cup Holder

Remote Stop

Rear View Mirror

Ram Phone Mount

Control 5 Switch Box -
Button Controls

Manual Recline

Control 5 Switch Box -
Toggle Controls

Folding Back Canes

Stroller Handle

Colour Accents

Black

Blue

Green

Orange

Red

White