

Acheeva Learning Station

The Acheeva is most commonly found in schools and rehabs, promoting safe supported lying that enables students with complex and continuing healthcare needs to be comfortably and easily included in classroom activities. The versatile learning station provides full support in supine, prone or side-lying to protect body shape and respiratory function.

Safe and Supported Activity Chair



Full support in supine, prone or side-lying to protect body shape and respiratory function

Locking side rails and 3 adjustable safety belts means it's very secure

Large locking castors on all units make them easy to manoeuvre

3 way power: Height, tilt and back angle. Making it easy to find a comfortable position



Inbuilt Inclinometers for easy re-set



Height adjustment (60-100cm) for both Student and Graduate



Tilt range from -10 to +50 degrees



Recline range 20 degrees in either direction



ACCESSORIES

Acheeva Airmantle



Student: 7460-8820-000

Graduate: 7460-8830-000

Support Pack

Student: 7460-3150-000

Graduate: 7460-3152-000



Bumper Accessories

Small 7460-3501-000

Long 7460-3502-000

Chunky 7460-3503-000




Side Table

7460-1850-000



Prone Table

7460-1750-000



Size	Code	Overall Length	Overall Width	Mattress Length	Mattress Width	Power Tilt in Space	Power Back Support	Unloaded Unit Weight	Safe Working Load
Student	8460-0020-000	165cm	70cm	155cm	60cm	+/-20°	+50°/-10°	73kg	100kg
Graduate	8460-0030-000	218cm		208cm				75kg	120kg