


Bath Corner Chair

Compact, lightweight and adjustable, the Smirthwaite Bath Corner Seat is ideal for children requiring a little extra help with sitting balance in the bath.

Balanced in the Bath

Upright Seating

Basic upright support option for low to moderate level needs.



Compact Design

The small footprint is perfect for most baths, and makes cleaning and storage a breeze.

Adjustable Support

Both the medial thigh support and the pelvic support can be adjusted to suit each user.

Tailored Fit


Removable PU pads make it simple to tailor both the fit, and the feel, of the seat.

Ultimate Stability

Suction cups on the seat enhance stability.



Designed to deliver low-level postural support, the removable padding makes it a breeze to tailor comfort and fit, while the angled back provides extra safety for user and carer alike. The pommel and pelvic strap can also be adjusted to suit each child's needs.



The compact design makes the Bath Corner Seat quick to clean and easy to store, and means it can fit in a wide range of baths. It's size also allows it to sit low in the tub, amplifying comfort by enabling children to fully extend their legs.

Base package includes:

- Bath corner chair complete with removable padding
- Medial thigh support
- Pelvic support belt



Size	Code	Seat Height	Seat Width	Seat Depth	Seat Depth with Pads	Back Support Height	Weight Limit	Product Weight
1	6101-4501-000	75cm	34cm	26cm	24.5cm	35cm	30kg	3kg
2	6101-4502-000			31cm	29.5cm	46cm	35kg	4kg