

Combi Bathroom Chair

Designed as a starter seat for children requiring low to moderate levels of postural support, the Smirthwaite Combi toileting chair helps nurture a sense of independence and confidence.

The Gift of Independence

Beginner's potty

The moulded plastic potty is perfect for younger children and those just embarking on their toilet training.



Adjustable toileting support

Tailor seat height, seat depth and arm support positioning to each individual user's unique needs.

Advanced users' commode

The commode, with upholstered seat, quickly transforms the Combi into a seat to suit more advanced toilet users.

Hygienic frame

The anti-microbial properties of the wooden frame help reduce the risk of cross-infection.

Pelvic strap

The velcro pelvic strap enhances stability and safety.

Base package includes:

- Upholstered Seat Unit with Adjustable Armrests
- Velcro Pelvic Strap
- Plastic Commode Seat with Lid
- Moulded Potty Seat



Seat height, seat depth and arm support position can all be quickly and easily adjusted, providing optimal support for each individual user. This incredibly versatile product can even evolve from a potty for younger children to a commode seat for more advanced users, ensuring years of value.

The Combi also features a wooden frame with anti-microbial properties which ensures surfaces remain as hygienic as possible. This helps reduce the risk of cross-infection and is particularly advantageous in multi-user environments like nurseries and schools.

ACCESSORIES

Hand Rail Toggle

Size 1: 6115-7201-000

Size 2: 6115-7202-000

Size 3: 6115-7203-000



Splash Guard

Used in conjunction with Hand Rail Toggle.

Size 1: 6115-7601-000

Size 2: 6115-7602-000



Reversible Footboard

Size 1: 6115-5021-000

Size 2: 6115-5022-000

Size 3: 6115-5023-000



Size	Code	Seat Height	Seat Width	Seat Depth	Back Support Height	Weight Limit	Product Weight
1	6115-1101-000	215-29cm	29cm	20.5-29cm	32.5cm	35kg	11kg
2	6115-1102-000	275-37.5cm	32.5cm	24.5-35.5cm	37cm	45kg	13kg
3	6115-1103-000	30-40cm	36.5cm	31-38.5cm	40cm	55kg	14kg