

Indoor Seating Comparison Chart



SPECIFICATIONS	JUNI CHAIR				ZOOMI+ HIGH CHAIR	BINGO				BINGO SPEX		GRAVITY CHAIR					JENX ATOM	JENX JUNIOR+ SEAT			SYMMETRIKIT CHAIR				
	1	2	3	4		Mini	1	2	3	1	2	Sml	Med	Lrg	XLrg	XL-Wide		Size 1	Size 2	Size 3	Low	Med	High		
DESCRIPTION	Exceptional adjustment, along with a great range of optional accessories, your child's Juni can be tailored to them.				Ideal feeding chair for users who require extra postural support.	High & Low level seating for the BINGO.				High & Low level seating for the BINGO Spex.		Supportive alternative seating and foam positioning system to augment postural control for a wide range of individuals. Also used as a 'hugging' chair for individuals who benefit from this type of sensory input.					A modular seating system for users 6 months-5 years.	Comfortable, contoured and stable modular seating system, with outstanding levels of adjustment and a wide variety of support options.			Tilt-in-space chair provides the solution for time out of regular positioning systems. Adaptable postural management for sitting and lying.				
APPROX AGE OF USER*	18mths-12 years				4 mths-adult	4mths-4yrs	1-5yrs	3-12yrs	8-15yrs	2-5yrs	3-7yrs	2-6yrs	6-10yrs	10-13yrs	13+ yrs	Adult	6mths-5yrs	18mths-5yrs	3-8yrs	5-12yrs	6-12 yrs	10-15 yrs	14-Adult		
USER WEIGHT (kg)	30	40	50	60	85	25	50	50	50	50	50	50	50	100		25		50				140			
SEAT WIDTH (cm)	24	28	32.5	36.5	43	13-25	18-30	23-35	23-35	17-28	22-33	28 front 21 back	30 front 23 back	37 front 27 back	39 front 30 back	50 front 41 back	18-27		15-37.5		17-32 cm Narrow	31.5-45.5 cm Standard	45-60 cm Wide		
SEAT DEPTH (cm)	18-24	22-30	26-36	32-42	15-26	12-27	16-34	22-40	27-45	18-30	23-38	24	30	38	43	43	9-34		21-30.5	24-35.5	26-44		31.5 - 52		
BACK SUPPORT HEIGHT (cm)	26	29	32.5	40	27	37-51	45-63	54-75	62-83	36-44	40-48	60 total back	69 total back	74 total back	81 total back	81 total back	27-38		31.5-51	34.5-51	41-57.5		Total backrest: 75-95 cm		
FOOT SUPPORT HEIGHT (cm)		x			26-60	16-34	16-37	16-40	16-40	25-40	25-45	20 to floor	23 to floor	26 to floor	29 to floor	29 to floor	11.5-36		19.5-27	23-35	25-40.5	12-25	24-56	41-52	
TILT IN SPACE		x			x		45° posterior					Use wedge underneath					10° anterior-30° posterior			10° anterior-30° posterior			-20°- +30°	-10° - +30°	
RECLINE		x			x		70°								x			-15° - +15°		-15° - +15°				✓	
GROWTH ADJUST SEAT		✓			✓	✓	✓	✓	✓	✓	✓						x	✓		✓			x		
SEAT HEIGHT (cm)	16 - 22	20 - 28	24-34	30-40	40-65	27-61	27-61	27-61	27-61	27-61	27-61						N/A	28.5-75	29-76	33-76	36.5-76	32	58	71	
HI LO OPERATION		x			x		Hand Operated Gas Assist											N/A	Hydraulic Foot Pump			Hydraulic Foot Pump			N/A
TRAY		✓			✓	✓	✓	✓	✓	✓	✓						✓	✓		✓			Mobile Table		
POSTURAL SUPPORT LEVEL	Mild/Moderate				Mild/Moderate	Mild/Moderate/Complex						Mild/Moderate					Mild/Moderate/Complex			Mild/Moderate/Complex			Mild/Moderate/Complex		

* Indicator only, for accurate product selection client measurements are the best guide.

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.