

Hoggi Cruiser Walker

Delivering perfect support for users from 18-months to 16 years' old, the Hoggi Cruiser Walker really is the ultimate anterior walking frame.

For a range of reasons, some children simply aren't suited to posterior walking aids. The Hoggi Cruiser Walker caters to those users, providing the stability and maneuverability they need to navigate the world with confidence.

Cruising with Confidence!



The large wheels are great for both indoor and outdoor use.

Tailor the height and width of the handles in accordance with your child's unique needs.

Rear-wheel drag brakes enhance safety and are standard with all Hoggi Cruiser Walkers.

Compact design and an innovative folding mechanism make it a breeze to store.

A wide range of accessories are also available, including swivel locks, storage baskets and parking brakes.

Package Includes:

- Anterior walking frame
- Large diameter wheels
- Height adjustable and width adjustable handlebar
- Lightweight frame complete with folding mechanism and upholstered sling seat
- Drag brakes for rear wheels included as standard



ACCESSORIES

Driving and parking brakes – size 1
5660-8450-000



Driving and parking brakes – size 2
5660-8452-000



Storage basket – size 1
5660-8011-000



Storage basket – size 2
5660-8012-000



Swivel lock
2880-8440-000



Available in multiple sizes, the Hoggi Cruiser Walker has large diameter wheels and height/width adjustable handles. The lightweight frame, which includes an upholstered sling seat, features a folding mechanism that makes it more compact and easier to store than any other anterior walker on the market.



Size	Code	Colour	Handle Height	Handle Distance	Seat Height	Total Width	Total Length	Folding Size	Frame Weight	Max User Weight
1	5660-0001-060	Pastel Blue	32-59 cm	22-36 cm	33 cm	51 cm	55 cm	50x55x27 cm	5.6 kg	50 kg
2	5660-0002-054	Pastel Green	44-83 cm	36-50 cm	45 cm	65 cm	58 cm	68x58x29 cm	5.9 kg	65 kg