

Hoggi CRUISER



ANTERIOR WALKER

SCRIPTING AND ORDERING GUIDE

Distributed by
medifab

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		



Measurements	Size 1 (Pastel Blue)	Size 2 (Pastel Green)
Order Code	5660-0001-060	5660-0002-054

Mobility Base with Swivelling Wheels		
Handle height (cm)	32-59	44-83
Handle distance (cm)	22-36	36-50
Seat height (cm)	33	45
Total width (cm)	51	65
Total length (cm)	55	58
Folded size (cm)	50 x 55 x 27	68 x 58 x 29
Frame weight (kg)	5.6	5.9
Load capacity (kg)	50	65

Package Includes:

- Anterior walking frame with large diameter wheels
- Height adjustable and width adjustable handlebar
- Lightweight frame complete with folding mechanism and upholstered sling seat
- Drag brakes for rear wheels

Accessories

Driving and parking brakes		Storage Basket	
Size 1 5660-8450-000		Size 1 5660-8011-000	
Size 2 5660-8452-000		Size 2 5660-8012-000	

Swivel Lock	
2880-8440-000	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: CRUISER_20240722