

Symmetrikit Tiltrite Chair

Rest, work, play. An essential part of any 24 hour postural management plan, this tilt-in-space chair provides the solution for time out of regular positioning systems such as wheelchairs or task seating. The Symmetrikit Tiltrite Chair can be quickly repositioned to provide support in functional upright positions or for relaxing and lying.

Adaptable postural management for sitting & lying



Push handles for ease of moving

Washable blue suede covers with CoolOver™ liner

Fingerform™ back included as standard

Locks into various positions at a **range of adjustment points**

Chairs are calibrated for easy resetting of angles

Multi-adjustable, conforming backrest with upper back which can fold forward for Kyphotic backs

Table accessory option: 3830-7190-000

Multi-position armrests fold down to facilitate side transfers

Two rear braked castors

Elevating height/angle adjustable footplate

ONE CHAIR FITS ALL



Both Kerry and Courtney are using exactly the same Medium chair. Both have had their chairs adjusted to fit perfectly. Settings can be replicated for them exactly, time and time again.

LIE FLAT



Use as a bed or for postural drainage.

PRONE WEDGE



This multifunctional chair can be used as a positioning wedge.

POWERED ASSIST

Powered control option for tilt, recline and leg elevation.




Fingerform backrest. Complex trunk positioning and support.

Optimal tilt-in-space angle can result in up to 40% reduction in peak & mean pressures compared to neutral.

Reclining / side lying also possible. Use chair head down for postural drainage or 'back to front' for prone stretch.

	Code	Approx Age Range	Seat Depth	Seat Width	Back Height	Leg Length	Seat Tilt	Foot Plate Angle	Seat-To-Floor	Max Weight
Powered	3830-1000-366	7yrs - Adult	315-52 cm	315-45.5 cm	75-95 cm	235-46 cm	10° / -30°	180°	58 cm	140 kg
Manual	3830-1500-066									