

FLUX Walker



VERSATILE POSTERIOR WALKING AID

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE

Distributed by
medifab

Client Details

Client/User Name		Notes
Date		
NDIS/Funding		

Measurements

	Size 1	Size 2	Size 3	Size 4
Standard grab rail grip height	48-61cm	57-81cm	69-87cm	82-102cm
Standard grab rail grip width	41cm	46cm	51cm	51cm
Universal grip grab rail grip height	43-68cm	51-88cm	62-94cm	75-102cm
Universal grip grab rail grip width	29-41cm	34-46cm	39-51cm	39-51cm
Total frame width	64cm	69cm	74cm	78cm
Total frame length	65cm	80cm	88cm	91cm
Total load capacity	60kg	75kg	75kg	100kg
Front / rear wheel diameter	15cm / 20cm	15cm / 20cm	15cm / 20cm	20cm / 30cm

Package Options

Package includes:	Size 1	Size 2	Size 3	Size 4
<ul style="list-style-type: none"> • Standard Handlebar • Front Swivel Wheels • Polyurethan Tyres • Bumper Wheels • Friction Brakes • Anti-Reverse Mechanism 	<p>Luminous Red 5641-0001-000</p> 	<p>Sky Blue 5642-0001-000</p>  <p>(Shown with universal handlebar)</p>	<p>Gentian Blue 5643-0001-000</p>  <p>(Shown with universal handlebar)</p>	<p>Grey Aluminium 5644-0001-000</p> 

Accessories

Anti-Tippers without Castors

Size 1-3

5640-8421-000



Size 4

5640-8424-000



Anti-Tippers with Castors

Size 1-3

5640-8425-000



Size 4

5640-8428-000



Fold Up Seat

Size 1-2

5640-0201-000

Size 3-4

5640-0204-000



Hip Pad

Multi-Adjustable

5640-1900-000



Universal Handlebar

Size 1

5641-5751-000

Size 2

5642-5752-000

Size 3

5643-5753-000

Size 4

5644-5754-000



Driving Brake

Size 1-3

Includes wheel lock.

5640-8451-000

Size 4

Includes wheel lock.

5640-8454-000





Swivel Lock

To hold castors straight when brisk walking.



2880-8440-000




Accessories Continued

Dynamic Hip Control (Hip Belt Required)		Hip Belt (Dynamic Hip Control Required)	
Suspension, soft 5640-1952-000		Small Hip Size: 53-83cm 5640-6201-000	
Suspension, hard 5640-1953-000		Large Hip Size: 65-115cm 5640-6202-000	
Kit for orthoses 5640-1951-000			

Forearm Supports – 2 Axial Adjust

Small 14 x 9.5cm 5640-5500-000		Standard 22 x 10.5cm 5640-5501-000	
---	---	---	---

Hand Disc

5640-5550-000	
---------------	---

Notes

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Flux_20241024_AU