

# Jenx Atom



## MODULAR SEATING SYSTEM

NEW ZEALAND  
SCRIPTING AND ORDERING GUIDE

Distributed by  
**medifab**

# Client Details

Client/User Name:

Date:


NDIS/Funding:

Measurements	User Measurements
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	

## Product Specifications

Average age	6 Months-5 Years
Weight capacity	25kg
Seat width	18cm-27cm
Seat depth	9cm-34cm
Chest support width	20cm-28cm
Seat to top of backpad height	27cm-38cm
Foot support angle	-20° to 15°
Foot support height	11cm-36cm
Back support angle	-75° to 105°
Tilt range	30° posterior, 10° anterior

## Notes



Therapists Normally Choose

## Seat Unit/Colour Choice

## Base Options


Red	Pink	Y-Base
<p><b>Includes:</b></p> <ul style="list-style-type: none"> <li>Padded seat with memory foam</li> <li>Shaped back pad</li> <li>Arm supports</li> <li>Nursery tray</li> <li>Adjustable flip away hip support pads</li> <li>Push handle</li> <li>Upper back pad</li> <li>Flip up angle adjustable footplate</li> </ul> <p>3431-1001-022</p>	<p><b>Includes:</b></p> <ul style="list-style-type: none"> <li>Padded seat with memory foam</li> <li>Shaped back pad</li> <li>Adjustable flip away hip support pads</li> <li>Push handle</li> <li>Upper back pad</li> <li>Arm supports</li> <li>Nursery tray</li> <li>Shoulder protractors</li> <li>High rise flexi thoracic pads</li> <li>Thoracic pad covers</li> <li>2 Link multigrip head support</li> <li>Flip up angle adjustable footplate</li> </ul> <p>3431-0010-028</p>	<p>Black Y-base with lockable swivel castors            Floor to seat height: 30-75cm            Tilt range: 30° posterior 10° anterior            3615-8740-000</p> <p></p>

# Head Supports

1 Link Multigrip	2 Link Multigrip	
Single link wings for some degree of lateral support. 33cm (W) x 5cm (D) x 7.5cm (H) <b>For Red Atom:</b> 3431-4601-022 <b>For Pink Atom:</b> 3431-4601-028	2 Links in each wing for extra support and positioning.  47.5cm (W) x 5cm (D) x 12cm (H) <b>For Red Atom:</b> 3431-4602-022 <b>For Pink Atom:</b> 3431-4602-028	
Wipeable Cover (Optional) 5356-4671-000	Wipeable Cover (Optional) 5356-4672-000	

# Back Pad Options





Lumbar Pad	Back Pad Infill
Provides additional lumbar support. 3431-2501-000	
	A 50mm back pad to reduce the seat depth for smaller children. 3431-2370-000
	

Washable Cover	
For Seat and Back Support. Comes in Red. 3431-2470-022	

# Shoulder & Trunk Supports

Shoulder Protractors	High-Rise Short Thoracic Pad	High-Rise Flexible Thoracic Pads
Flip away when not required. Independently adjustable for width and angle. 13cm (H) x 15.5cm (D) 3431-2750-022	Short, slim trunk support pads with single-handed flip away action. No strap. 3430-3151-000	Slimline and super flexible trunk support pads with straps and single-handed flip away action.  8cm (H) x 18cm (D) 3431-3251-022
		
	Cover (Optional) 3615-3171-022	Cover (Required) 3431-3271-022

# Harnesses & Hip Belt

Waistcoat Harness	Wonderseat Harness	3 Point Lap Strap	4-Point Hip Belt
<p>Provides additional trunk support to children with low tone. 3615-6701-000</p>	<p>Size 0, Pad height 22cm Padded Centrepoint Shoulder Harness. 1409-5660-017</p>	<p>Dual lapstrap with additional safety of a groin harness. 3431-6251-000</p>	<p>Wonderseat Padded Hip Belt, Pad width 24cm. 1406-5420-027</p>
			



# Leg Supports

Leg Positioning Pads	Long Leg Sitting Board	Abduction Block	Abduction Block
<p>Ideal for controlling windsweeping or creating an abducted seating position. Required fixed angle foot support. 3431-1751-000</p>	<p>For very young children or those needing hamstring stretch. Used in place of foot support. Comes in red. 3431-0281-022</p>	<p><b>Size 0</b> Length: 9cm Height above Seat: 8.5cm Overall height: 14.5cm Width Front: 9cm Width Rear: 7cm 3615-1650-000</p>	<p><b>Size 1</b> Length: 18cm Height above Seat: 7.5cm Overall height: 15.5cm Width Front: 10cm Width Rear: 8.5cm 3615-1651-000</p>
			

# Foot Supports

Foot Support Pad
<p>Attaches to foot support for softer surface. Can only be used on the flip-up foot support. 3cm thickness. 3431-5701-000</p>


# Foot Supports (cont.)

Sandals	Sandal Raising Blocks	Ankle Snugs
<p>Ideal for positioning feet where a firmer hold on the foot is required. 8cm (W) x 17cm (L) 5356-6081-000</p> 	<p>Easily accommodate leg length discrepancies by giving sandals extra height where needed. 4 x 1.5cm Blocks 5356-6091-000</p> 	<p>Features unique magnetic buckle fastening mechanism. Sandals required. Inner Circumference: 12cm-22.5cm 5356-6131-000</p> 



# Arm Supports

Covers for Arm Supports	
<p>(Red) Provide soft top for arm supports when tray is not used. 3431-5510-022</p>	

# Tray & Accessories

Nursery Tray Infill	
<p>Designed to work with the Nursery Tray to reduce the cut-out for smaller children. 3431-7182-000</p>	

Nursery Tray Padding	Elbow Blocks
<p>Half Padding 3431-7442-000</p>	<p>Comes as a pair. Mounts to Nursery Tray corners for elbow support. 9.5cm (H) x 10cm (W) 3430-7350-000</p>
<p>Full Padding 3431-7452-000</p>	

Activity Frame	Grab Bar
<p>Great for encouraging head control and hand/eye co-ordination. Requires tray. 54cm (W) x 16cm (D) x 35cm (H) 3615-7475-000</p>	<p>Nursery Tray required. Bar dimensions: 54cm (W) x 9.5cm (D) x 7cm (H) 3615-7200-000</p>
	

# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Jenx Atom\_20241030\_NZ