

# Zoomi+



## INDOOR SEATING

AUSTRALIAN  
SCRIPTING AND ORDERING GUIDE

Distributed by  
**medifab**

# Client Details

|                   |  |
|-------------------|--|
| Client/User Name: |  |
| Date:             |  |
| NDIS/Funding:     |  |

Get it Faster

9999-ASMB-002

## Medifab Assembly Service for Easy Delivery

Medifab offers our popular assembly service on this product, meaning the product arrives ready to be delivered, straight from the box, saving time and speeding up the delivery process.

| Order Code | Size 1        | Notes |
|------------|---------------|-------|
| Black      | 3212-1001-099 |       |
| Blue       | 3212-1001-060 |       |
| Red        | 3212-1001-022 |       |
| Lime       | 3212-1001-054 |       |

| Seat Unit Measurements          |   |
|---------------------------------|---|
| Seat height (cm)                | 81  |
| Seat width (cm)                 | 43  |
| Seat depth (cm)                 | 15-26<br>(increased adjustment possible with accessories) |
| Foot support height (cm)        | 26-60   |
| Back support height (cm)        | 27  |
| Load capacity of seat unit (kg) | 85  |
| Product weight (kg)             | 16  |

1 Ranges given are recommendations only - further adjustment range is possible.  
 2 Foot board can be removed, leaving option of foot flat to floor.

### Base package includes:

- Seat Unit
- Contoured Padded Seat
- High Back Support
- Arm Supports
- 5 Point Harness
- Foot Support
- Quick Release Tray
- Mobile Ski's



# Accessories

## Transfer Step

Add on accessory to facilitate a standing transfer.

**Black**

3212-5901-099

**Blue**

3212-5901-060

**Red**

3212-5901-022

**Lime**

3212-5901-054



## Extended Slotted Footboard

**Black**

3212-5021-099

**Blue**

3212-5021-060

**Red**

3212-5021-022

**Lime**

3212-5021-054



## Extended Contoured Padded Seat

Upgraded seat depth for bigger users, addition 10cm seat depth.

**Black**

3212-8991-099

**Blue**

3212-8991-060

**Red**

3212-8991-022

**Lime**

3212-8991-054



## Wonderseat Padded Centrepoint Shoulder Harness

**Size 0**

Pad height 22cm

1409-5660-017

**Size 1**

Pad height 28cm

1409-5661-017



## Spex 2-Point Hip Belt

**X-small**

2 point slim line, dual pull with neoprene sleeve padding.

1405-6220-027



# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Zoomin+, 20241030\_AU