

# Standing Frame Comparison Chart



SPECIFICATIONS	Prone only			Supine or prone			Supine only		
	JENX MONKEY	JENX PRONE STANDER		JENX MULTISTANDER		JENX STANDZ	JENX SUPINE		
		Size 1	Size 2	Size 1	Size 2		Size 1	Size 2	Size 3
DESCRIPTION	Gentle introduction to standing with easy to adjust positioning.	Dedicated to assisted prone standing with the child supported at the front		Highly supportive standing frame in prone, supine, or upright.		Abducted standing systems offering prone and supine abducted standing.	Robust and easy to use standing solution.		
APPROX AGE OF USER*	9mths-4yrs	3-7yrs	6-14yrs	9mths-6yrs	3-13yrs	1-9yrs	9mths-4yrs	3-11yrs	10-Adult
USER WEIGHT (kg)	25	55kg	70kg	35	60	45	25	50	100
MAX USER HEIGHT (cm)	-	-	-	116	175	138	105	145	190
FOOT TO AXILLA (cm)	(Top of chest pad) 51-75	66-100	82-127	Prone 49-86 Supine 57-86	Prone 70-117 Supine 77-117	58-100	50-107	90-147	154-164
CHEST WIDTH (cm)	15-23	19-29	23-35	15-26	17-32	18-30	13-31	13-31	22-44
HIP WIDTH (cm)	15-23	19-29	23-35	12-26	17-32	18-30	13-31	13-31	22-44
SUPINE TILT	N/A	N/A		90-10°		90°-0°	0-80°		
PRONE TILT	20°	0-40°		10-90°		0-90°	N/A		
MULTI ADJUST KNEE SUPPORT	✗	✗		✓		✓	✓		
MULTI ADJUST FOOT SUPPORTS	✓	✓		✓		✓	✓		
TRAY TABLE	✓	✓		✓		✓	✓		
HEAD SUPPORT	✗	✗		✓ (Supine)		✓ (Supine)	✓		
POSTURAL SUPPORT LEVEL	Mild/Mod.	Mild/Mod.		Mild/Mod.		Mild/Mod./Complex	Mod./Complex		
FRAME FOOTPRINT (cm)	57x88	90x65	99.5x65	56x71	65x95	63.5x90	57x66	65x95	65x95
LEG ABDUCTION RANGE**	~10°	~10°		~10°		30°	Neutral		
OPERATION	Gas Assist	Gas Assist		Gas Assist (size 2 electric option)		Gas Assist	Battery/Electric		
KNEE CONTRACTURE ACCOMMODATE	Nil	Nil		Supine: Extensive Prone: Mild		Supine: Extensive Prone: Mild	Mild	Mild	Extensive

\* Indicator only, for accurate product selection client measurements are the best guide.  
\*\* Measured from vertical (0°) leg alignment to abducted angle

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.