

Vehicle Seating Comparison Chart



SPECIFICATIONS	HERO-NXT		STARLIGHT-NXT		HERO-XL	CARROT 3000	CARROT XL	KIDSFLEX		
	With Std Cushion	Without Std Cushion	Standard	With XL Kit				1	2	3
DESCRIPTION	Swivelling special purpose car seat for forward or rear facing.		Excellent posture support option with tilt and swivel functions.		Innovative extended rear facing support for those that need it. Up to 125cm tall.	Extension accessories and multi-adjustable features allow longer use.	Larger version of Carrot 3000 for older or heavier children and small adults.	Fully adjustable specialised seating with excellent posture support options.		
USER WEIGHT RANGE (KG)	0 - 18		9 - 36		36 (max)	13.6 - 49	36 - 75	0 - 22	0 - 36	0 - 75
SEAT WIDTH (CM)	21	25	30		24	31	51	31	35	40
SEAT DEPTH (CM)	25	28	26	32	25 - 34	25 - 40	38 or 42	28 - 33	34 - 39	40 - 46
SHOULDER HEIGHT (CM)	24 - 34	27 - 35	32 - 42	44 - 53	28 - 45	32 - 57	49.5 - 62	37 - 51	41 - 57	47 - 65
APPROX PRODUCT WEIGHT (KG)	16		13.5		9	8.7	10.5	8.2	8.9	12.3
TILT IN SPACE	✓		✓		✓	Using wedge	Using wedge	✓		
SWIVEL BASE	✓		✓		✗	✓	✗	✓		
GROWTH ADJUSTABLE	✓		✓		✓	✓	✓	✓		
FOOT SUPPORT	✗		✗		✗	✓	✗	✓		
TRAY/FOREARM SUPPORT	✓		✓		✓	✓	✓	✓		
POSTURAL SUPPORT LEVEL	Mild/Mod/Complex		Mild/Mod/Complex		Mild/Mod/Complex	Mild/Mod/Complex	Mild/Mod/Complex	Mild/Mod/Complex		
IS ISOFIX REQUIRED IN VEHICLE?	Essential		Optional		No	Only required for Swivel Accessory	Optional	Optional		
REAR FACING CAPACITY	✓		✗		✓ (Rear Facing Only)	✗	✗	✗		
SUPPORT/PAD OPTIONS	✓		✓		✓	✓	✓	✓		
STANDARD/REGULATION COMPLIANCE	ECE r129		ECE r44		ECE r129	ECE R44 & FMVSS 213	Medical Device Regulation (EU) 2017/745	Medical Device Regulation (EU) 2017/745		
ANTI-ESCAPE OPTIONS	✗		✗		✗	✓	✓	✗		
RECLINE	✗		✗		✗	✓	✓	✓		

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.