



Walking Aids Comparison Chart



SPECIFICATIONS	GRILLO				FLUX WALKER			
	Mini	Small	Medium	Large	Size 1	Size 2	Size 3	Size 4
DESCRIPTION	Versatile aid that can be either posterior or anterior. Interchangeable supports make it suitable for wide range of users.				Lightweight yet sturdy walking aid with unique design that promotes natural walking patterns.			
APPROX AGE OF USER*	1-6yrs	3-8yrs	8-18yrs	Adults	1-6yrs	3-8yrs	8-18yrs	Adults
MAX USER WEIGHT (kg)	35	40	90	110	60	75	75	100
TOTAL BODY HEIGHT (cm)	75-90	90-135	10-160	190	65-115	95-130	110-170	Adult
FOOT - AXILLA (cm)	52-77	10-160	10-160	81-143	N/A			
FOOT - HIP (cm)	41-63	10-160	10-160	81-143	N/A			
PRODUCT WEIGHT (kg)	13.5	Posterior: 15.7 Anterior: 16.2	Posterior: 17.7 Anterior: 19.5	Posterior: 22 Anterior: 23	6.5	6.8	7.1	9.1
ANTI-REVERSE MECHANISM	✓				✓			
CHEST SUPPORT	✓				✗			
HIP POSITIONER	✓				✓			
FOREARM SUPPORTS	✓				✓			
LEG DIVIDER/SEPARATOR	✓				✗			
FOOTPRINT (cm)	57x72	65x75	76x91	80x98	64x65	69x81	74x88	78x91
FOOT TO HAND GRIP HEIGHT(cm)	Variable				48-61 (Std..) 43-68 (uni.)	57-81 (Std..) 51-88 (uni.)	69-87 (Std..) 62-94 (uni.)	82-102 (Std..) 75-102 (uni.)
SADDLE SEAT	✓				✗			
SLING SEAT	✓				✓			
ATTENDANT HANDLE	✓				✗			
FRICTION/SPEED BREAK	✓				✓			
FOLDING SEAT OPTION	Posterior Only				✓			

* Indicator only, for accurate product selection client measurements are the best guide.

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.