

Indoor Seating Comparison Chart



SPECIFICATIONS	JUNI CHAIR				ZOOMI+ HIGH CHAIR	BINGO				BINGO SPEX		GRAVITY CHAIR					JENX ATOM	JENX JUNIOR+ SEAT		
	1	2	3	4		Mini	1	2	3	1	2	Sml	Med	Lrg	XLrg	XL-Wide		Size 1	Size 2	Size 3
DESCRIPTION	Exceptional adjustment, along with a great range of optional accessories, your child's Juni can be tailored to them.				Ideal feeding chair for users who require extra postural support.	High & Low level seating for the BINGO.				High & Low level seating for the BINGO Spex.		Supportive alternative seating and foam positioning system to augment postural control for a wide range of individuals. Also used as a 'hugging' chair for individuals who benefit from this type of sensory input.					A modular seating system for users 6 months-5 years.	Comfortable, contoured and stable modular seating system, with outstanding levels of adjustment and a wide variety of support options.		
APPROX AGE OF USER*	18mths-12 years				4 mnths-adult	4mths-4yrs	1-5yrs	3-12yrs	8-15yrs	2-5yrs	3-7yrs	2-6yrs	6-10yrs	10-13yrs	13+ yrs	Adult	6mths-5yrs	18mths-5yrs	3-8yrs	5-12yrs
USER WEIGHT (kg)	30	40	50	60	85	25	50	50	50	50	50	50	50	100		25		50		
SEAT WIDTH (cm)	24	28	32.5	36.5	43	13-25	18-30	23-35	23-35	17-28	22-33	28 front 21 back	30 front 23 back	37 front 27 back	39 front 30 back	50 front 41 back	18-27		15-37.5	
SEAT DEPTH (cm)	18-24	22-30	26-36	32-42	15-26	12-27	16-34	22-40	27-45	18-30	23-38	24	30	38	43	43	9-34	21-30.5	24-35.5	26-44
BACK SUPPORT HEIGHT (cm)	26	29	32.5	40	27	37-51	45-63	54-75	62-83	36-44	40-48	60 total back	69 total back	74 total back	81 total back	81 total back	27-38	31.5-51	34.5-51	41-57.5
FOOT SUPPORT HEIGHT (cm)		x			26-60	16-34	16-37	16-40	16-40	25-40	25-45	20 to floor	23 to floor	26 to floor	29 to floor	29 to floor	16.5-32.5 (Fixed Angle Foot Support) 11.5-36 (Flip Up Foot Support)	19.5-27	23-35	25-40.5
TILT IN SPACE		x			x		45° posterior					Use wedge underneath					Wooden Base: 20° posterior Y Base: 10° anterior-30° posterior			
RECLINE		x			x		70°								x		-15° - +15°			-15° - +15°
GROWTH ADJUST SEAT		✓			✓	✓	✓	✓	✓	✓	✓						✓		✓	
SEAT HEIGHT (cm)	16 - 22	20 - 28	24-34	30-40	40-65	27-61	27-61	27-61	27-61	27-61	27-61				N/A		28.5-75	29-76	33-76	36.5-76
HI LO OPERATION		x			x		Hand Operated Gas Assist											Hydraulic Foot Pump		Hydraulic Foot Pump
TRAY		✓			✓	✓	✓	✓	✓	✓	✓				✓		✓		✓	
POSTURAL SUPPORT LEVEL	Mild/Moderate				Mild/Moderate	Mild/Moderate/Complex					Mild/Moderate					Mild/Moderate/Complex	Mild/Moderate/Complex			

* Indicator only, for accurate product selection client measurements are the best guide.

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.